

# OFIT<sup>®</sup> Femoral Stem System

## Restore Activity

This documentation is intended for health care professionals and the b-ONE ORTHO sales force.  
Distribution to any other recipient is prohibited.

Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. The treating surgeon is responsible for determining the appropriate treatment, technique(s), and product(s) for each individual patient. b-ONE ORTHO does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

Contact your b-ONE representative if you have questions about the availability of b-ONE products in your area.

All content herein is protected by copyright, trademarks and other intellectual property rights owned or licensed to b-ONE ORTHO Corp. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of b-ONE ORTHO.

**Please refer to the following for compatible surgical techniques:**

Tapered Wedge Femoral Stem System Surgical Technique  
Direct Anterior Approach Surgical Technique


**BIOLOX<sup>®</sup> delta is a trademark of CeramTec AG.**

BJLIT-00012-V01 OFIT Surgical Technique EN  
©2023 b-ONE ORTHO



**b-ONE ORTHO Corp.**  
3 Wing Drive, Suite 259  
Cedar Knolls, NJ 07927  
**USA**  
[www.b1.co](http://www.b1.co)

**b-ONE**  
O R T H O



# TABLE OF CONTENTS

**01** | System Introduction

**03** | Indications and  
Contraindications

**05** | Pre-Operative Templating

**07** | Femoral Neck Resection

**08** | Preparing the Femoral Canal

**09** | Broaching

**11** | Trial Reduction

**12** | Head Trial

**13** | Cement Plug Insertion

**14** | Implanting The Stem

**15** | Implanting The Stem

**16** | Implant Removal

# SYSTEM INTRODUCTION



The OFIT® Femoral Hip System is a bone compacting femoral stem that is intended for Cementless and Cemented fixation. This hip stem is evolved from a traditional bone compacting design to optimize fit throughout a wide demographic of patient anatomies with a size-specific medial curvature geometry and fixed neck length design across the stem families.

The versatile stem design incorporates a reduced proximal lateral shoulder, a lateral distal relief, and an overall shorter stem length compared to traditional designs to allow ease of insertion in various surgical approaches.

The OFIT® Femoral Hip System is a broach-only system comprised of simple, intuitive instruments designed to optimize operative workflow and efficiency. Instrumentation carriers are designed to seamlessly stack into Steril Containers™ by Aesculap® or can be processed individually to accommodate hospital protocol or preference.

## The OFIT® Hip System Includes

- Three Stem Options in One Family  
OFIT® family offers OFIT® HA, Basic and Cemented options to meet various clinical and economical needs. OFIT® instrumentation are compatible amongst all three options.
- Fully HA Coated  
Fully HA coated stem induces rapid osteointegration which improves early and long term stability of the stem.
- Grit Blasted Surface for Osteointegration  
Grit blast in various levels of roughness to enhance overall bone on-growth performance.
- High Nitrogen Stainless Steel  
High nitrogen stainless steel is a nickel free solution that offers lightweight, durability and biocompatibility compared to other commonly used materials.

# INDICATIONS AND CONTRAINDICATIONS

## INDICATIONS

The b-ONE® Total Hip System is intended for primary or revision total hip replacement in skeletally mature patients with a severely disabled hip joint and/or hip damage due to the following conditions:

Osteoarthritis, traumatic arthritis, avascular necrosis of the femoral head, noninflammatory degenerative joint disease (NIDJD), slipped capital epiphysis, fused hip, fracture of the pelvis, and diastrophic variant. Hip components are also indicated for inflammatory degenerative joint disease including rheumatoid arthritis and congenital dysplasia; treatments of nonunion, acute traumatic fracture of the femoral head or neck; failed endoprosthesis, femoral osteotomy, or Girdlestone resection; and fracture-dislocation of the hip.

OFIT® HA and OFIT® Basic stems are intended for cementless use only.

OFIT® Cemented stainless steel stems are intended for cemented use only.


b-ONE® Total Hip System components are not intended for use with other total hip systems.

## CONTRAINDICATIONS

- Use of this system is contraindicated when there is active systemic infection, infection localized to the site of the proposed implantation, or when the patient has demonstrated allergy or foreign body sensitivity to any of the implants materials.
- Severe osteoporosis or osteopenia may prevent adequate fixation and thus preclude the use of these or any other orthopedic implants.
- Conditions that may place excessive stresses on bone and implants, such as severe obesity, or pregnancy, are relative contraindications. The decision to use these devices in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
- Use of these implants is relatively contraindicated in patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation, or lifestyle may interfere with their ability to follow postoperative restrictions and who may place undue stresses on the implant during body healing and may be at a higher risk of implant failure.
- Using a BIOLOX® delta head in combination with a prosthesis stem left in situ in a revision surgery is contraindicated. A BIOLOX® delta head must only be used with a brand-new, unused, and undamaged stem taper.



# PRE-OPERATIVE TEMPLATING



Preoperative planning supports the determination of appropriate stem style, size, level of femoral neck

Qualitative and quantitative factors such as patient bone quality, density, and morphology should be considered to select the appropriate implant system for the patient. Preoperative templating should serve only as a guide.

Preoperative templating requires quality radiographs with known and correct magnification. The desired magnification for all imaging should be 20% magnification to correspond with the b-ONE® x-ray templates,


with x-ray magnification calibration used whenever possible. Generally, proper radiographs include a single anteroposterior (A/P) radiograph of the pelvis as well as A/P and lateral radiographs of the affected hip to show the proximal one-third of the femur. A/P views with the limbs in 15 degrees of internal rotation are preferred.

The following templating technique assumes the patient has a normal, symmetrical pelvis

## Template Acetabulum

Overlay the acetabular template on the x-ray, ensuring the medial border of the cup approximates the ilioischial line, and the inferior border of the cup is at the inferior aspect of the teardrop. The cup should be positioned with an abduction angle of 40-45 degrees. Mark the center of rotation of the acetabular component.

**Check b1LIT-00002 for b-ONE Primary Acetabular Technique.**



## Assess Leg Length



First note any possible hip flexion contracture which could make the leg appear short on x-ray. Use clinical evaluation with radiographic analysis to determine intraoperative leg length management.

Beginning with the A/P of the pelvis, draw a reference along the inferior border of the ischial tuberosities, ensuring the line extends beyond the medial cortices of the femurs. Alternatively, a reference line through the inferior aspect of the teardrop landmarks

can be used. Then mark a reference point on each femur, such as the most proximal aspect of each lesser trochanter. Measure the distance between the reference line and each femoral reference point. Often, a line parallel to the reference line is drawn through each femoral reference point to assist with this measurement. The difference between the two measurements will indicate leg length discrepancy.

## Template Femur

The OFIT® femoral stem fixed neck angle of 132° with multiple neck lengths proportional to the stem size. The OFIT® template has markings that indicate the center of the femoral head for the range of head options for each femoral neck offset option. Choose the appropriate stem size that achieves mediolateral cortical engagement at the proximal two-thirds of the stem and recreates the desired leg length and offset.

The relative positioning of the head center of rotation markings on the femoral template with respect to the acetabular center of rotation previously marked on the x-ray will predict the change in leg length and offset. For example, a given head center of rotation marking superior to the acetabular center of rotation mark will lengthen the limb, while a head center of rotation inferior to the acetabular center of rotation will shorten the limb. The desired change in leg length is determined by the radiographic leg length

inequality and clinical evaluation previously determined. The predicted change in offset is also considered by comparing the relative medial/lateral position of the center of rotation markings of the femoral and acetabular components.

Mark the anticipated neck resection level. This will be used as a reference during neck resection.




# FEMORAL NECK RESECTION

The neck resection level affects the final fit and placement of the stem.

The Neck Resection Guide can aid in marking the appropriate neck resection level by placing it on the anterior/posterior aspect of the exposed femur, with the centerline aligned with the axis of the femoral canal. Care should be taken to reference the anatomic landmarks determined during preoperative templating, as well as visual inspection in relation to the lesser trochanter prior to making the cut. After the femoral resection is marked, the resection is made with an oscillating saw.

To remove the femoral head, a Modular Corkscrew is available and can be connected to the modular T-Handle or power. After the femoral head is removed, typically the acetabulum is prepared for the acetabular component.




# PREPARING THE FEMORAL CANAL

## Access the Femoral Canal

Position the leg to provide the best exposure for preparation of the femoral canal. Use the Box Osteotome and Mallet to initiate entry into the femoral canal. Ensure the orientation of the Box Osteotome reflects the desired anteversion, which is typically 10-15 degrees.

Use the starter reamer to create a pathway into the medullary canal. To minimize the risk of varus placement or under-sizing of the femoral prosthesis, remove adequate bone from the lateral aspect of the canal with the Starter Reamer, Box Osteotome, or a rongeur.



# BROACHING

Begin broaching with the smallest available broach. The broach size can be identified on the broach.

The Broach Handle is designed for easy attachment to the broach by extending and closing the lever handle. Be sure to orient the broach in the correct version and pay special attention to the varus/valgus and anterior/posterior placement of the broach. With the Mallet, deliver solid impacts to the strike plate on the Broach Handle to advance the broach.

Sequentially increase the size of the Broach until adequate fill is achieved.


There are markings on the broach to identify its size. Note that all OFIT® family is designed to be used with same broaching system.



The final Broach should sit firmly on medial and lateral cortical bone. Increased resistance to advancement and change in pitch during impaction serve as clues to achieving adequate size.

**Broach can be advanced so that the medial aspect of broach is in line with resected bone.**

Leave the final Broach in the femoral canal and remove the broach handle to proceed with calcar planning (optional) or trialing.




If calcar planning is desired, attach the Calcar Planer to power, ensuring the power setting is set to ream. Advance the Calcar Planer over the broach post, confirming alignment and stability. Power should be initiated prior to contacting bone. Slowly advance the Calcar Planer on continuous power until the stop engages the broach post and adequate bone is removed.

If the Calcar Planer cannot fully engage the broach post, remove the broach and recut the neck resection or consider the next larger broach.




# TRIAL REDUCTION



## Neck Trial


The OFIT® hip stem is designed with fixed neck angle across all stem sizes, neck length is proportional to the stem size. See the attached Implant Information table for more information.

Choose the desired neck trial and assemble onto the broach.



Size	Neck Angle	A	B	C
		Stem Length (mm)	Neck Length (mm)	Offset
8.5	132°	120	30.0	33.0
9		130	30.0	34.0
10		140	30.0	36.9
11		145	32.5	37.4
12		150	32.5	38.4
13		155	32.5	40.7
14		160	35.0	41.2
15		165	35.0	42.2
16		170	35.0	44.6
17		180	37.5	45.6
18	190	37.5	46.6	

# HEAD TRIAL



Insert the appropriate b-ONE® 12/14 taper head trial onto the neck trial. Note the femoral head offsets differ depending on the femoral head size as shown in the chart below:

HEAD DIAMETER	OFFSET			
	S	M	L	XL
22mm;22.2mm	0	2	4	N/A
28mm	-3.5	0	3.5	7
32mm	-4	0	4	7
36mm	-4	0	4	8

When performing the trial reduction, it is recommended to perform the following:

- Inspect the reduction of the femoral head in the acetabular cup. The reduction should be concentric and the appropriate amount of coverage of the femoral head achieved.
- Appropriate tissue tension should be assessed. Pulling the leg in a neutral position is important to obtain a true assessment of tissue tension.
- Assess stability through a full functional range of motion, checking any maneuvers that lead to instability. Note any acetabular osteophytes that may cause the hip to sublux out of the cup.
- Assess the leg lengths.


## CEMENT PLUG INSERTION

Prior to implanting cement plug, it is recommended to clean the femoral canal of debris and to open the interstices of the bone.

Select the appropriate size of cement plug identified during pre-operative templating to fit the distal canal. Attach it to the cement plug inserter and insert the trial cement plug to the planned depth. Check that it is firmly seated in the canal.

The planned level should be 10mm below the distal tip of the implant. Reference the Implant Information table attached at the end of the surgical technique to determine the desired depth.

Once the desired depth is reached and the Cement Plug is firmly seated in the canal, remove the Cement Plug Inserter, clean the femoral canal of the debris.



## IMPLANTING THE STEM

### Cemented Application

OFIT<sup>®</sup> Cemented Stem is used for cemented applications.

For the cemented application, femoral canal must be first cleaned and dried. Cement is then applied and pressurized to ensure the cement is firmly fixed into the trabecular bone. Continually inject cement during the period of pressurization.

The OFIT<sup>®</sup> Cemented Stem is recommended to be placed into the prepared femoral canal with the Stem Remover. Introduce the implant in line with the long axis of the femur in one slow movement. Its entry point should be lateral, close to the greater trochanter.

During stem insertion maintain thumb pressure on the cement at the medial femoral neck. Insert the stem up to the resection level. If necessary, a few light taps on the stem inserter will bring the stem to the right level.



## IMPLANTING THE STEM


Once the final construct is determined, remove all trial components from the body. Ensure that the selected stem size matches the broach and neck trial combination that was determined during the trial procedure.

### Cementless Application

OFIT® Basic and OFIT® HA Stems are used for cementless applications. Implantation technique is the same for both of these stem types, the following shows the example of OFIT HA Stem.



Threaded Insertion



Stem Impactor/Insertor

The OFIT® Basic and HA Femoral Stems can be placed into the prepared canal by hand and seated with the femoral stem impactor.

Place the prosthesis in the femoral canal, ensuring the correct version prepared during femoral preparation. Position the impactor/insertor into the impaction slot on the superior lateral aspect of the prosthesis. Impact the Femoral Impactor/insertor with the Mallet until the stem is fully seated.

**NOTE:** The Straight Stem Impactor do not lock to the prosthesis. If a threaded insertion tool is preferred, the Stem Remover can be used. Thread the Stem Remover into the threaded hole contained in the impaction slot of the prosthesis. Proceed with implantation as previously described.

## IMPLANT REMOVAL

If the stem must be removed, the impaction slot of the femoral prosthesis contains a threaded hole that mates with the Stem Remover. Clean the impaction slot and threaded hole thoroughly to ensure debris does not prevent engagement of the Stem Remover threads. Attach the Stem Remover to the Modular Axial Handle. Thread the assembled Stem Remover into the threaded hole of the prosthesis, being careful to avoid cross-threading.

Proceed to reverse impact the stem remover handle with a mallet to remove the stem.

