







7.0ChLP femoral periprosthetic plate

3.7220; 3.7221;
3.7222; 3.7223;
3.7224; 3.7225

- *IMPLANTS*
- *INSTRUMENT SET*
- *SURGICAL TECHNIQUE*



SYMBOLS DESCRIPTION

	Titanium or titanium alloy		H length [mm]
	Cobalt		Angle
	Left		available lengths
	Right		Available number of holes
	Available versions: left/right		Thickness [mm]
	Length		Scale 1:1
	Torx drive		Number of threaded holes in the shaft part of the plate
	Torx drive cannulated		Number of locking holes in the plate
	Hexagonal drive		Variable angle
	Hexagonal drive cannulated		Cortical
	Cannulated		Cancellous
	Locking		Available in sterile/ non- sterile condition
	Diameter [mm]		Refer to surgical technique
<hr/>			
	Caution - pay attention to a special procedure.		
	Perform the activity under X-Ray control.		
	Information about the next stages of a procedure.		
	Proceed to the next stage.		
	Return to the specified stage and repeat the activity.		
	Before using the product, carefully read the Instructions for Use. It contains, among others, indications, contraindications, side effects, recommendations and warnings related to the use of the product.		
	The above description is not a detailed instruction of conduct. The surgeon decides about choosing the operating procedure.		

www.chm.eu

Document No ST/80-703
Date of issue 12.06.2019
Review date P-002-18.01.2021

The manufacturer reserves the right to introduce design changes.

Updated INSTRUCTIONS FOR USE are available at the following website: ifu.chm.eu

1. INTRODUCTION	5
2. IMPLANT DESCRIPTION	6
3. SURGICAL TECHNIQUE - 7.0ChLP FEMORAL PERIPROSTHETIC PLATE 3.7220.XXX	9
3.1. PATIENT'S POSITIONING	9
3.2. SURGICAL APPROACH	9
3.3. FRACTURE REDUCTION	9
3.4. Implant selection	9
3.5. PLATE INSERTION	9
3.6. TEMPORARY PLATE STABILIZATION	10
3.7. SCREWS INSERTION	10
3.8. WOUND CLOSURE	10
4. SURGICAL PROCEDURES	11
4a. PROCEDURE OF TEMPORARY IMPLANT STABILIZATION	11
4b. PROCEDURE OF CORTICAL SELF-TAPPING SCREW 4.5 [3.1471] INSERTION	12
4c. PROCEDURE OF 7.0ChLP SELF-TAPPING SCREW 5.0 [3.5210] INSERTION	13
5. POSTOPERATIVE PROCEDURE	14
6. IMPLANT REMOVAL	14
7. CATALOGUE PAGES	15
7a. INSTRUMENT SET	15
7b. Implants	17
7c. Screws	19

1. INTRODUCTION

This surgical technique applies to 7.0ChLP locked plating system used for stabilization of proximal femur fractures. The plates are a part of the ChLP locked plating system developed by **ChM**. The presented range of implants is made of materials in accordance with ISO 5832 standards.

The system includes:

- implants (*plates and screws*),
- instrument set used in the surgery,
- surgical technique.

Indications

- trochanteric osteotomies,
- trochanteric fractures,
- periprosthetic femur fractures.

Plate selection and shaping

The plates are available in various lengths and for left and right limb separately. This allows for optimal selection of the implant to the fracture type. Shaping of the plates in their epiphyseal part is not allowed.









Before using the product, carefully read the Instructions for Use. It contains, among others, indications, contraindications, side effects, recommendations and warnings related to the use of the product.




The above description is not a detailed instruction of conduct. The surgeon decides about choosing the operating procedure.

2. IMPLANT DESCRIPTION

Femoral periprosthetic plates are a part of 7.0ChLP system. This system includes also compatible locking screws. To facilitate their identification, both titanium plate and screws are blue anodized.


Femoral periprosthetic plates left		Femoral periprosthetic plates right	
Short epiphysis		Short epiphysis	
			
Long epiphysis		Long epiphysis	
			
Trochanteric periprosthetic plates			
Short epiphysis		Long epiphysis	
			

7.0ChLP femoral periprosthetic plate




• New locking hole design:

- the screws heads do not protrude above the surface of the plate what significantly reduces the irritation of periimplant tissues,
- cortical screws use possible.



• K-wire holes in shaft and epiphysis:

- facilitated plate positioning over the prosthesis
- provisional plate stabilization



• Hooks in epiphyseal part:


- additional stabilization of the greater trochanter.

• Plate design:

- design aligned with anatomy of a bone
- narrowed epiphyseal part
- gentle upper edges

• Transverse holes in the epiphyseal and shaft parts:


- for cables to attach the plate to the bone.



3. SURGICAL TECHNIQUE - 7.0ChLP FEMORAL PERIPROSTHETIC PLATE 3.7220.XXX


3.1. PATIENT'S POSITIONING

Place a patient supine. Ensure lateral and AP visualization.



3.2. SURGICAL APPROACH

Lateral access. Perform a more or less extensive incision of the skin (*depending on the implant used*). The incision shall start from the top of the greater trochanter to the lateral condyle of the femur.



3.3. FRACTURE REDUCTION


Perform fracture reduction. If need be, temporarily stabilize the bone fragments with Kirschner wires and/or reduction pliers.

3.4. IMPLANT SELECTION

Select the right size of an implant to the type of fracture, bone size and structure.

3.5. PLATE INSERTION

Position the implant correctly on the bone.




3.6. TEMPORARY PLATE STABILIZATION

Stabilize the position of the implant inserting Kirschner wires into appropriate holes (acc. to procedure 4a).



The Kirschner wires in the epiphyseal part indicate (in the transverse view) the plane of the screws to be inserted. They facilitate positioning of the plate over the prosthesis.



3.7. SCREWS INSERTION

Insert screws of a suitable length into the holes of the plate.

- Cortical self-tapping screw 4.5 **[3.1471]** (acc. to procedure 4b).
- 7.0ChLP self-tapping screw 5.0 **[3.5210]** (acc. to procedure 4c).



Insert the cortical screws 4.5 into a bone fragment before inserting the locking screws.



The doctor decides about the order and number of locking and cortical screws to be inserted.

3.8. WOUND CLOSURE


Before closing the wound, take an X-Ray image in at least two projections to confirm implant position and fracture reduction. Make sure all the screws are properly tightened and do not penetrate the joint surface. Use appropriate surgical technique to close the wound.

4. SURGICAL PROCEDURES

4a. PROCEDURE OF TEMPORARY IMPLANT STABILIZATION



Stabilization using Kirschner wires

- Stabilize temporary the implant inserting Kirschner wires 2.0/210 [40.4815.210] into dedicated holes in the plate.





Stabilization in locking holes using Kirschner wires

- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.
- Insert Kirschner wire [40.4815.210] through the guide sleeve 7.0/4.0 [40.5705.740].



Stabilization using setting-compressing screw

- Insert guide sleeve 7.0/4.0 [40.5705.740] into the locking hole of the plate.
- Insert setting-compressing screw 4.0/180 [40.5706.740] through the guide sleeve 7.0/4.0 [40.5705.740].
- Tighten the nut of the setting-compressing screw 4.0/180 [40.5706.740] and push the plate to the bone.



4b. PROCEDURE OF CORTICAL SELF-TAPPING SCREW 4.5 [3.1471] INSERTION

Compression guide positioning

Position the compression guide 3.2 [40.4802.732] in a desired position:



NEUTRAL POSITION: Push the guide to the plate. It will position itself so as neutral insertion of the screw is allowed.

COMPRESSION POSITION: Do not push the guide and move it to the edge of the compression hole. The hole drilled in this position allows compressive insertion of the screw.

ANGULAR POSITION: Angular position of the guide may also be applied.

Hole drilling

Perform a hole through both cortices for a cortical screw 4.5 insertion. For drilling, use drill with scale 3.2/210 [40.5650.212] and compression guide in a desired position.



Measurement of hole depth

Insert depth measure [40.4639.550] into drilled hole until the hook of the measure rests against the outer surface of the second cortex.



Screw insertion

Insert cortical screw using torque limiting ratchet handle T 4Nm [40.6660.000] and screwdriver tip T25 [40.5684.200].



NEUTRAL POSITION

COMPRESSION POSITION

ANGULAR POSITION

4c. PROCEDURE OF 7.0ChLP SELF-TAPPING SCREW 5.0 [3.5210] INSERTION

Guide sleeve insertion

- Insert guide sleeve 7.0/4.0 **[40.5705.740]** into a locking hole of the plate.



40.5705.740

Hole drilling

Drill using drill with scale 4.0/210 **[40.5651.212]** until desired depth is reached.



40.5651.212

Measurement of hole depth

OPTION I: Read the length of the screw from the drill measure **[40.5651.212]**



40.5651.212

OPTION II: or use screw length measure **[40.5675.500]**.



40.5675.500


OPTION III: Having removed the guide sleeve 7.0/4.0 **[40.5705.740]**, use depth measure **[40.4639.550]** to determine the length of a screw.



40.4639.550

Screw insertion




Remove the guide sleeve 7.0/4.0 **[40.5705.740]**. Use torque limiting ratchet handle T 4Nm **[40.6660.000]** and screwdriver tip T25 **[40.5684.200]** to insert the locking screw.




40.6660.000




40.5684.200




OPTION I:



OPTION II:



OPTION III:



5. POSTOPERATIVE PROCEDURE

Introduce appropriate postoperative treatment. The physician decides on the post-operative treatment and its conduct. In order to avoid patient's movement limitations, introduce exercises as soon after surgery as possible. However, make sure that the limb is not fully loaded before fragments osteosynthesis is complete.

6. IMPLANT REMOVAL

















The physician decides about implant removal. In order to remove the implants from the body, unlock all the locking screws first and then remove them from the bone. This will prevent any rotation of the plate when removing the last locking screw.

7. CATALOGUE PAGES

7a. INSTRUMENT SET










Instrument set for 7.0ChLP 4x4 1/2H

15.0207.201

No.		Name	Catalogue No.	Pcs
1		Tray for 7.0ChLP instrument set 4x4 1/2H	14.0207.201	1
2		Kirschner wire 2.0/210	40.4815.210	4
3		Drill with scale 3.2/210	40.5650.212	2
4		Drill with scale 4.0/210	40.5651.212	2
5		Cannulated drill with scale 5.0/2.2/210	40.5652.212	1
6		Screwdriver tip T25-1/4	40.5684.200	1
7		Cannulated screwdriver tip T30-1/4	40.5685.200	1
7		Compression guide 3.2	40.4802.732	1
8		Guide sleeve 9.0/5.0	40.5654.750	1
9		Guide sleeve 9.0/3.2	40.5707.732	1
10		Guide sleeve 5.0/2.0	40.5654.120	1
11		Torque limiting ratchet handle T 4Nm	40.6660.000	1
12		Protective guide 9/7	40.5708.000	2
13		Guide sleeve 7.0/3.2	40.5705.732	2
14		Guide sleeve 7.0/4.0	40.5705.740	4
15		Depth measure	40.4639.550	1


Instrument set for 7.0ChLP 4x4 1/2H

15.0207.202

No.		Name	Catalogue No.	Pcs
1		Tray for 7.0ChLP instrument set 4x4 1/2H	14.0207.202	1
2		Plates bender 7.0	40.4643.700	2
3		Setting-compressing screw 4.0/180	40.5706.740	2
4		Tripod screwdriver tip 7.0ChLP	40.6271.700	1
5		Screwdriver tip T25 with holder	40.6259.200	1
6		Tap 7.0ChLP - 5.0	40.5646.200	1
7		Cortical tap HA 4.5	40.5647.200	1
8		Screw length measure	40.5675.500	1
Optional instrument				
9		Torque connector 4Nm	40.5927.040	1

Ti


7.0ChLP femoral periprosthetic plate



L R
Ster
Non
Ster

		Ti		SHORT	
①	Len	L	R		
6	222	3.7221.606	3.7220.606		
8	274	3.7221.608	3.7220.608		
10	326	3.7221.610	3.7220.610		
12	378	3.7221.612	3.7220.612		


		LONG	
②			
6	229	3.7223.606	3.7222.606
8	281	3.7223.608	3.7222.608
10	333	3.7223.610	3.7222.610
12	384	3.7223.612	3.7222.612



	Ti					
3.5210.xxx	✓	✓	✓	✓		5.0
3.1471.xxx	✓	✓		✓		4.5

Ti


7.0ChLP trochanteric periprosthetic plate



①			Ti	SHORT
①	Len			
2	42	3.7224.600		
②				LONG
2	50	3.7225.600		




		Ti							
	3.5210.xxx	✓	✓	✓	✓				5.0
	3.1471.xxx	✓	✓		✓				4.5





7.0ChLP self-tapping screw 5.0



Len	Ti
16	3.5210.016
18	3.5210.018
20	3.5210.020
22	3.5210.022
24	3.5210.024
26	3.5210.026
28	3.5210.028
30	3.5210.030
32	3.5210.032
34	3.5210.034
36	3.5210.036
38	3.5210.038
40	3.5210.040
42	3.5210.042
44	3.5210.044
46	3.5210.046
48	3.5210.048
50	3.5210.050
52	3.5210.052
54	3.5210.054
56	3.5210.056
58	3.5210.058
60	3.5210.060
65	3.5210.065
70	3.5210.070
75	3.5210.075
80	3.5210.080
85	3.5210.085
90	3.5210.090
95	3.5210.095
100	3.5210.100
105	3.5210.105
110	3.5210.110

Cortical self-tapping screw 4.5



Len	Ti
16	3.1471.016
18	3.1471.018
20	3.1471.020
22	3.1471.022
24	3.1471.024
26	3.1471.026
28	3.1471.028
30	3.1471.030
32	3.1471.032
34	3.1471.034
36	3.1471.036
38	3.1471.038
40	3.1471.040
42	3.1471.042
44	3.1471.044
46	3.1471.046
48	3.1471.048
50	3.1471.050
52	3.1471.052
54	3.1471.054
56	3.1471.056
58	3.1471.058
60	3.1471.060
65	3.1471.065
70	3.1471.070
75	3.1471.075
80	3.1471.080
85	3.1471.085
90	3.1471.090
95	3.1471.095
100	3.1471.100
105	3.1471.105
110	3.1471.110

ChM sp. z o.o.

Lewickie 3b
16-061 Juchnowiec Kościelny
Polska
tel. +48 85 86 86 100
fax +48 85 86 86 101
chm@chm.eu
www.chm.eu



CE 0197