

# Total Hip System

## b-ONE Bi-polar Heads

Adding Variety to Surgical Options

### Reduces Wear

b-ONE Bi-polar heads reduces wear through proper distribution of the load stress.

### Easy to Assemble

Femoral heads can be easily snapped into the Bi-polar insert during the surgery.

### Accurate Fitting

b-ONE Bi-polar heads offer multiple size option for wide variety of patients.



In combination with the b-ONE femoral components and the surgeon's preferred Surgical Technique, the b-ONE Bi-polar can be a cost effective and efficient solution not only to patients with femoral neck fractures and avascular necrosis, but also to surgeons looking for good patient outcome.

**b-ONE**  
O R T H O

7715 NW 46th Street Doral, FL 33166, USA

+1 (305) 675-8012

sales@novutechusa.com

**NOVUTECH**  
www.novutechusa.com

## 28.0mm Bi-polar System

Code	Product Description
3815001011	Bipolar Liner, ID 28mm, Size 1
3815002011	Bipolar Liner, ID 28mm, Size 2
3815003011	Bipolar Liner, ID 28mm, Size 3
3815004011	Bipolar Liner, ID 28mm, Size 4

Code	Product Description
3815390001	Bipolar Shell, OD 39mm, ID Size 1
3815400001	Bipolar Shell, OD 40mm, ID Size 1
3815410001	Bipolar Shell, OD 41mm, ID Size 1
3815420001	Bipolar Shell, OD 42mm, ID Size 1
3815430001	Bipolar Shell, OD 43mm, ID Size 2
3815440001	Bipolar Shell, OD 44mm, ID Size 2
3815450001	Bipolar Shell, OD 45mm, ID Size 2
3815460001	Bipolar Shell, OD 46mm, ID Size 3
3815480001	Bipolar Shell, OD 48mm, ID Size 3
3815500001	Bipolar Shell, OD 50mm, ID Size 3
3815520001	Bipolar Shell, OD 52mm, ID Size 4
3815540001	Bipolar Shell, OD 54mm, ID Size 4
3815560001	Bipolar Shell, OD 56mm, ID Size 4


Code	Product Description
7811128035	CoCr Head, 28mm-3.5
7811128000	CoCr Head, 28mm+0
7811128350	CoCr Head, 28mm+3.5
781112870S	CoCr Head, 28mm+7

# Total Hip System

## b-ONE Bi-polar Heads

Product Catalogue

 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

  
www.novutechusa.com

# Total Hip System

## OFIT™ Hip System

Cost-Effective Solution for All



### **132° Neck Angle**

Catered for Asian Anatomy

### **Modified Neck Geometry**

Improves Neck Strength and Range of Motion

### **Proportional Neck Length**

Logical Approach to Restoring Leg Length

### **Small Distal Tip**

Easy to Insert, Reduces Thigh Pain

### **Constant Insertion Path**

Easy to Insert and Easy to Assess Rotational Stability

### **Proximal Cross Section**

Provides Rotational Stability and Self Locking


### **Distal Vertical Grooves**

Improves Contact Area and Rotational Stability

### **Sand Blasting Technology**

Improves Bony On-Growth

 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

**NOVUTECH**  
www.novutechusa.com

# Total Hip System

## OFIT® Femoral Stem

Product Catalogue

### Code


7810110085  
7810110009  
7810110010  
7810110011  
7810110012  
7810110013  
7810110014  
7810110015  
7810110016

### Product Description

Femoral Stem, Proportional, 132°, Size 8.5  
Femoral Stem, Proportional, 132°, Size 9  
Femoral Stem, Proportional, 132°, Size 10  
Femoral Stem, Proportional, 132°, Size 11  
Femoral Stem, Proportional, 132°, Size 12  
Femoral Stem, Proportional, 132°, Size 13  
Femoral Stem, Proportional, 132°, Size 14  
Femoral Stem, Proportional, 132°, Size 15  
Femoral Stem, Proportional, 132°, Size 16



 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

NOVUTECH  
[www.novutechusa.com](http://www.novutechusa.com)

# Total Hip System

## FUGU™ Acetabular System

Cost-Effective Solution for All

### Firm Fixation

Macro Structures on the FUGU™ acetabular cup are oriented in two different directions for secure primary and secondary fixations.

### Easy to Insert

Hemispherical design of the FUGU™ acetabular cup's main surface is design line to line with the reamer. Which facilitates the insertion of the cup.

### Patented Locking

Patented locking mechanism of the FUGU™ cup is constitute of 12 asymmetric locking ear and tabs that provides secondary locking.

FUGUTM aims to improve clinically proven fixation technology, and incorporates patented secure locking mechanism to provide higher value for patients and surgeons world wide.



**b-one**  
ORTH O

7715 NW 46th Street Doral, FL 33166, USA

X: +1 (305) 675-8012

sales@novutechusa.com

**NOVUTECH**  
[www.novutechusa.com](http://www.novutechusa.com)

# Total Hip System

## FUGU® Acetabular System

Product Catalogue



### Code

781331441C  
781331461D  
781331481E  
781331501F  
781331521G  
781331541H  
781331561I  
781331581J  
781331601K  
781331621L  
781331641L

### Product Description

Acetabular Shell, Code C, 44mm  
Acetabular Shell, Code D, 46mm  
Acetabular Shell, Code E, 48mm  
Acetabular Shell, Code F, 50mm  
Acetabular Shell, Code G, 52mm  
Acetabular Shell, Code H, 54mm  
Acetabular Shell, Code I, 56mm  
Acetabular Shell, Code J, 58mm  
Acetabular Shell, Code K, 60mm  
Acetabular Shell, Code L, 62mm  
Acetabular Shell, Code L, 64mm

### Code

781231281C  
781231281D  
781231281E  
781231281F  
781231281G  
781231281H  
781231281I  
781231281J  
781231281K  
781231281L

### Product Description

UHMWPE Liner, Hooded, Code C, 28mm  
UHMWPE Liner, Hooded, Code D, 28mm  
UHMWPE Liner, Hooded, Code E, 28mm  
UHMWPE Liner, Hooded, Code F, 28mm  
UHMWPE Liner, Hooded, Code G, 28mm  
UHMWPE Liner, Hooded, Code H, 28mm  
UHMWPE Liner, Hooded, Code I, 28mm  
UHMWPE Liner, Hooded, Code J, 28mm  
UHMWPE Liner, Hooded, Code K, 28mm  
UHMWPE Liner, Hooded, Code L, 28mm


### Code

7814106520  
7814106525  
7814106530  
7814106535  
7814106540

### Product Description

Shell Screw, 6.5×20mm  
Shell Screw, 6.5×25mm  
Shell Screw, 6.5×30mm  
Shell Screw, 6.5×35mm  
Shell Screw, 6.5×40mm

 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

 **NOVUTECH**  
www.novutechusa.com