

Total Hip System

b-ONE Bi-polar Heads

Adding Variety to Surgical Options

Reduces Wear

b-ONE Bi-polar heads reduces wear through proper distribution of the load stress.

Easy to Assemble

Femoral heads can be easily snapped into the Bi-polar insert during the surgery.

Accurate Fitting

b-ONE Bi-polar heads offer multiple size option for wide variety of patients.



In combination with the b-ONE femoral components and the surgeon's preferred Surgical Technique, the b-ONE Bi-polar can be a cost effective and efficient solution not only to patients with femoral neck fractures and avascular necrosis, but also to surgeons looking for good patient outcome.

b-ONE
O R T H O

7715 NW 46th Street Doral, FL 33166, USA

+1 (305) 675-8012

sales@novutechusa.com

NOVUTECH
www.novutechusa.com

28.0mm Bi-polar System

Code	Product Description
3815001011	Bipolar Liner, ID 28mm, Size 1
3815002011	Bipolar Liner, ID 28mm, Size 2
3815003011	Bipolar Liner, ID 28mm, Size 3
3815004011	Bipolar Liner, ID 28mm, Size 4

Code	Product Description
3815390001	Bipolar Shell, OD 39mm, ID Size 1
3815400001	Bipolar Shell, OD 40mm, ID Size 1
3815410001	Bipolar Shell, OD 41mm, ID Size 1
3815420001	Bipolar Shell, OD 42mm, ID Size 1
3815430001	Bipolar Shell, OD 43mm, ID Size 2
3815440001	Bipolar Shell, OD 44mm, ID Size 2
3815450001	Bipolar Shell, OD 45mm, ID Size 2
3815460001	Bipolar Shell, OD 46mm, ID Size 3
3815480001	Bipolar Shell, OD 48mm, ID Size 3
3815500001	Bipolar Shell, OD 50mm, ID Size 3
3815520001	Bipolar Shell, OD 52mm, ID Size 4
3815540001	Bipolar Shell, OD 54mm, ID Size 4
3815560001	Bipolar Shell, OD 56mm, ID Size 4


Code	Product Description
7811128035	CoCr Head, 28mm-3.5
7811128000	CoCr Head, 28mm+0
7811128350	CoCr Head, 28mm+3.5
781112870S	CoCr Head, 28mm+7

Total Hip System

b-ONE Bi-polar Heads

Product Catalogue

 7715 NW 46th Street Doral, FL 33166, USA

 +1 (305) 675-8012

 sales@novutechusa.com


www.novutechusa.com

Total Hip System

OFIT™ Hip System

Cost-Effective Solution for All



132° Neck Angle

Catered for Asian Anatomy

Modified Neck Geometry

Improves Neck Strength and Range of Motion

Proportional Neck Length

Logical Approach to Restoring Leg Length

Small Distal Tip

Easy to Insert, Reduces Thigh Pain

Constant Insertion Path

Easy to Insert and Easy to Assess Rotational Stability

Proximal Cross Section

Provides Rotational Stability and Self Locking


Distal Vertical Grooves


Improves Contact Area and Rotational Stability

Sand Blasting Technology

Improves Bony On-Growth

 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

NOVUTECH
www.novutechusa.com

Total Hip System

OFIT® Femoral Stem

Product Catalogue

Code


7810110085
7810110009
7810110010
7810110011
7810110012
7810110013
7810110014
7810110015
7810110016

Product Description

Femoral Stem, Proportional, 132°, Size 8.5
Femoral Stem, Proportional, 132°, Size 9
Femoral Stem, Proportional, 132°, Size 10
Femoral Stem, Proportional, 132°, Size 11
Femoral Stem, Proportional, 132°, Size 12
Femoral Stem, Proportional, 132°, Size 13
Femoral Stem, Proportional, 132°, Size 14
Femoral Stem, Proportional, 132°, Size 15
Femoral Stem, Proportional, 132°, Size 16



 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

 **NOVUTECH**
www.novutechusa.com

Total Hip System

FUGU™ Acetabular System

Cost-Effective Solution for All

Firm Fixation

Macro Structures on the FUGU™ acetabular cup are oriented in two different directions for secure primary and secondary fixations.

Easy to Insert

Hemispherical design of the FUGU™ acetabular cup's main surface is design line to line with the reamer. Which facilitates the insertion of the cup.

Patented Locking

Patented locking mechanism of the FUGU™ cup is constitute of 12 asymmetric locking ear and tabs that provides secondary locking.

FUGUTM aims to improve clinically proven fixation technology, and incorporates patented secure locking mechanism to provide higher value for patients and surgeons world wide.



b-one
O R T H O

7715 NW 46th Street Doral, FL 33166, USA

+1 (305) 675-8012

sales@novutechusa.com

NOVUTECH
www.novutechusa.com

Total Hip System

FUGU® Acetabular System

Product Catalogue



Code

781331441C
781331461D
781331481E
781331501F
781331521G
781331541H
781331561I
781331581J
781331601K
781331621L
781331641L

Product Description

Acetabular Shell, Code C, 44mm
Acetabular Shell, Code D, 46mm
Acetabular Shell, Code E, 48mm
Acetabular Shell, Code F, 50mm
Acetabular Shell, Code G, 52mm
Acetabular Shell, Code H, 54mm
Acetabular Shell, Code I, 56mm
Acetabular Shell, Code J, 58mm
Acetabular Shell, Code K, 60mm
Acetabular Shell, Code L, 62mm
Acetabular Shell, Code L, 64mm

Code

781231281C
781231281D
781231281E
781231281F
781231281G
781231281H
781231281I
781231281J
781231281K
781231281L

Product Description

UHMWPE Liner, Hooded, Code C, 28mm
UHMWPE Liner, Hooded, Code D, 28mm
UHMWPE Liner, Hooded, Code E, 28mm
UHMWPE Liner, Hooded, Code F, 28mm
UHMWPE Liner, Hooded, Code G, 28mm
UHMWPE Liner, Hooded, Code H, 28mm
UHMWPE Liner, Hooded, Code I, 28mm
UHMWPE Liner, Hooded, Code J, 28mm
UHMWPE Liner, Hooded, Code K, 28mm
UHMWPE Liner, Hooded, Code L, 28mm


Code

7814106520
7814106525
7814106530
7814106535
7814106540

Product Description

Shell Screw, 6.5×20mm
Shell Screw, 6.5×25mm
Shell Screw, 6.5×30mm
Shell Screw, 6.5×35mm
Shell Screw, 6.5×40mm

 7715 NW 46th Street Doral, FL 33166, USA

 X: +1 (305) 675-8012

 sales@novutechusa.com

 novUTECH
www.novutechusa.com