



4,0 ChM Locked Plating
ChLPsystem

5,0 ChM Locked Plating
ChLPsystem

7,0 ChM Locked Plating
ChLPsystem

WINKELSTABILE PLATTEN
LOCKING PLATES
PLACAS DE BLOQUEO



www.chm.eu

ZEICHENERKLÄRUNG	7	4,0ChLP TALUSPLATTE	30
SYMBOL DESCRIPTION		4.0ChLP TALUS PLATE	
DESCRIPCIÓN DE SÍMBOLOS		4,0ChLP PLACA PARA ASTRÁGALO	
WIE KATALOG LESEN	8	4,0ChLP RECHTECKIGE PLATTE	31
HOW TO READ THIS CATALOGUE		4.0ChLP RECTANGULAR PLATE	
CÓMO LEER ESTE CATÁLOGO		4,0ChLP PLACA RECTAGULAR	
4,0ChLP VOLARE RADIUSPLATTE	9	4,0ChLP 3D-PLATTE	32
4.0ChLP RADIUS VOLAR PLATE		4.0ChLP 3D PLATE	
4,0ChLP PLACA VOLAR DE RADIO		4,0ChLP PLACA 3D	
4,0ChLP ERWEITERTE VOLARE RADIUSPLATTE	10	4,0ChLP GEBOGENE PLATTE	33
4.0ChLP EXPANDED RADIUS VOLAR PLATE		4.0ChLP CURVED PLATE	
4,0ChLP PLACA RADIAL DE MANO EXTENDIDA		4,0ChLP PLACA CURVADA	
4,0ChLP DORSALE RADIUSPLATTE	11	4,0ChLP KUBOIDPLATTE	34
4.0ChLP DORSAL-RADIAL PLATE		4.0ChLP PLATE FOR CUBOID BONE	
4,0ChLP PLACA RADIAL DORSAL		4,0ChLP PLACA PARA HUESO CUBOIDES	
4,0ChLP DORSALE RADIUSPLATTE	12	4,0ChLP NAVIKULARPLATTE	35
4.0ChLP DORSAL-RADIAL PLATE		4.0ChLP PLATE FOR NAVICULAR BONE	
4,0ChLP PLACA RADIAL DORSAL		4,0ChLP PLACA PARA HUESO NAVICULAR	
PLATTEN FÜR DISTAL RADIUS – DORSAL	13	4,0ChLP PLATTE TYP-6	36
PLATES FOR DISTAL RADIUS – DORSAL SIDE		4.0ChLP 6-TYPE PLATE	
PLACAS DEL RADIO DISTAL – LADO DORSAL		4,0ChLP PLACA TIPO-6	
4,0ChLP GERADE PLATTE	14	4,0ChLP SCHRÄGE T-PLATTE	37
4.0ChLP STRAIGHT PLATE		4.0ChLP T PLATE OBLIQUE	
4,0ChLP PLACA RECTA		4,0ChLP PLACA T OBLICUA	
4,0ChLP L-PLATTE	15	4,0ChLP L-PLATTE	38
4.0ChLP L PLATE		4.0ChLP L PLATE	
4,0ChLP PLACA L		4,0ChLP PLACA L	
4,0ChLP SCHRÄGE L-PLATTE	16	4,0ChLP GERADE BREITE PLATTE	39
4.0ChLP OBLIQUE L PLATE		4.0ChLP STRAIGHT WIDE PLATE	
4,0ChLP PLACA L OBLICUA		4,0ChLP PLACA RECTA ANCHA	
4,0ChLP SCHRÄGE L-PLATTE	17	4,0ChLP MITTELFUß- UND ZEHEKNOCHEPLATTE	40
4.0ChLP OBLIQUE L PLATE		4.0ChLP LISFRANC PLATE	
4,0ChLP PLACA L OBLICUA		4,0ChLP PLACA METATARSOFALENGICA	
4,0ChLP T-PLATTE	18	4,0ChLP REKONSTRUKTIONSPLATTE	41
4.0ChLP T-PLATE		4.0ChLP RECONSTRUCTION PLATE	
4,0ChLP PLACA T		4,0ChLP PLACA DE RECONSTRUCCIÓN	
4,0ChLP L-PLATTE	19	INSTRUMENTARIUM FÜR 4,0ChLP - FUSS	42
4.0ChLP L PLATE		SET FOR 4.0ChLP - FOOT	
4,0ChLP PLACA L		SET PARA 4,0ChLP - PIE	
4,0ChLP PLATTE FÜR DISTALE ULNA	20	INSTRUMENTE 4,0ChLP	43
4.0ChLP DISTAL ULNAR PLATE		INSTRUMENT SET FOR 4.0ChLP	
4,0ChLP PLACA PARA EL CÚBITO DISTAL		SET DE INSTRUMENTAL PARA 4,0ChLP	
INSTRUMENTARIUM FÜR 4,0ChLP - DISTALER RADIUS	21	4,0ChLP 2-LOCH-FUSSPLATTE	44
SET 4.0ChLP - DISTAL RADIUS		4.0ChLP 2-HOLES FOOT PLATE	
SET PARA 4,0ChLP - RADIO DISTAL		4,0ChLP PLACA PARA PIE DE 2 ORIFICIOS	
PALETTE FÜR 4,0 ChLP IMPLANTATE - DISTALER RADIUS	22	4,0ChLP SCHRÄGE T-PLATTE FÜR KLEINKNOCHEN	45
PALETTE FOR 4.0ChLP IMPLANTS - DISTAL RADIUS		4.0ChLP OBLIQUE T PLATE FOR SMALL BONES	
PALETA PARA IMPLANTES 4,0ChLP - RADIO DISTAL		4,0ChLP PLACA T OBLICUA PARA HUESOS PEQUEÑOS	
INSTRUMENTE 4,0ChLP	23	4,0ChLP X-PLATTE	46
INSTRUMENT SET FOR 4.0ChLP		4.0ChLP X PLATE	
SET DE INSTRUMENTAL PARA 4,0ChLP		4,0ChLP PLACA X	
4,0ChLP HANDGELENKPLATTE FÜR ARTHRODESE	24	4,0ChLP T-PLATTE	47
4.0ChLP WRIST ARTHRODESIS PLATE		4.0ChLP T PLATE	
4,0ChLP PLACA PARA ARTRODESIS DE MUÑECA		4,0ChLP PLACA T	
BEHÄLTER FÜR IMPLANTATE 4,0ChLP 3.7206 4X2 1/2H	25	INSTRUMENTARIUM FÜR 4,0ChLP - HALLUX	48
STAND FOR 4.0ChLP IMPLANTS 3.7206 4X2 1/2H		SET FOR 4.0ChLP - HALLUX	
CONTENEDOR PARA IMPLANTES 4,0ChLP 3.7206 4X2 1/2H		SET PARA 4,0ChLP - HALLUX	
INSTRUMENTE 4,0ChLP	26	4,0ChLP T-KEILPLATTE	49
INSTRUMENT SET FOR 4.0ChLP		4.0ChLP WEDGE T PLATE	
SET DE INSTRUMENTAL PARA 4,0ChLP		4,0ChLP PLACA EN T TRAPEZOIDAL	
4,0ChLP PLATTE FÜR RADIUSKOPF	27	4,0ChLP X-KEILPLATTE	50
4.0ChLP RADIAL HEAD PLATE		4.0ChLP WEDGE X PLATE	
4,0ChLP PLACA DE LA CABEZA DE RADIO		4,0ChLP PLACA EN X TRAPEZOIDAL	
4,0ChLP PLATTE FÜR RADIUSHALS	28	INSTRUMENTARIUM FÜR 4,0ChLP - KEILOSTEOTOMIE	51
4.0ChLP RADIAL HEAD NECK PLATE		SET 4.0ChLP – WEDGE OSTEOTOMY	
4,0ChLP PLACA DEL CUELLO DE LA CABEZA DE RADIO		EL SET DE 4,0ChLP - OSTEOTOMÍA EN CUÑA	
INSTRUMENTARIUM 4,0ChLP - RADIUSKOPF	29	4,0ChLP INTRAOSSEÄRE PLATTE	52
SET FOR 4.0ChLP - RADIAL HEAD BONE		4.0ChLP ENDOSTEAL PLATE	
SET PARA 4,0ChLP - CABEZA DE RADIO		4,0ChLP PLACA INTRAÓSEA	

INSTRUMENTARIUM FÜR 4,0ChLP 3.7060/3.7061 4x4 H	53	5,0ChLP DISTALE DORSOLATERALE HUMERUSPLATTE	78
SET FOR 4.0ChLP 3.7060/3.7061 4x4 H		5.0ChLP DISTAL DORSOLATERAL HUMERUS PLATE	
SET PARA 4,0ChLP 3.7060/3.7061 4x4 H		5,0ChLP PLACA DISTAL DORSOLATERAL DEL HÚMERO	
INSTRUMENTARIUM FÜR 4,0ChLP 3.7060/3.7061 4x4 1/2H	54	5,0ChLP GERADE PLATTE FÜR UNTERARM	79
INSTRUMENT SET FOR 4.0ChLP 3.7060/3.7061 4x4 1/2H		5,0ChLP FOREARM PLATE STRAIGHT	
SET DE INSTRUMENTAL PARA 4,0ChLP 3.7060/3.7061 4x4 1/2H		5,0ChLP PLACA DE ANTEBRAZO RECTA	
INSTRUMENTARIUM FÜR 4,0ChLP 3.7060/3.7061 4x2 1/2H	55	5,0ChLP OLECRANONPLATTE	80
INSTRUMENT SET FOR 4.0ChLP 3.7060/3.7061 4x2 1/2H		5.0ChLP OLECRANON PLATE	
SET DE INSTRUMENTAL PARA 4,0ChLP 3.7060/3.7061 4x2 1/2H		5,0ChLP PLACA DE OLÉCRANON	
VERRIEGELUNGSSCHRAUBEN	56	5,0 ChLP HANDGELENKPLATTEN FÜR ARTHRODESE	81
LOCKING ELEMENTS		5.0ChLP WRIST ARTHRODESIS PLATES	
ELEMENTOS DE BLOQUEO		5,0ChLP PLACA PARA ARTRODESIS DE MUÑECA	
4,5ChLP SCHMALE GERADE REKONSTRUKTIONSPLATTE	59	5,0ChLP VOLARE RADIUSPLATTE	82
4.5ChLP NARROW STRAIGHT RECONSTRUCTION PLATE		5.0ChLP RADIUS VOLAR PLATE	
4,5ChLP PLACA DE RECONSTRUCCIÓN RECTA ESTRECHA		5,0ChLP PLACA RADIAL DE MANO	
4,5ChLP GERADE REKONSTRUKTIONSPLATTE	60	5,0ChLP PROXIMALE TIBIAPLATTE T	83
4.5ChLP STRAIGHT RECONSTRUCTION PLATE		5.0ChLP PROXIMAL TIBIAL T PLATE	
4,5ChLP PLACA DE RECONSTRUCCIÓN RECTA		5,0ChLP PLACA TIBIAL PROXIMAL T	
5,0 ChLP SCHMALE KOMPRESSIONSPLATTE LOW CONTACT	61	5,0ChLP PROXIMALE LATERALE TIBIAPLATTE	84
5.0ChLP NARROW COMPRESSION LIMITED CONTACT PLATE		5.0ChLP PROXIMAL LATERAL TIBIA PLATE	
5,0ChLP PLACA ESTRECHA DE COMPRESIÓN DE CONTACTO LIMITADO		5,0ChLP PLACA DE LA TIBIA PROXIMAL LATERAL	
5,0ChLP GERADE DRITTELROHRPLATTE 1/3	62	5,0ChLP SCHMALE TIBIAPLATTE	85
5.0ChLP TUBULAR STRAIGHT PLATE 1/3		5.0ChLP NARROW TIBIAL PLATE	
5,0ChLP PLACA RECTA TUBULAR 1/3		5,0ChLP PLACA TIBIAL ESTRECHA	
5,0ChLP GERADE DRITTELROHRPLATTE 1/3	63	5,0 ChLP SCHMALE POSTERIORE TIBIAPLATTE	86
5.0ChLP TUBULAR STRAIGHT PLATES 1/3		5.0ChLP POSTERIOR NARROW TIBIA PLATE	
5,0ChLP PLACA RECTA TUBULAR 1/3		5,0ChLP PLACA POSTERIOR DE LA TIBIA ESTRECHA	
5,0ChLP AC-GELENK-PLATTE	64	5,0 ChLP BREITE POSTERIORE TIBIAPLATTE	87
5.0ChLP ACJ PLATE		5.0ChLP POSTERIOR WIDE TIBIA PLATE	
5,0ChLP placa ACJ		5,0ChLP PLACA TIBIAL POSTERIOR ANCHA	
INSTRUMENTARIUM 5,0ChLP 3.7242 4x2 1/2H	65	5,0ChLP DISTALE MEDIALE TIBIAPLATTE	88
INSTRUMENT SET FOR 5.0ChLP 3.7242 4x2 1/2H		5.0ChLP TIBIAL DISTAL MEDIAL PLATE	
INSTRUMENTAL PARA 5,0ChLP 3.7242 4x2 1/2H		5,0ChLP PLACA TIBIAL DISTAL MEDIAL	
4,5ChLP KLAVIKULA-HAKENPLATTE	66	5,0ChLP DISTALE MEDIALE TIBIAPLATTE	89
4.5ChLP CLAVICULAR HOOK PLATE		5.0ChLP DISTAL MEDIAL TIBIAL PLATE	
4,5ChLP PLACA CLAVÍCULA CON GANCHO		5,0ChLP PLACA TIBIAL DISTAL MEDIAL	
5,0ChLP KLAVIKULA S-PLATTE	67	5,0ChLP DISTALE MEDIALE PLATTE	90
5.0ChLP CLAVICLE S PLATE		5.0ChLP DISTAL MEDIAL TIBIAL PLATE	
5,0ChLP PLACA CLAVICULAR EN S		5,0ChLP PLACA PARA TIBIA DISTAL MEDIAL	
5,0ChLP S-PLATTE FÜR KLAVIKULASCHAFT	68	5,0ChLP DISTALE L-TIBIAPLATTE	91
5.0ChLP CLAVICLE SHAFT S PLATE		5.0ChLP DISTAL TIBIAL L PLATE	
5,0ChLP PLACA EN S PARA EL EJE DE LA CLAVÍCULA		5,0ChLP PLACA TIBIAL DISTAL L	
5,0ChLP MEDIALE SCHULTERBLATTPLATTE	69	5,0 ChLP DISTALE LATERALE FIBULAPLATTE	92
5.0ChLP SCAPULAR MEDIAL PLATE		5.0ChLP DISTAL LATERAL FIBULA PLATE	
5,0ChLP PLACA MEDIAL DE ESCÁPULA		5,0ChLP PLACA DISTAL LATERAL DE PERONÉ	
5,0ChLP LATERALE SCHULTERBLATTPLATTE	70	5,0ChLP KALKANEUSPLATTE	93
5.0ChLP SCAPULAR LATERAL PLATE		5.0ChLP CALCANEAL PLATE	
5,0ChLP PLACA LATERAL DE ESCÁPULA		5,0ChLP PLACA PARA CALCÁNEO	
5,0ChLP PFANNEN-SCHULTERBLATTPLATTE	71	5,0ChLP KALKANEUSPLATTE	94
5.0ChLP SCAPULAR GLENOID PLATE		5.0ChLP CALCANEAL PLATE	
5,0ChLP PLACA DE GLENOIDE		5,0ChLP PLACA PARA CALCÁNEO	
5,0ChLP HUMERUSPLATTE	72	4,5ChLP KALKANEUS-HAKENPLATTE	95
5.0ChLP HUMERUS PLATE		4.5ChLP CALCANEAL PLATE WITH HOOK	
5,0ChLP PLACA HUMERAL		4,5ChLP PLACA PARA CALCÁNEO CON GANCHO	
5,0ChLP PŁYTKA RAMIENNA DALSZA Y	73	INSTRUMENTARIUM FÜR 5,0ChLP PLATTEN	96
5.0ChLP DISTAL HUMERUS Y PLATE		INSTRUMENT SET FOR 5.0ChLP PLATES	
5,0ChLP PLACA HUMERAL DISTAL Y		SET DE INSTRUMENTAL PARA PLACAS 5,0ChLP	
5,0 ChLP DISTALE LATERALE HUMERUSPLATTE	74	VERRIEGELUNGSSCHRAUBEN	101
5.0ChLP DISTAL LATERAL HUMERUS PLATE		LOCKING ELEMENTS	
5,0ChLP PLACA DISTAL LATERAL DEL HÚMERO		ELEMENTOS DE BLOQUEO	
5,0 ChLP DISTALE POSTERIORLATERALE HUMERUSPLATTE	75	PALETEN FÜR IMPLANTATE	106
5.0ChLP DISTAL POSTERIOR LATERAL HUMERUS PLATE		PALETTES FOR IMPLANTS	
5,0ChLP PLACA DISTAL POSTEROLATERAL DEL HÚMERO		PALETAS PARA IMPLANTES	
5,0 ChLP DISTALE POSTERIORMEDIALE HUMERUSPLATTE	76	7,0 ChLP LOW-CONTACT KOMPRESSIONSPLATTE, BREIT	118
5.0ChLP DISTAL POSTERIOR MEDIAL HUMERUS PLATE		7.0ChLP WIDE COMPRESSION LIMITED CONTACT PLATE	
5,0ChLP PLACA DISTAL POSTEROMEDIAL DEL HÚMERO		7,0ChLP PLACA ANCHA DE COMPRESIÓN DE CONTACTO LIMITADO	
5,0ChLP DISTALE MEDIALE HUMERUSPLATTE	77	7,0ChLP FEMURKONDYLENPLATTE	119
5.0ChLP DISTAL MEDIAL HUMERUS PLATE		7.0ChLP CONDYLAR FEMORAL PLATE	
5,0ChLP PLACA DISTAL MEDIAL DEL HÚMERO		7,0ChLP PLACA CONDILAR PARA FÉMUR	

7,0ChLP TIBIAKONDYLENPLATTE 120	7,0ChLP DYNAMISCHE HÜFTPLATTE 143
7,0ChLP CONDYLAR TIBIAL PLATE	7,0ChLP DYNAMIC HIP PLATE
7,0ChLP PLACA CONDILAR PARA TIBIA	7,0ChLP PLACA DINÁMICA DE CADERA
7,0ChLP BREITE TIBIAPLATTE 121	7,0ChLP DYNAMISCHE HÜFTPLATTE 144
7,0ChLP WIDE TIBIAL PLATE	7,0ChLP DYNAMIC HIP PLATE
7,0ChLP PLACA TIBIAL ANCHA	7,0ChLP PLACA DINÁMICA DE CADERA
7,0 ChLP PROXIMALE LATERALE TIBIAPLATTE 122	INSTRUMENTARIUM FÜR DIE DYNAMISCHE HÜFTPLATTE 145
7,0ChLP PROXIMAL LATERAL TIBIAL PLATE	INSTRUMENT SET FOR DYNAMIC HIP PLATE
7,0ChLP PLACA TIBIAL PROXIMAL LATERAL	SET DE INSTRUMENTAL PARA PLACA DINÁMICA DE CADERA
7,0ChLP LATERALE MEDIALE TIBIAPLATTE 123	7,0ChLP DISTANZPLATTE FÜR TIBIA 147
7,0ChLP MEDIAL PROXIMAL TIBIAL PLATE	7,0ChLP DISTANCE TIBIAL PLATE
7,0ChLP PLACA TIBIAL MEDIAL PROXIMAL	7,0ChLP PLACA DISTANCIA TIBIAL
7,0ChLP ANTEROLATERALE PLATTE FÜR ARTHRODESE TT 124	7,0ChLP DISTANZPLATTE FÜR TIBIA 148
7,0ChLP ANTEROLATERAL ARTHRODESIS PLATE TT	7,0ChLP DISTANCE TIBIAL PLATE
7,0ChLP PLACA PARA ARTHRODESIS ANTEROLATERAL TT	7,0ChLP PLACA DISTANCIA TIBIAL
7,0ChLP LATERALE PLATTE FÜR ARTHRODESE TTC 125	7,0ChLP DISTANZPLATTE FÜR FEMUR 149
7,0ChLP LATERAL ARTHRODESIS PLATE TTC	7,0ChLP DISTANCE FEMUR PLATE
7,0ChLP PLACA PARA ARTHRODESIS LATERAL TTC	7,0ChLP PLACA DISTANCIA FEMORAL
7,0ChLP LATERALE PLATTE FÜR ARTHRODESE TT 126	INSTRUMENTARIUM FÜR OSTEOTOMIE 150
7,0ChLP LATERAL ARTHRODESIS PLATE TT	INSTRUMENT SET FOR OSTEOTOMY
7,0ChLP PLACA PARA ARTHRODESIS LATERAL TT	SET DE INSTRUMENTAL PARA OSTEOTOMÍA
7,0ChLP POSTERIORE PLATTE FÜR ARTHRODESE TTC 127	VERRIEGELUNGSSCHRAUBEN 152
7,0ChLP POSTERIOR ARTHRODESIS PLATE TTC	LOCKING ELEMENTS
7,0ChLP PLACA PARA ARTHRODESIS POSTERIOR TTC	ELEMENTOS DE BLOQUEO
7,0ChLP POSTERIORE PLATTE FÜR ARTHRODESE TT 128	PALETEN FÜR IMPLANTATE 157
7,0ChLP POSTERIOR ARTHRODESIS PLATE TT	PALETES FOR IMPLANTS
7,0ChLP PLACA PARA ARTHRODESIS POSTERIOR TT	PALETAS PARA IMPLANTES
7,0ChLP ANTERIORE PLATTE FÜR ARTHRODESE TT 129	PROBEIMPLANTATE FÜR 4,0 ChLP PLATTEN 160
7,0ChLP ANTERIOR ARTHRODESIS PLATE TT	4,0ChLP PLATE TRIAL
7,0ChLP PLACA PARA ARTHRODESIS ANTERIOR TT	PLANTILLAS PARA PLACAS 4,0ChLP
INSTRUMENTARIUM FÜR 7,0ChLP 130	PROBEIMPLANTATE FÜR 5,0 ChLP PLATTEN 162
INSTRUMENT SET FOR 7,0ChLP PLATES	5,0ChLP PLATE TRIAL
SET DE INSTRUMENTAL PARA PLACAS 7,0ChLP	PLANTILLAS PARA PLACAS 5,0ChLP
INSTRUMENTARIUM FÜR 7,0ChLP 131	PROBEIMPLANTATE FÜR 7,0 ChLP PLATTEN 164
INSTRUMENT SET FOR 7,0ChLP PLATES	7,0ChLP PLATE TRIAL
SET DE INSTRUMENTAL PARA PLACAS 7,0ChLP	PLANTILLAS PARA PLACAS 7,0ChLP
7,0ChLP PERIPROTHETISCHE FEMURPLATTE 132	ZUSATZINSTRUMENTE FÜR PLATTEN 165
7,0ChLP PROXIMAL FEMORAL PERIPROTHETIC PLATE	ADDITIONAL INSTRUMENTS FOR PLATES
7,0ChLP PLACA FEM. PERIPROTÉSICA	INSTRUMENTOS ADICIONALES PARA PLACAS
7,0ChLP PERIPROTHETISCHE TROCHANTERPLATTE 133	
7,0ChLP TROCHANTERIC PERIPROTHETIC PLATE	
7,0ChLP PLACA PERIPROTÉSICA TROCANTÉREA	
CERKLAGE DRAHT MIT KLEMME - GEFLEHT - 2,0X600 134	
WIRE FOR CERCLAGE WITH CLAMP-WEAVE-2.0X600	
ALAMBRE PARA CERCLAJE CON GRAPA - TEJIDO - 2.0X600	
INSTRUMENTARIUM - CERCLAGE 9X4 H 135	
INSTRUMENT SET - CERCLAGE 9X4H	
SET DE INSTRUMENTAL - CERCLAJE 9X4H	
INSTRUMENTARIUM - CERCLAGE 9X4 H 136	
INSTRUMENT SET - CERCLAGE 9X4H	
SET DE INSTRUMENTAL - CERCLAJE 9X4H	
7,0ChLP PROXIMALE FEMURPLATTE 137	
7,0ChLP PROXIMAL FEMORAL PLATE	
7,0ChLP PLACA FEMORAL PROXIMAL	
7,0ChLP PLATTE FÜR FEMURSCHAFT 138	
7,0ChLP FEMORAL SHAFT PLATE	
7,0ChLP PLACA DEL EJE FEMORAL	
7,0 ChLP SCHMALE KOMPRESSIONSPLATTE LOW CONTACT 139	
7,0ChLP NARROW COMPRESSION PLATE WITH LIMITED CONTACT	
7,0ChLP PLACA ESTRECHA DE COMPRESIÓN DE CONTACTO LIMITADO	
INSTRUMENTARIUM FÜR 7,0ChLP PLATTEN 140	
INSTRUMENT SET FOR 7,0ChLP PLATES	
SET DE INSTRUMENTAL PARA PLACAS 7,0ChLP	



Sprawdź również nasz katalog ORTOPEDIA DZIECIĘCA, gdzie znajdziesz płytki Heki i inne implanty

ChM®



Innovative Solutions for *Pediatric Orthopaedics*

	Reines Titan Pure titanium Titanio puro		Dicke [mm] Thickness [mm] Espesor [mm]
	Titan-Legierung Titanium alloy Aleación de titanio		Maßstab 1:1 Scale 1:1 Escala 1:1
	Kobalt-Legierung Cobalt alloy Aleación de cobalto		Zahl der Gewindelöcher im Plattenschaft Number of threaded holes in the shaft part of the plate Número de orificios roscados en la parte de eje de la placa
	Links Left Izquierdo		Anzahl der Gewinde- und Kompressionslöcher im Plattenschaft Number of threaded and compression holes in the shaft part of the plate Número de orificios roscados y de compresión en la parte de eje de la placa
	Rechts Right Derecho		Gesamtzahl der Gewindelöcher im Plattenschaft Number of locking holes in the plate Número de orificios bloqueados en la placa
	Erhältlich in der linken/rechten Version Available versions: left/right Disponible en versión izquierda/derecha		Variabler Winkel Variable angle Ángulo variable
	Länge Length Longitud		Kortikalis Cortical Cortical
	Torx-Schraubenkopf Torx drive Cabeza torx		Spongiosa Cancellous Esponjoso
	Kanülierter Torx-Schraubenkopf Torx drive cannulated Cabeza torx canulada		Erhältlich in der sterilen/nicht sterilen Version Available in sterile/ non- sterile condition Disponible en condición estéril/no-estéril
	Hexagonaler Schraubenkopf Hexagonal drive Cabeza hexagonal		Siehe chirurgische Anleitung See surgery technique Vea técnica quirúrgica
	Hexagonaler kanülierter Schraubenkopf Hexagonal drive cannulated Cabeza hexagonal canulada		
	Kanüliert Cannulated Canulado		
	Verriegelung Locking Bloqueado		
	Durchmesser [mm] Diameter [mm] Diámetro [mm]		
	Länge H [mm] H length [mm] Longitud H [mm]		
	Winkel Angle Ángulo		
	Erhältliche Längen Available lengths Longitudes disponibles		
	Erhältliche Zahl der Löcher Available number of holes Número disponible de orificios		

Die Implantate beschreibenden Piktogramme
Icons describing the implant
Iconos que describen el implante

**Kat.-Nr.,
Catalogue no.
Número de catálogo**

Die Implantate beschreibenden Piktogramme
Icons describing the implant
Iconos que describen el implante



7.0ChLP PROXIMALE FEMURPLATTE
7.0ChLP PROXIMAL FEMORAL PLATE
7.0ChLP PLACA FEMORAL PROXIMAL

7.0ChLP system

	2	4	6	8
132	3.7023.602	3.7022.602		
174	3.7023.604	3.7022.604		
216	3.7023.606	3.7022.606		
258	3.7023.608	3.7022.608		



7.0ChLP PLATTE FÜR FEMURSCHAFT
7.0ChLP FEMORAL SHAFT PLATE
7.0ChLP PLACA DEL EJE FEMORAL

7.0ChLP system

	12	14	16
209	3.7062.610	3.7062.612	3.7062.614
251		3.7062.614	3.7062.616
292			
333			

	3.5210.xxx	3.5232.xxx	3.5216.xxx	3.5221.xxx	3.5224.xxx	3.5228.xxx	3.1471.xxx	3.5219.008	3.1221.070
5.0	✓	✓	✓					✓	
5.4		✓							
5.0	✓								
7.3				✓					
7.3					✓				
6.5						✓			
4.5							✓		
5.0								✓	
									✓

	3.5210.xxx	3.5232.xxx	3.5216.xxx	3.1471.xxx	3.5219.008	3.1221.070
5.0	✓	✓	✓		✓	
5.4		✓				
5.0			✓			
4.5				✓		
5.0					✓	
						✓



7.0ChLP PROXIMALE FEMURPLATTE
7.0ChLP PROXIMAL FEMORAL PLATE
7.0ChLP PLACA FEMORAL PROXIMAL

7.0ChLP system



7.0ChLP PLATTE FÜR FEMURSCHAFT
7.0ChLP FEMORAL SHAFT PLATE
7.0ChLP PLACA DEL EJE FEMORAL

7.0ChLP system

Paletten für 7.0ChLP Platten 3.7022/7023 mit Inlay
Paletten für 7.0ChLP plates 3.7022/3.7023
Paletas para placas 7.0ChLP 3.7022/7023 con inlay

40.5704.900

Paletten für 7.0ChLP Platten 3.7062
Paletten für 7.0ChLP plates 3.7062
Paletas para placas 7.0ChLP 3.7062

40.5704.900

Zurücklock 1 (3.7023)
Averaj block 1 (3.7023)
Bloque guía 1 (3.7023)

40.5702.100

Zurücklock K (3.7022)
Averaj block K (3.7022)
Bloque guía K (3.7022)

40.5702.200

Platte 3.7022.604/3.7023.604
Proximale Implantat
Placa 3.7022.604/3.7023.604
Proximale para placa 3.7022.604/3.7023.604

40.7022.004

Platte 3.7062.610/3.7062.612
Femurshaft Implantat
Placa 3.7062.610/3.7062.612
Femurshaft para placa 3.7062.610/3.7062.612

40.7062.004





7,0CHLP SPONGIOSASCHRAUBE 5,4

7,0CHLP CANCELLOUS SCREW 5.4

7,0CHLP TORNILLO ESPONJOSA 5,4

VERRIEGELUNGSSCHRAUBEN
LOCKING ELEMENTS
ELEMENTOS DE BLOQUEO

Die Implantate beschreibenden
Piktogramme
Icons describing the implant
Iconos que describen el implante





EMPFOHLEN
SUGGESTED
RECOMENDADOS

35	3.5232.035	35	3.1380.035
40	3.5232.040	40	3.1380.040
45	3.5232.045	45	3.1380.045
50	3.5232.050	50	3.1380.050
55	3.5232.055	55	3.1380.055
60	3.5232.060	60	3.1380.060
65	3.5232.065	65	3.1380.065
70	3.5232.070	70	3.1380.070
75	3.5232.075	75	3.1380.075
80	3.5232.080	80	3.1380.080







EMPFOHLEN
SUGGESTED
RECOMENDADOS

50	3.5224.050	50	3.1665.050
55	3.5224.055	55	3.1665.055
60	3.5224.060	60	3.1665.060
65	3.5224.065	65	3.1665.065
70	3.5224.070	70	3.1665.070
75	3.5224.075	75	3.1665.075
80	3.5224.080	80	3.1665.080
85	3.5224.085	85	3.1665.085
90	3.5224.090	90	3.1665.090
95	3.5224.095	95	3.1665.095


Länge
Length
Longitud

Kat.-Nr.,
Catalogue no.
Número de catálogo





Zielblock R (3.7226/3.7232)	
Aiming block R (3.7226/3.7232)	40.8236.000
Blque guia R (3.7226/3.7232)	
Zielblock L (3.7227/3.7233)	
Aiming block L (3.7227/3.7233)	40.8237.000
Blque guia L (3.7227/3.7233)	
Zielblock R (3.7228/3.7234)	
Aiming block R (3.7228/3.7234)	40.8238.000
Blque guia R (3.7228/3.7234)	
Zielblock L (3.7229/3.7235)	
Aiming block L (3.7229/3.7235)	40.8239.000
Blque guia L (3.7229/3.7235)	
Zielblock R (3.7230/3.7236)	
Aiming block L (3.7230/3.7236)	40.8240.000
Blque guia R (3.7230/3.7236)	
Zielblock L (3.7231/3.7237)	
Aiming block L (3.7231/3.7237)	40.8241.000
Blque guia L (3.7231/3.7237)	

02.9



				Ti
①	4	56	3.7233.604	3.7232.604
	6	72	3.7233.606	3.7232.606
②	4	56	3.7235.604	3.7234.604
	6	72	3.7235.606	3.7234.606
③	4	56	3.7237.604	3.7236.604
	6	72	3.7237.606	3.7236.606

		TiA	Co			VA		
3.5164.xxx	✓		✓	✓	✓	✓	2.4	
3.5165.xxx	✓		✓	✓	✓	✓	2.7	
3.1220.xxx	✓		✓			✓	2.7	
4.5235.xxx		✓	✓	✓	✓	✓	2.4	




Palette für 4,0ChLP Implantate - distal Radius
Palette for 4.0ChLP implants - distal radius
Paleta implantes 4,0ChLP-radio distal


40.5713.170



Zielblock R (3.7226/3.7232)	40.8236.000	
Aiming block R (3.7226/3.7232)		
Bloque guía R (3.7226/3.7232)		
Zielblock L (3.7227/3.7233)	40.8237.000	
Aiming block L (3.7227/3.7233)		
Bloque guía L (3.7227/3.7233)		
Zielblock R (3.7228/3.7234)	40.8238.000	Platte 3.7234.606 Probeimplantat
Aiming block R (3.7228/3.7234)		Plate 3.7234.606 trial
Bloque guía R (3.7228/3.7234)		Plantilla para placa 3.7234.606
Zielblock L (3.7229/3.7235)	40.8239.000	Platte 3.7235.606 Probeimplantat
Aiming block L (3.7229/3.7235)		Plate 3.7235.606 trial
Bloque guía L (3.7229/3.7235)		Plantilla para placa 3.7235.606
Zielblock R (3.7230/3.7236)	40.8240.000	
Aiming block R (3.7230/3.7236)		
Bloque guía R (3.7230/3.7236)		
Zielblock L (3.7231/3.7237)	40.8241.000	
Aiming block L (3.7231/3.7237)		
Bloque guía L (3.7231/3.7237)		



	Len	L	R	Ti
①	4,0ChLP Dorsale Radiusplatte H schmal 4.0ChLP dorsal-radial h plate narrow 4,0ChLP placa radial dorsal en H estrecha			
	63	3.7079.000	3.7080.000	
②	4,0ChLP Dorsale Radiusplatte H breite 4.0ChLP dorsal-radial h plate wide 4,0ChLP placa radial dorsal en H ancha			
	64	3.7081.000	3.7082.000	
③	4,0ChLP Dorsale Radiusplatte Pi schmal 4.0ChLP dorsal-radial pi plate narrow 4,0ChLP placa radial dorsal en Pi estrecha			
	63	3.7075.000	3.7076.000	
④	4,0ChLP Dorsale Radiusplatte Pi breite 4.0ChLP dorsal-radial pi plate wide 4,0ChLP placa radial dorsal en Pi ancha			
	64	3.7077.000	3.7078.000	



		TiA	Co		VA		
	3.5164.xxx	✓		✓	✓	✓	2.4
	3.5165.xxx	✓		✓	✓	✓	2.7
	3.1220.xxx	✓		✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	2.4



Palette für 4,0 ChLP Implantate - distaler Radius
Palette for 4.0ChLP implants - distal radius
Paleta para implantes 4,0ChLP - radio distal






40.5713.150

Austauschbare Platten zur Verwendung mit gerader Platte in der 2-Platten-Technik
Interchangeable plates used with straight plate in 2-plates technique.
Placas intercambiables utilizadas con placas rectas en la técnica de 2 placas.

1:1



1:1



3.4047.6xx



3.4041.6xx

3.4043.6xx

3.4037.6xx

3.4045.6xx

3.4048.6xx





3.4042.6xx

3.4044.6xx

3.4038.6xx

3.4046.6xx





			
	Len		Ti
3	47	3.4047.603	
4	58	3.4047.604	

								
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4




Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)

40.5713.040



Platte 3.4047.603 Probeimplantat
Plate 3.4047.603 trial
Plantilla para placa 3.4047.603

43.4047.603



			Ti
	Len	L	R
2	29	3.4041.602	3.4042.602
3	40	3.4041.603	3.4042.603
4	51	3.4041.604	3.4042.604

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)

40.5713.040

Platte 3.4041.603 Probeimplantat
Plate 3.4041.603 trial
Plantilla para placa 3.4041.603
43.4041.603
Platte 3.4042.603 Probeimplantat
Plate 3.4042.603 trial
Plantilla para placa 3.4042.603
43.4042.603



1:1

			Ti
	Len	L	R
2	30	3.4037.602	3.4038.602
3	41	3.4037.603	3.4038.603
4	52	3.4037.604	3.4038.604

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)

40.5713.040




Platte 3.4037.603 Probeimplantat
Plate 3.4037.603 trial
Plantilla para placa 3.4037.603
Platte 3.4038.603 Probeimplantat
Plate 3.4038.603 trial
Plantilla para placa 3.4038.603

43.4037.603
43.4038.603



1:1



Ti

			
2	32	3.4045.602	3.4046.602
3	43	3.4045.603	3.4046.603
4	54	3.4045.604	3.4046.604

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4




Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)

40.5713.040



Platte 3.4045.603 Probeimplanta	
Plate 3.4045.603 trial	43.4045.603
Plantilla para placa 3.4045.603	
Platte 3.4046.603 Probeimplantat	
Plate 3.4046.603 trial	43.4046.603
Plantilla para placa 3.4046.603	

1:1




3	40	3.4048.603	
4	51	3.4048.604	

3.5164.xxx	✓	✓	✓	✓	✓	2.4		
3.5165.xxx	✓	✓	✓	✓	✓	2.7		
3.1220.xxx	✓	✓	✓	✓	✓	2.7		
4.5235.xxx	✓	✓	✓	✓	✓	2.4		



Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)


40.5713.040



Platte 3.4048.603 Probeimplantat
Plate 3.4048.603 trial
Plantilla para placa 3.4048.603

43.4048.603

1:1



			Ti	
	Len	L		R
2	29	3.4043.602		3.4044.602
3	40	3.4043.603		3.4044.603
4	51	3.4043.604		3.4044.604

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4




Palette für 4,0 ChLP Implantate - distaler Radius
(mit Instrumenten)
Palette for 4.0ChLP implants - distal radius
(with instruments)
Paleta para implantes 4,0ChLP - radio distal (con
instrumentos)


40.5713.040




Platte 3.4043.603 Probeimplantat
Plate 3.4043.603 trial
Plantilla para placa 3.4043.603
43.4043.603
Platte 3.4044.603 Probeimplantat
Plate 3.4044.603 trial
Plantilla para placa 3.4044.603
43.4044.603



	Len	Ti
3	45	3.4099.603
4	53	3.4099.604
5	61	3.4099..605

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



40.5710.171	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Instrumente 4,0ChLP Instrument set for 4.0ChLP Set de instrumental para 4,0ChLP	1	40.5711.200
	Palette für 4,0ChLP Implantate - distal Radius* Palette for 4.0ChLP implants - distal radius* Paleta implantes 4,0ChLP-radio distal*	1	40.5713.170
	Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum lid 1/2 306x272x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	1	12.0751.200
	Behälter mit solidem Boden 1/2 306x272x85mm Container with solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	1	12.0751.100

* Behälter ohne Implantaten / Tray does not include implants / Bandeja no incluye implantes

40.5713.050

Name
Name
Nombre

Stk. Kat.-Nr.,
Pcs Catalogue no.
Uds. Número de catálogo



Palette für 4,0 ChLP Implantate - distaler Radius
Palette for 4.0ChLP implants - distal radial
Paleta para implantes 4,0ChLP - radio distal

1

40.5713.150



Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm
Perforated aluminum cover 1/2 306x272x15mm gray
Tapa de aluminio perforada 1/2 306x272x15mm Gris

1

12.0751.200

Behälter mit solidem Boden 1/2 306x272x85mm
Container with solid bottom 1/2 306x272x85mm
Contenedor de fondo sólido 1/2 306x272x85mm


1

12.0751.100


40.5711.200	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4,0 Guía roscada M3,5/1,8 - 4,0	4	40.4896.018
	Kompressionsführung 1,8 Compression guide 1.8 Guía de compresión 1,8	1	40.4897.018
	Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	1	40.5928.018
	Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja de Kirschner 1,0/180	5	40.4814.000
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	2	40.2063.181
	Messlehre für Verriegelungsschrauben- länge Locking screw length measure Medidor de tornillos de bloqueo	1	40.4818.100
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
	Schraubendrehereinsatz T8,0 Screwdriver tip T8.0 Punta del destornillador T8,0	1	40.5682.000
	Einsatz T8 mit Greifer Screwdriver tip T8 with holder Punta del destornillador T8 con sujeción	1	40.5989.000
	Gewindebohrer HA 2,7 Cortical tap HA 2.7 Macho de roscar HA 2,7	1	40.5988.000
	Gewindebohrer 4,0ChLP - 2,4 Tap 4.0ChLP - 2.4 Macho de roscar 4,0ChLP - 2,4	1	40.5987.024
	Plattenrepositionsschraube 1,8/120 Setting-compressing screw 1.8/120 Tornillo de ajuste-compresión 1,8/120	2	40.5678.000
	Drehmomentgriff mit Kupplung 1Nm Torque limiting ratchet handle 1.0Nm Mango dinamométrico con carraca 1Nm	1	40.6650.000
	Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	1	40.0669.100
	Biegegerät für Platten 4,0 Plates bender 4,0 Grifa para doblar placas 4,0	2	40.4643.000
	Anatomische Pinzette Standard, 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas Standard 14,5cm	1	30.3303.000
	Palette für Instrumente 4,0ChLP Palette for instruments 4.0ChLP Paleta para instrumentos 4,0ChLP	1	40.5712.100




①



②



①	6	14	3.7206.014
②	12	21	3.7206.021




	TiA	Co	Locking	VA	Locking	Locking
3.5164.xxx	✓		✓	✓	✓	2.4
4.5235.xxx		✓	✓	✓	✓	2.4

Behälter für Implantate 4,0ChLP 3.7206 4x2 1/2H
Stand for 4.0ChLP implants 3.7206 4x2 1/2H
Contenedor para implantes 4,0ChLP 3.7206 4x2 1/2H

15.0204.602

	Name Name Nombre	Kat.-Nr Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Fräser 15 Cutter 15 Fresa 15	40.8262.015	1
	Fräser 22 Cutter 22 Fresa 22	40.8262.022	1
	Kirschnerdraht 1,5/180 Kirschner wire 1.5/180 Aguja Kirschner 1,5/180	40.4592.180	5







Behälter für Implantate 4,0ChLP 3.7206 4x2 1/2H
Stand for 4.0ChLP implants 3.7206 4x2 1/2H
Contenedor para implantes 4,0ChLP 3.7206 4x2 1/2H

14.0204.602

1



Behälter für Implantate 4,0ChLP 3.7206 4x2 1/2H
Stand for 4.0ChLP implants 3.7206 4x2 1/2H
Contenedor para implantes 4,0ChLP 3.7206 4x2 1/2H

14.0204.602

	4,0ChLP Handgelenkplatte für Arthrodese 4.0ChLP wrist arthrodesis plate 4,0ChLP placa para artrodesis de muñeca			
	Stk. Pcs Uds.		2	2
	4,0ChLP Schraube 2,4 4.0ChLP Screw 2.4 4,0ChLP tornillo 2,4			
	L [mm]	6 8 10 12 14 16 18 20		
	Stk. Pcs Uds.	5 5 5 5 5 5 5 5		
	4,0ChLP VA Schraube 2,4 4.0ChLP screw VA 2.4 4,0ChLP tornillo VA 2,4			
	L [mm]	6 8 10 12 14 16 18 20		
	Stk. Pcs Uds.	5 5 5 5 5 5 5 5		

* Behälter ohne Implantaten / Trav does not include implants / Bandeia no incluye implantes

* Behälter ohne Implantaten / Tray does not include implants / Bandeja no incluye implantes

	Name Name Nombre	Kat.-Nr. Catalogue no. Número de catálogo	Stk. Pcs Uds.	
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4,0 Guía roscada M3,5/1,8 - 4,0	40.4896.018	4	40.5711.200
	Kompressionsführung 1,8 Compression guide 1.8 Guía de compresión 1,8	40.4897.018	1	
	Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	40.5928.018	1	
	Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja Kirschner 1,0/180	40.4814.000	5	
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	40.2063.181	2	
	Längenmesslehre für Verriegelungsschraubenlänge Length measure of locking screw Medidor de tornillos de bloqueo	40.4818.100	1	
	Tiefenmesslehre Depth measure Medidor de profundidad	40.4640.000	1	
	Schraubendrehereinsatz T 8,0 Screwdriver tip T8.0 Punta de destornillador T8,0	40.5682.000	1	
	Einsatz T8 mit Greifer T8 screwdriver tip with holder Punta de destornillador T8 con sujeción	40.5989.000	1	
	Gewindebohrer HA 2,7 Cortical tap HA 2.7 Macho de roscar HA 2,7	40.5988.000	1	
	Gewindebohrer 4,0ChLP - 2,4 Tap 4.0ChLP - 2.4 Macho de roscar 4,0ChLP - 2,4	40.5987.024	1	
	Plattenrepositionsschraube 1,8/120 Setting-compressing screw 1.8/120 Tornillo de ajuste-compresión 1,8/120	40.5678.000	2	
	Handgriff mit Drehmomentbegrenzung 1,0Nm Torque limiting ratchet handle 1.0Nm Mango dinamométrico 1,0Nm	40.6650.000	1	
	Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	40.0669.100	1	
	Plattenbiegegerät 4,0 Plates bender 4.0 Grifa para doblar placas 4,0	40.4643.000	2	40.5711.300
	Anatomische Pinzette Standard 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas STANDARD, 14,5cm	30.3303.000	1	
	Palette für Instrumente 4,0ChLP Palette for instruments 4.0ChLP Paleta para instrumentos 4,0ChLP	40.5712.100	1	
	Behälter mit solidem Boden 1/2 306x272x85mm Container with solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	12.0751.100	1	
	Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum lid 1/1 595x275x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	12.0751.200	1	




		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4

			Ti	
1	Len	L	R	
1	36	3.7025.501	3.7024.501	
2	47	3.7025.502	3.7024.502	
2				
1	36	3.7027.501	3.7026.501	
2	47	3.7027.502	3.7026.502	
1-3	36-58			


KLEINE
SMALL
PEQUEÑA

GROßE
LARGE
LARGA




Palette für 4,0ChLP Instrumente und Implantate -
-Radiuskopf
Palette for 4.0ChLP instruments and implants-radial
head bone
Paleta para instrumentos e implantes de 4,0ChLP
para la cabeza de radio

40.5713.330





Platte 3.7024.502 Probeimplantat Plate 3.7024.502 trial Plantilla para placa 3.7024.502	43.7024.502
Platte 3.7025.502 Probeimplantat Plate 3.7025.502 trial Plantilla para placa 3.7025.502	43.7025.502
Platte 3.7026.502 Probeimplantat Plate 3.7026.502 trial Plantilla para placa 3.7026.502	43.7026.502
Platte 3.7027.502 Probeimplantat Plate 3.7027.502 trial Plantilla para placa 3.7027.502	43.7027.502



1:1

			
1	32		3.7028.501
2	43		3.7028.502
			
1-3	32	54	

								
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	3.1220.xxx	✓		✓			✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4





Palette für 4,0ChLP Instrumente und Implantate -
-Radiuskopf
Palette for 4.0ChLP instruments and implants-radial
head bone
Paleta para instrumentos e implantes de 4,0ChLP
para la cabeza de radio

40.5713.330



Platte 3.7028.502 Probeimplantat
Plate 3.7028.502 trial
Plantilla para placa 3.7028.502

43.7028.502


40.5713.230	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4,0 Guía roscada M3,5/1,8 - 4,0	4	40.4896.018
	Kompressionsführung 1,8 Compression guide 1.8 Guía de compresión 1,8	1	40.4897.018
	Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	1	40.5928.018
	Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja de Kirschner 1,0/180	5	40.4814.000
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	2	40.2063.181
	Messlehre für Verriegelungsschraubenlänge Locking screw length measure Medidor de tornillos de bloqueo	1	40.4818.100
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
	Schraubendrehereinsatz T8,0 Screwdriver tip T8.0 Punta del destornillador T8,0	1	40.5682.000
	Drehmomentgriff mit Kupplung 1Nm Torque limiting ratchet handle 1.0Nm Mango dinamométrico con carraca 1Nm	1	40.6650.000
	Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	1	40.0669.100
	Biegegerät für Platten 4,0 Plates bender 4,0 Grifa para doblar placas 4,0	2	40.4643.000
	Anatomische Pinzette Standard, 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas Standard 14,5cm	1	30.3303.000
	Palette für 4,0ChLP Instrumente und Implantate - Radiuskopf Palette for 4.0ChLP instruments and implants-radial head bone Paleta para instrumentos e implantes de 4,0ChLP para la cabeza de radio	1	40.5713.330
	Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum cover 1/2 306x272x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	1	12.0751.200
	Behälter mit solidem Boden 1/2 306x272x85mm Container with solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	1	12.0751.100



1:1

Len	Ti
41	3.7021.041

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie



40.5714.020

1:1



Len	Ti
22	3.4061.022
27	3.4061.027

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7
	3.1220.xxx	✓			✓	2.7



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie



40.5714.020



1:1

Ti		
Len	L	R
40	3.4063.040	3.4064.040
45	3.4063.045	3.4064.045
50	3.4063.050	3.4064.050

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




		Ti	
	Len	L	R
4	39	3.4066.004	
5	46	3.4073.005	3.4066.005
6	54	3.4066.006	

	TiA			VA		
3.5164.xxx	✓	✓		✓	2.4	
3.5165.xxx	✓	✓		✓	2.7	





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




	Ti	
L ₁ xL ₂	L	R
15x26	3.4067.026	3.4068.026
18x29	3.4067.029	3.4068.029

	TiA			VA		
	3.5164.xxx	✓	✓	✓	✓	2.4
	3.5165.xxx	✓	✓	✓	✓	2.7





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




Len	Ti
53	3.4069.053

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie


40.5714.020



1:1

Ti		
Len	L	R
46	3.4071.046	3.4070.046
50	3.4071.050	3.4070.050
54	3.4071.054	3.4070.054

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020



Ti		
Len	L	R
33	3.4085.033	3.4086.033

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




		Ti		
	Len	L	R	
2	28	3.4087.028	3.4088.028	
3	38	3.4087.038	3.4088.038	
4	48	3.4087.048	3.4088.048	

	TiA			VA		
3.5164.xxx	✓	✓		✓	2.4	
3.5165.xxx	✓	✓		✓	2.7	
3.1220.xxx	✓			✓	2.7	





Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




Len	Ti
51	3.4084.051
60	3.4084.060
66	3.4084.066

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7
	3.1220.xxx	✓			✓	2.7




Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020




Ti		
Len	L	R
39	3.4075.601	3.4074.601
48	3.4075.602	3.4074.602
59	3.4075.603	3.4074.603

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7
	3.1220.xxx	✓			✓	2.7



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020





1:1



	Len	Ti
3	54	3.4072.503
4	64	3.4072.504
5	74	3.4072.505
6	84	3.4072.506
7	94	3.4072.507
8	104	3.4072.508
9	114	3.4072.509
10	124	3.4072.510

	TiA			VA		
	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7
	3.1220.xxx	✓			✓	2.7



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

40.5714.020

40.5710.020


Name
Name
Nombre

Stk. Kat.-Nr.
Pcs Catalogue no.
Uds. Número de catálogo



Instrumente 4,0ChLP
Instruments set for 4.0ChLP
Set de instrumental para 4,0ChLP

1 40.5711.000



Palette für 4,0ChLP Implantate - Fuß
Palette for 4.0ChLP implants - foot
Paleta para implantes 4,0ChLP - pie

1 40.5714.020




Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm
Perforated aluminum cover 1/2 306x272x15mm gray
Tapa de aluminio perforada 1/2 306x272x15mm Gris

1 12.0751.200

Behälter mit solidem Boden 1/2 306x272x85mm
Container with solid bottom 1/2 306x272x85mm
Contenedor de fondo sólido 1/2 306x272x85mm


1 12.0751.100

40.5711.000	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4.0 Guía roscada M3,5/1,8 - 4,0	4	40.4896.018
	Kompressionsführung 1,8 Compression guide 1.8 Guía de compresión 1,8	1	40.4897.018
	Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja de Kirschner 1,0/180	5	40.4814.000
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	2	40.2063.181
	Messlehre für Verriegelungsschraubenlänge Locking screw length measure Medidor de tornillos de bloqueo	1	40.4818.000
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
	Schraubendrehereinsatz T8,0 Screwdriver tip T8.0 Punta del destornillador T8,0	1	40.5682.000
	Drehmomentgriff mit Kupplung 1Nm Torque limiting ratchet handle 1.0Nm Mango dinamométrico con carraca 1Nm	1	40.6650.000
	Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	1	40.0669.100
	Biegegerät für Platten 4,0 Plates bender 4,0 Grifa para doblar placas 4,0	2	40.4643.000
	Anatomische Pinzette Standard, 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas Standard 14,5cm	1	30.3303.000
	Palette für Instrumente 4,0ChLP Palette for 4.0ChLP instruments Paleta para instrumentos de 4,0ChLP	1	40.5712.000



A	Len	Ti
10	16	3.7020.016
14	20	3.7020.020
18	24	3.7020.024

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für Instrumente und Implantate 4,0ChLP – Hallux
Palette for 4.0ChLP implants and instruments - hallux
Paleta para instrumentos e implantes de 4,0ChLP - hallux

40.5741.300

1.8

Ster
Non
Ster

L R



1:1

3.4031.032

Ti

Len	L	R
28	3.4031.028	3.4032.028
30	3.4031.030	3.4032.030
32	3.4031.032	3.4032.032

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für Instrumente und Implantate 4,0ChLP – Hallux
Palette for 4.0ChLP implants and instruments - hallux
Paleta para instrumentos e implantes de 4,0ChLP - hallux



40.5741.300



1:1

Len	Ti
23	3.4065.723
25	3.4065.725
30	3.4065.730
35	3.4065.735

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4




Palette für Instrumente und Implantate 4,0ChLP – Hallux
Palette for 4,0ChLP implants and instruments - hallux
Paleta para instrumentos e implantes de 4,0ChLP - hallux

40.5741.300



Len	Ti
38	3.4062.038
40	3.4062.040
42	3.4062.042
44	3.4062.044
46	3.4062.046


		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für Instrumente und Implantate 4,0ChLP – Hallux
Palette for 4.0ChLP implants and instruments - hallux
Paleta para instrumentos e implantes de 4,0ChLP - hallux


40.5741.300

40.5741.200	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4.0 Guía roscada M3,5/1,8 - 4,0	4	40.4896.018
	Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	1	40.5928.018
	Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja de Kirschner 1,0/180	5	40.4814.000
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	2	40.2063.181
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
	Schraubendrehereinsatz T8,0 Screwdriver tip T8.0 Punta del destornillador T8,0	1	40.5682.000
	Drehmomentgriff mit Kupplung 1Nm Torque handle 1.0Nm Mango dinamométrico con carraca 1Nm	1	40.6650.000
	Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	1	40.0669.100
	Anatomische Pinzette Standard, 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas Standard 14,5cm	1	30.3303.000
	Repositionszange, schräg 175mm Bone holding forceps oblique 175mm Pinzas para huesos, oblicuas 175mm	1	40.4146.000
	Palette für Instrumente und Implantate 4,0ChLP - Hallux Palette for 4.0ChLP instruments and implants - hallux Paleta para instrumentos e implantes de 4,0ChLP - hallux	1	40.5741.300
	Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum cover 1/2 306x272x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	1	12.0751.200
	Behälter mit solidem Boden 1/2 306x272x85mm Container with solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	1	12.0751.100




		
2	30	3.7056.002
3	30	3.7056.003
4	30	3.7056.004
5	30	3.7056.005
6	30	3.7056.006
7	30	3.7056.007

								
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4




Palette für Instrumente und Implantate 4,0ChLP-Keilosteotomie
Pallette for 4.0ChLP implants and instruments- opening wedge osteotomy
Paleta para instrumental e implantes 4,0ChLP - osteotomía en cuña

40.6298.000



	Len	Ti
2	23	3.7057.002
3	23	3.7057.003
4	23	3.7057.004
5	23	3.7057.005
6	23	3.7057.006
7	23	3.7057.007

		TiA	Co			VA		
	3.5164.xxx	✓		✓	✓		✓	2.4
	3.5165.xxx	✓		✓	✓		✓	2.7
	4.5235.xxx		✓	✓	✓	✓	✓	2.4



Palette für Instrumente und Implantate 4,0ChLP-Keilosteotomie
Pallette for 4.0ChLP implants and instruments- opening wedge osteotomy
Paleta para instrumental e implantes 4,0ChLP - osteotomía en cuña

40.6298.000


40.6297.000		Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
		Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3,5/1,8 - 4,0 Guía roscada M3,5/1,8 - 4,0	4	40.4896.018
		Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	1	40.5928.018
		Kirschnerdraht 1,0/180 Kirschner wire 1.0/180 Aguja Kirschner 1,0/180	5	40.4814.000
		Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	2	40.2063.181
		Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
		Schraubendrehereinsatz T8,0 Screwdriver tip T8,0 Punta de destornillador T8,0	1	40.5682.000
		Handgriff mit Drehmomentbegrenzung 1,0Nm Torque limiting ratchet handle 1.0Nm Mango de torsión 1,0Nm	1	40.6650.000
		Sechskant-Schraubendreher T8 Star screwdriver T8 Destornillador de estrella T8	1	40.0669.100
		Anatomische Pinzette Standard, 14,5cm Dissecting forceps Standard 14.5cm Pinzas anatómicas STANDARD, 14,5cm	1	30.3303.000
		Repositionszange, schräg 175mm Sarf bone holding forceps 175 mm Pinzas para huesos oblicuas 175mm	1	40.4146.000
		Palette für Instrumente und Implantate 4,0ChLP-Keilosteotomie Pallette for 4.0ChLP implants and instruments- opening wedge osteotomy Paleta para instrumental e implantes 4,0ChLP - osteotomía en cuña	1	40.6298.000
		Behälter mit solidem Boden 1/2 306x272x85mm Container solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	1	12.0751.100
		Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum lid 1/2 306x272x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	1	12.0751.200







Ti		
① Len	L	R
40	3.7061.140	3.7060.140
45	3.7061.145	3.7060.145
50	3.7061.150	3.7060.150


② Len	L	R
40	3.7061.240	3.7060.240
45	3.7061.245	3.7060.245
50	3.7061.250	3.7060.250

	TiA	Co	⊗	⊙	VA	⌚	↗
3.5164.xxx	✓		✓	✓		✓	2.4
4.5235.xxx		✓	✓	✓	✓	✓	2.4



15.0204.002	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Instrumentarium für 4,0ChLP 3.7060/3.7061 4x4 1/2H Instrument set for 4.0ChLP 3.7060/3.7061 4x4 1/2H Set de instrumental para 4,0ChLP 3.7060/3.7061 4x4 1/2H	1	15.0204.201
	Instrumentarium für 4,0ChLP 3.7060/3.7061 4x2 1/2H Instrument set for 4.0ChLP 3.7060/3.7061 4x2 1/2H Set de instrumental para 4,0ChLP 3.7060/3.7061 4x2 1/2H	1	15.0204.203
	Behälter für Implantate 4,0ChLP 3.7060/3.7061 4x2 1/2H Stand for 4.0ChLP implants 3.7060/3.7061 4x2 1/2H Contenedor para implantes 4,0ChLP 3.7060/3.7061 4x2 1/2H	1	14.0204.601
	Behälterdeckel 4,0ChLP 3.7060/3.7061 4x4 4.0ChLP container lid 3.7060/3.7061 4x4 Tapa de contenedor 4,0ChLP 3.7060/3.7061 4x4	1	14.0204.104
	Behälter 4,0ChLP 3.7060/3.7061 4x4 4.0ChLP container 3.7060/3.7061 4x4 Contenedor 4,0ChLP 3.7060/3.7061 4x4	1	14.0204.103

15.0204.201	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Gewindeführung M3,5/1,8 - 4,0 Threaded guide M3.5/1.8 - 4.0 Guía roscada M3,5/1,8 - 4,0	2	40.4896.018
	Bohrer 1,8/180 Drill 1.8/180 Broca 1,8/180	1	40.2063.181
	Schraubendrehereinsatz T8,0 Screwdriver tip T8.0 Punta del destornillador T8,0	1	40.5682.000
	Einsatz T8 mit Greifer Screwdriver tip T8 with holder Punta del destornillador T8 con sujeción	1	40.5989.000
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4640.000
	Messlehre für Verriegelungsschraubenlänge Locking screw length measure Medidor de tornillos de bloqueo	1	40.4818.100
	Buck-Gramcko Hebel 7,5 BUCK-GRAMCKO elevator 7.5 Elevador BUCK-GRAMCKO 7,5	2	40.2185.000
	Schlitzhammer Mallet Martillo	1	40.6284.000
	Ausschlaginstrument Extractor Extractor	1	40.6283.000
	Drehmomentgriff mit Kupplung 1Nm Torque limiting ratchet handle 1Nm Mango dinamométrico con carraca 1Nm	1	40.6650.000
	Führung VA 1,8 Guide VA 1.8 Guía VA 1,8	1	40.5928.018
	Palette für Instrumente 4,0ChLP 3.7060/3.7061 4x4 1/2H Tray for 4.0ChLP instrument set 3.7060/3.7061 4x4 1/2H Paleta para instrumental 4,0ChLP 3.7060/3.7061 4x4 1/2H	1	14.0204.201

15.0204.203	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Raspatorium Raspatory Periostótomo	1	40.6285.000
	Zielgerät für Interknochenplatte - links Targeter for endosteal plate - left Guía para la placa intraósea - izquierda	1	40.6281.100*
	Zielgerät für Interknochenplatte - rechts Targeter for endosteal plate - right Guía para la placa intraósea - derecha	1	40.6282.100**
	Bohrer 4,0 Wiertło 4,0 Broca 4,0	1	40.6278.000
	Palette für Instrumentarium für 4,0ChLP 3.7060/3.7061 4x2 1/2H Tray for 4.0ChLP instrument set 3.7060/3.7061 4x2 1/2H Paleta para instrumental 4,0ChLP 3.7060/3.7061 4x2 1/2H	1	14.0204.203

* verwenden mit Platten/ use with plates/ usar con las placas 3.7061.1xx; 3.7061.2xx

** verwenden mit Platten/ use with plates/ usar con las placas 3.7060.1xx;
3.7060.2xx



	3.5164.xxx	✓	✓		✓	2.4
	3.5165.xxx	✓	✓		✓	2.7
	3.1220.xxx	✓			✓	2.7



	4.5235.xxx	✓	✓	✓	✓	2.4



4,0ChLP SCHRAUBE 2,4
4.0ChLP SCREW 2.4
4,0ChLP TORNILLO 2,4



6	3.5164.006
8	3.5164.008
10	3.5164.010
12	3.5164.012
14	3.5164.014
16	3.5164.016
18	3.5164.018
20	3.5164.020
22	3.5164.022
24	3.5164.024
26	3.5164.026
28	3.5164.028
30	3.5164.030
32	3.5164.032
34	3.5164.034
36	3.5164.036
38	3.5164.038
40	3.5164.040

Kerndurchmesser core diameter diámetro de núcleo	1.8	
Bohrer drill broca	1.8	40.2063.181
Gewindeführung threaded guide guía roscada	M3.5/1.8	40.4896.018
Schraubendrehereinsatz screwdriver tip punta de destornillador	T8	40.5682.000

4,0ChLP SCHRAUBE 2,7
4.0ChLP SCREW 2.7
4,0ChLP TORNILLO 2,7



6	3.5165.006
8	3.5165.008
10	3.5165.010
12	3.5165.012
14	3.5165.014
16	3.5165.016
18	3.5165.018
20	3.5165.020
22	3.5165.022
24	3.5165.024
26	3.5165.026
28	3.5165.028
30	3.5165.030
32	3.5165.032
34	3.5165.034
36	3.5165.036
38	3.5165.038
40	3.5165.040

Kerndurchmesser core diameter diámetro de núcleo	1.8	
Bohrer drill broca	1.8	40.2063.181
Gewindeführung threaded guide guía roscada	M3.5/1.8	40.4896.018
Schraubendrehereinsatz screwdriver tip punta de destornillador	T8	40.5682.000



SELBSTSCHNEIDENDE KORTIKALISSCHRAUBE 2,7
CORTICAL SELF-TAPPING SCREW 2.7
TORNILLO CORTICAL AUTO-ROSCANTE 2,7



6	3.1220.006
8	3.1220.008
10	3.1220.010
12	3.1220.012
14	3.1220.014
16	3.1220.016
18	3.1220.018
20	3.1220.020
22	3.1220.022
24	3.1220.024
26	3.1220.026
28	3.1220.028
30	3.1220.030
32	3.1220.032
34	3.1220.034
36	3.1220.036
38	3.1220.038
40	3.1220.040

Kerndurchmesser core diameter diámetro de núcleo	1.8	
Bohrer drill broca	1.8	40.2063.181
Kompressionsführung compression guide guía de compresión	1.8	40.4897.018
Sechskant-Schraubendreher star screwdriver destornillador de estrella	T8	40.0669.000

4,0ChLP VA SCHRAUBE 2,4
4.0ChLP SCREW VA 2.4
4,0ChLP TORNILLO VA 2,4



6	4.5235.006
8	4.5235.008
10	4.5235.010
12	4.5235.012
14	4.5235.014
16	4.5235.016
18	4.5235.018
20	4.5235.020
22	4.5235.022
24	4.5235.024
26	4.5235.026
28	4.5235.028
30	4.5235.030
32	4.5235.032
34	4.5235.034
36	4.5235.036
38	4.5235.038
40	4.5235.040

Kerndurchmesser core diameter diámetro de núcleo	1.8	
Bohrer drill broca	1.8	40.2063.181
Führung VA guide VA guía VA	1.8	40.5928.018
Schraubendrehereinsatz screwdriver tip punta de destornillador	T8	40.5682.000



	4	88	3.4080.504
	5	102	3.4080.505
	6	116	3.4080.506
	7	130	3.4080.507
	8	144	3.4080.508



		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				




3.4022.806



4	88	3.4022.804
5	102	3.4022.805
6	116	3.4022.806
7	130	3.4022.807
8	144	3.4022.808
10	172	3.4022.810



	3.5200.xxx	✓	✓	✓	✓	3.5
	3.5205.xxx	✓	✓		✓	3.5
	3.5231.xxx	✓	✓	✓		3.9
	3.5225.xxx	✓	✓	✓	✓	2.4
	3.1306.xxx	✓	✓		✓	3.5
	4.5236.xxx		✓	✓	✓	3.5
	3.5207.006	✓	✓	✓		3.5
	3.1221.050	✓	✓	✓		



Palette für 5,0ChLP Platten 3.4022/3.4080 4x4 1/2H
Tray for 5.0ChLP plates 3.4022/3.4080 4x4 1/2H
Paleta para placas 5,0ChLP 3.4022/3.4080 4x4 1/2H


14.0205.402



	Len	Ti
5	88	3.3157.505
6	103	3.3157.506
7	118	3.3157.507
8	133	3.3157.508
10	163	3.3157.510
12	193	3.3157.512



	3.5200.xxx	✓	✓	✓	✓	3.5
	3.5205.xxx	✓	✓		✓	3.5
	3.5231.xxx	✓	✓	✓	✓	3.9
	3.5225.xxx	✓	✓	✓	✓	2.4
	3.1306.xxx	✓	✓		✓	3.5
	4.5236.xxx	✓	✓	✓	✓	3.5
	3.5207.006	✓	✓	✓		3.5
	3.1221.050	✓	✓	✓		



Palette für 5,0ChLP Platten 3.3157-3.3159 4x4 1/2H
Tray for 5.0ChLP plates 3.3157-3.3159 4x4 1/2H
Paleta para placas 5,0ChLP 3.3157-3.3159 4x4 1/2H

14.0205.404



	Len	Ti
4	45	3.4081.004
5	57	3.4081.005
6	69	3.4081.006
7	81	3.4081.007
8	93	3.4081.008
10	117	3.4081.010
3-12	33-141	


		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				






	Len	Ti
4	46	3.4081.604
5	58	3.4081.605
6	70	3.4081.606
7	82	3.4081.607
8	94	3.4081.608
9	106	3.4081.609
10	118	3.4081.610

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				





10	3.7242.110
13	3.7242.113
16	3.7242.116

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5



15.0205.601	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Skaliertes Bohrer 2,8/210 Drill with scale 2.8/210 Broca graduada 2,8/210	1	40.5653.212
	Schraubendrehereinsatz T15 Screwdriver tip T15 Punta de destornillador T15	1	40.5677.000
	Knebelgriff mit Kupplung Handle ratchet device Mango con embrague	1	40.6654.000
	Führungshülse 5,0/2,8 Guide sleeve 5.0/2.8 Guía de la broca 5,0/2,8	1	40.5673.728
	Messlehre Measure Plantilla	1	40.8263.000
	Behälter f. Instrumentarium 5,0ChLP 3.7242 4x2 1/2H Stand for 5.0ChLP implants 3.7242 4x2 1/2 Contenedor para instr. 5,0ChLP 3.7242 4x2 1/2H	1	14.0205.601


INSTRUMENTARIUM 5,0CHLP 3.7242 4X1 1/2H* STAND FOR 5.0CHLP IMPLANTS 3.7242 4X1 1/2H* CONTENEDOR PARA TORNILLOS 5,0CHLP 3.7242 4X1 1/2H*

15.0205.430	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Messlehre Measure Plantilla	1	40.8263.000
	Behälter für Implantate 5,0ChLP 3.7242 4x1 1/2H Stand for 5.0ChLP implants 3.7242 4x1 1/2 Contenedor para implantes 5,0ChLP 3.7242 4x1 1/2H	1	14.0205.430

*Ergänzendes Instrumentarium. Verwendung mit Instrumentarium 5.0ChLP 4x4 1/2H [15.0205.206].


*Complementary instruments. Use with instrument set for 5.0ChLP 4x1 1/2H [15.0205.206].

*Instrumentación complementaria. Utilizar con Instrumentarium 5.0ChLP 4x4 1/2H [15.0205.206].



			(Ti)	
	(H)	(R)	(L)	(R)
12	5		3.4025.505	3.4026.505
	6		3.4025.506	3.4026.506
	7		3.4025.507	3.4026.507
15	5		3.4027.505	3.4028.505
	6		3.4027.506	3.4028.506
	7		3.4027.507	3.4028.507

		(TiA)	(Co)			(VA)			
3.5200.xxx	✓		✓	✓		✓			3.5
3.5205.xxx	✓		✓			✓			3.5
3.5231.xxx	✓		✓	✓				✓	3.9
3.5225.xxx	✓		✓	✓		✓			2.4
3.1306.xxx	✓		✓			✓			3.5
4.5236.xxx		✓	✓	✓	✓	✓			3.5
3.5207.006	✓		✓	✓					3.5



Palette für 5,0ChLP Platten 3.4025-3.4030 4x4 1/2H
Tray for plates 5,0ChLP 3.4025-3.4030 4x4 1/2H
Paleta para placas 5,0ChLP 3.4025-3.4030 4x4 1/2H
14.0205.403





Platte 3.4025.505 Probeimplantat Plate 3.4025.505 trial Plantilla para placa 3.4025.505	43.4025.505
Platte 3.4026.505 Probeimplantat Plate 3.4026.505 trial Plantilla para placa 3.4026.505	43.4026.505
Platte 3.4027.505 Probeimplantat Plate 3.4027.505 trial Plantilla para placa 3.4027.505	43.4027.505
Platte 3.4028.505 Probeimplantat Plate 3.4028.505 trial Plantilla para placa 3.4028.505	43.4028.505



3.7015.508

		Ti	
		L	R
3	71	3.7015.503	3.7014.503
4	80	3.7015.504	3.7014.504
5	90	3.7015.505	3.7014.505
6	99	3.7015.506	3.7014.506
8	116	3.7015.508	3.7014.508
3-11	71-143		

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5



Palette für 5,0ChLP Platten 3.7014/3.7015 4x4 1/2H
Tray for 5.0ChLP plates 3.7014/3.7015 4x4 1/2H
Paleta para placas 5,0ChLP 3.7014/3.7015 4x4 1/2H

14.0205.405



Platte 3.7014.505 Probeimplantat
Plate 3.7014.505 trial
Plantilla para placa 3.7014.505
Platte 3.7015.505 Probeimplantat
Plate 3.7015.505 trial
Plantilla para placa 3.7015.505



43.7014.505

43.7015.505



			
6	66	3.7049.506	3.7048.506
8	84	3.7049.508	3.7048.508
10	104	3.7049.510	3.7048.510
			
5-11	57-114		

									
	3.5200.xxx	✓	✓	✓		✓			3.5
	3.5205.xxx	✓	✓			✓			3.5
	3.5231.xxx	✓	✓	✓			✓		3.9
	3.5225.xxx	✓	✓	✓		✓			2.4
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓	✓	✓					3.5



Palette für 5,0ChLP Platten 3.7048/3.7049 4x4 1/2H
Tray for 5.0ChLP plates 3.7048/3.7049 4x4 1/2H
Paleta para placas 5,0ChLP 3.7048/3.7049 4x4 1/2H

14.0205.406





Platte 3.7048.508 Probeimplantat
Plate 3.7048.508 trial
Plantilla para placa 3.7048.508
43.7048.508
Platte 3.7049.508 Probeimplantat
Plate 3.7049.508 trial
Plantilla para placa 3.7049.508
43.7049.508



3.7207.605

			Ti	
	Len	L		R
3	72	3.7207.603		3.7208.603
5	91	3.7207.605		3.7208.605

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5



Palette für 5,0ChLP Platten 3.7207/3.7212 4x4 1/2H
Tray for plates 5.0ChLP 3.7207-3.7212 4x4 1/2H
Paleta para placas 5,0ChLP 3.7207/3.7212 4x4 1/2H


14.0205.426



3.7209.607

			Ti	
	Len	L		R
5	91	3.7209.605		3.7210.605
7	111	3.7209.607		3.7210.607

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5



Palette für 5,0ChLP Platten 3.7207/3.7212 4x4 1/2H
Tray for plates 5,0ChLP 3.7207-3.7212 4x4 1/2H
Paleta para placas 5,0ChLP 3.7207/3.7212 4x4 1/2H

14.0205.426




39

3.7211.600

3.7212.600

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5





Palette für 5,0ChLP Platten 3.7207/3.7212 4x4 1/2H
 Tray for plates 5.0ChLP 3.7207-3.7212 4x4 1/2H
 Paleta para placas 5,0ChLP 3.7207/3.7212 4x4 1/2H

14.0205.426



		TiA	Co			VA				
	3.5200.xxx	✓		✓	✓		✓			3.5
	3.5205.xxx	✓		✓			✓			3.5
	3.5231.xxx	✓		✓	✓				✓	3.9
	3.5225.xxx	✓		✓	✓		✓			2.4
	3.1306.xxx	✓		✓			✓			3.5
	4.5236.xxx		✓	✓	✓	✓	✓			3.5
	3.5207.006	✓		✓	✓					3.5
	3.1221.050	✓		✓	✓					




40.5671.000



			Ti	
	!	Len	L	R
5		116	3.7243.605	3.7244.605
6		129	3.7243.606	3.7244.606
8		155	3.7243.608	3.7244.608
10		181	3.7243.610	3.7244.610
12		207	3.7243.612	3.7244.612

		TiA	Co			VA			
3.5200.xxx		✓		✓	✓		✓		3.5
3.5205.xxx		✓		✓			✓		3.5
3.5231.xxx		✓		✓	✓			✓	3.9
3.5225.xxx		✓		✓	✓		✓		2.4
3.1306.xxx		✓		✓			✓		3.5
4.5236.xxx			✓	✓	✓	✓	✓		3.5
3.5207.006		✓		✓	✓				3.5
3.1221.050		✓		✓	✓				



Palette für 5,0ChLP Platten 3.7243/3.7244 4x4 1/2H
Tray for plates 5.0ChLP 3.7243/3.7244 4x4 1/2H
Paleta para placas 5,0ChLP 3.7243/3.7244 4x4 1/2H

14.0205.429




Platte 3.7243.608 Probeimplantat
Plate 3.7243.608 trial
Plantilla para placa 3.7243.608
43.7243.608
Platte 3.7244.608 Probeimplantat
Plate 3.7244.608 trial
Plantilla para placa 3.7244.608
43.7244.608



		Ti	
	Len	L	R
4	91	3.7037.604	3.7036.604
6	111	3.7037.606	3.7036.606
8	131	3.7037.608	3.7036.608
10	151	3.7037.610	3.7036.610
12	171	3.7037.612	3.7036.612

									
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				



Palette für 5,0ChLP Platten 4x4 1/2H
Tray for 5.0ChLP plates 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419
14.0205.420

Zielblock R (3.7036)
Aiming block R [3.7036]
Bloque guía R (3.7036)
Zielblock L (3.7037)
Aiming block L [3.7037]
Bloque guía L (3.7037)

40.8210.000

40.8211.000

Platte 3.7036.606 Probeimplantat
Plate 3.7036.606 trial
Plantilla para placa 3.7036.606
Platte 3.7037.606 Probeimplantat
Plate 3.7037.606 trial
Plantilla para placa 3.7037.606

43.7036.606

43.7037.606



			Ti	
	Len	L	R	
4	75	3.7039.604	3.7038.604	
6	94	3.7039.606	3.7038.606	
8	113	3.7039.608	3.7038.608	
10	131	3.7039.610	3.7038.610	
12	150	3.7039.612	3.7038.612	


		TiA	Co			VA			
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	✓	3.5

mindestens 2 Löcher
min 2 holes
min. 2 orificios



90°

3.7071 / 3.7072
+
3.7038 / 3.7039



Palette für 5,0ChLP Platten 4x4 1/2H
Tray for 5.0ChLP plates 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419
14.0205.420

Zielblock R (3.7038)
Aiming block R [3.7038]
Bloque guía R (3.7038)
Zielblock L (3.7039)
Aiming block L [3.7039]
Bloque guía L (3.7039)

40.8212.000

40.8213.000

Platte 3.7038.606 Probeimplantat
Plate 3.7038.606 trial
Plantilla para placa 3.7038.606
Platte 3.7039.606 Probeimplantat
Plate 3.7039.606 trial
Plantilla para placa 3.7039.606


43.7038.606

43.7039.606



			Ti
	Len	L	R
4	84	3.7041.604	3.7040.604
6	104	3.7041.606	3.7040.606
8	124	3.7041.608	3.7040.608
10	144	3.7041.610	3.7040.610
12	164	3.7041.612	3.7040.612

									
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				



Paleta na płytki 5,0ChLP 4x4 1/2H
Tray for plates 5.0ChLP 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419
14.0205.420

Zielblock R (3.7040)
Aiming block R [3.7040]
Bloque guía R (3.7040)
Zielblock L (3.7041)
Aiming block L [3.7041]
Bloque guía L (3.7041)

40.8214.000

40.8215.000

Platte 3.7040.606 Probeimplantat
Plate 3.7040.606 trial
Plantilla para placa 3.7040.606
Platte 3.7041.606 Probeimplantat
Plate 3.7041.606 trial
Plantilla para placa 3.7041.606


43.7040.606

43.7041.606



			Ti	
	Len	L	R	
4	91	3.7071.604	3.7072.604	
6	111	3.7071.606	3.7072.606	
8	131	3.7071.608	3.7072.608	
10	151	3.7071.610	3.7072.610	
12	171	3.7071.612	3.7072.612	

		TiA	Co		VA			
3.5200.xxx	✓	✓	✓	✓	✓	3.5		
3.5205.xxx	✓	✓	✓	✓	✓	3.5		
3.5231.xxx	✓	✓	✓	✓	✓	3.9		
3.5225.xxx	✓	✓	✓	✓	✓	2.4		
3.1306.xxx	✓	✓	✓	✓	✓	3.5		
4.5236.xxx	✓	✓	✓	✓	✓	3.5		
3.5207.006	✓	✓	✓	✓	✓	3.5		
3.1221.050	✓	✓	✓	✓	✓	3.5		



Palette für 5,0ChLP Platten 4x4 1/2H
Tray for 5.0ChLP plates 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419
14.0205.420

Zielblock L (3.7071)
Aiming block L [3.7071]
Bloque guía L (3.7071)
Zielblock R (3.7072)
Aiming block R [3.7072]
Bloque guía R (3.7072)

40.8217.000

40.8216.000

Platte 3.7071.606 Probeimplantat
Plate 3.7071.606 trial
Plantilla para placa 3.7071.606
Platte 3.7072.606 Probeimplantat
Plate 3.7072.606 trial
Plantilla para placa 3.7072.606


43.7071.606

43.7072.606




			Ti
	Len	L	R
4	75	3.7073.604	3.7074.604
6	94	3.7073.606	3.7074.606
8	113	3.7073.608	3.7074.608
10	132	3.7073.610	3.7074.610
12	151	3.7073.612	3.7074.612

		TiA	Co			VA		
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	3.5



90°
3.7071 / 3.7072
+
3.7073 / 3.7074



Palette für 5,0ChLP Platten 4x4 1/2H
Tray for 5.0ChLP plates 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419
14.0205.420

Zielblock L (3.7073)
Aiming block L [3.7073]
Bloque guía L (3.7073)
Zielblock R (3.7074)
Aiming block R [3.7074]
Bloque guía R (3.7074)

40.8219.000

40.8218.000

Platte 3.7073.606 Probeimplantat
Plate 3.7073.606 trial
Plantilla para placa 3.7073.606
Platte 3.7074.606 Probeimplantat
Plate 3.7074.606 trial
Plantilla para placa 3.7074.606

43.7073.606

43.7074.606



6	82	3.7240.606
8	102	3.7240.608

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				







3.4097.505

			Ti
	Len	L	R
2	88	3.4097.502	3.4098.502
4	121	3.4097.504	3.4098.504
6	151	3.4097.506	3.4098.506
8	181	3.4097.508	3.4098.508
2-12	88 240		

		TiA	Co		VA			
	3.5200.xxx	✓		✓	✓		✓	3.5
	3.5205.xxx	✓		✓			✓	3.5
	3.5231.xxx	✓		✓	✓		✓	3.9
	3.5225.xxx	✓		✓	✓		✓	2.4
	3.1306.xxx	✓		✓			✓	3.5
	4.5236.xxx		✓	✓	✓	✓	✓	3.5
	3.5207.006	✓		✓	✓			3.5
	3.1221.050	✓		✓	✓			



Palette für 5,0ChLP Platten 3.4097/3.4098 4x4 1/2H
Tray for 5.0ChLP plates 3.4097/3.4098 4x4 1/2H
Paleta para placas 5,0ChLP 3.4097/3.4098 4x4 1/2H

14.0205.409

Zielblock L (3.4097)
Aiming block L [3.4097]
Bloque guía L (3.4097)
Zielblock R (3.4098)
Aiming block R [3.4098]
Bloque guía R (3.4098)


40.5731.100

40.5731.200

Platte 3.4097.504 Probeimplantat
Plate 3.4097.504 trial
Plantilla para placa 3.4097.504
Platte 3.4098.504 Probeimplantat
Plate 3.4098.504 trial
Plantilla para placa 3.4098.504


43.4097.504

43.4098.504



	Ti	
1 Len		GERADE STRAIGHT RECTA
130	3.7033.004	
2 Len		KURZE BIEGUNG SHORT BEND CURVATURA PEQUEÑA
131	3.7034.004	
3 Len		STANDARD BIEGUNG STANDARD BEND CURVATURA ESTÁNDAR
140	3.7035.004	

		TiA	Co			VA			
3.5200.xxx	✓		✓	✓		✓			3.5
3.5205.xxx	✓		✓			✓			3.5
3.5231.xxx	✓		✓	✓			✓		3.9
3.5225.xxx	✓		✓	✓		✓			2.4
3.1306.xxx	✓		✓			✓			3.5
3.1220.xxx	✓		✓			✓			2.7
4.5236.xxx		✓	✓	✓	✓	✓			3.5
3.5207.006	✓		✓	✓					3.5
3.1221.050				✓					




Palette für 5,0ChLP Platten 3.7033-3.7035 4x4 1/2H
Tray for 5,0ChLP plates 3.7033-3.7035 4x4 1/2H
Paleta para placas 5,0ChLP 3.7033-3.7035 4x4 1/2H

14.0205.410





Platte 3.7033.004 Probeimplantat Plate 3.7033.004 trial Plantilla para placa 3.7033.004	43.7033.004
Platte 3.7034.004 Probeimplantat Plate 3.7034.004 trial Plantilla para placa 3.7034.004	43.7034.004
Platte 3.7035.004 Probeimplantat Plate 3.7035.004 trial Plantilla para placa 3.7035.004	43.7035.004




		Ti		
	Len	L	R	
5	73	3.7067.605	3.7068.605	
7	97	3.7067.607	3.7068.607	
9	122	3.7067.609	3.7068.609	
11	148	3.7067.611	3.7068.611	
5-15	73-194			

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				




Palette für 5,0ChLP Platten 3.7067/3.7068 4x4 1/2H
Tray for 5.0ChLP plates 3.7067/3.7068 4x4 1/2H
Paleta para placas 5,0ChLP 3.7067/3.7068 4x4 1/2H

14.0205.411



Zielblock L (3.7067)	40.5742.100
Aiming block L [3.7067]	
Bloque guía L (3.7067)	
Zielblock R (3.7068)	40.5742.200
Aiming block R [3.7068]	
Bloque guía R (3.7068)	





3.7091.606

Ti

	Len		
4	110	3.7091.604	3.7092.604
5	125	3.7091.605	3.7092.605
6	140	3.7091.606	3.7092.606
8	170	3.7091.608	3.7092.608
3-15	95-275		

		TiA	Co			VA		
	3.5200.xxx	✓		✓	✓		✓	3.5
	3.5205.xxx	✓		✓			✓	3.5
	3.5231.xxx	✓		✓	✓			3.9
	3.5225.xxx	✓		✓	✓		✓	2.4
	3.1306.xxx	✓		✓			✓	3.5
	4.5236.xxx		✓	✓	✓	✓	✓	3.5
	3.5207.006	✓		✓	✓			3.5
	3.1221.050	✓		✓	✓			



Palette für 5,0ChLP Platten 3.7091/3.7092 4x4 1/2H
Tray for 5.0ChLP plates 3.7091/3.7092 4x4 1/2H
Paleta para placas 5,0ChLP 3.7091/3.7092 4x4 1/2H

14.0205.421

Zielblock L [3.7091]
Aiming block L [3.7091]
Bloque guía L [3.7091]
Zielblock R [3.7092]
Aiming block R [3.7092]
Bloque guía R [3.7092]

40.8223.000

40.8222.000




3.7203.608

Ti

	Len	L	R
4	121	3.7203.604	3.7202.604
5	134	3.7203.605	3.7202.605
6	147	3.7203.606	3.7202.606
8	173	3.7203.608	3.7202.608

4-16 121-277

	TiA	Co		VA		
3.5200.xxx	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	3.5



Palette für 5,0ChLP Platten 3.7202-3.7203 4x4 1/2H
Tray for plates 5.0ChLP 3.7202/3.7203 4x4 1/2H
Paleta para placas 5,0ChLP 3.7202/3.7203 4x4 1/2H

14.0205.424



Zielblock R [3.7202]
Aiming block R [3.7202]
Bloque guía R [3.7202]

40.8226.000

Zielblock L [3.7203]
Aiming block L [3.7203]
Bloque guía L [3.7203]

40.8227.000



Platte 3.7202.606 Probeimplantat
Plate 3.7202.606 trial
Plantilla para placa 3.7202.606

43.7202.606

Platte 3.7203.606 Probeimplantat
Plate 3.7203.606 trial
Plantilla para placa 3.7203.606


43.7203.606



3.4011.609

			Ti
	Len	L	R
3	100	3.4011.603	3.4012.603
5	126	3.4011.605	3.4012.605
7	152	3.4011.607	3.4012.607
9	178	3.4011.609	3.4012.609
1-17	74-282		

		TiA	Co			VA		
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	



Palette für 5,0ChLP Platten 3.4011/4012.6xx 4x4 1/2H
Tray for plates 5.0ChLP 3.4011/4012.6xx 4x4 1/2H
Paleta para placas 5,0ChLP 3.4011/4012.6xx 4x4 1/2H

14.0205.427



Zielblock R (3.4012.6xx)
Aiming block R (3.4012.6xx)
Bloque guía R (3.4012.6xx)
Zielblock L (3.4011.6xx)
Aiming block L (3.4011.6xx)
Bloque guía L (3.4011.6xx)

40.8232.000

40.8233.000



Platte 3.4012.605 Probeimplantat
Plate 3.4012.605 trial
Plantilla para placa 3.4012.605

Platte 3.4011.605 Probeimplantat
Plate 3.4011.605 trial
Plantilla para placa 3.4011.605

43.4012.605

43.4011.605





3.7094.604



	Len	
4	66	3.7094.604
6	86	3.7094.606



	3.5200.xxx	✓	✓	✓	✓	3.5
	3.5205.xxx	✓	✓		✓	3.5
	3.5231.xxx	✓	✓	✓		3.9
	3.5225.xxx	✓	✓	✓	✓	2.4
	3.1306.xxx	✓	✓		✓	3.5
	4.5236.xxx		✓	✓	✓	3.5
	3.5207.006	✓	✓	✓		3.5
	3.1221.050	✓	✓	✓		



Palette für 5,0ChLP Platten 3.7094/3.7095 4x4 1/2H
Tray for 5.0ChLP plates 3.7094/3.7095 4x4 1/2H
Paleta para placas 5,0ChLP 3.7094/3.7095 4x4 1/2H

14.0205.422





Zielblock [3.7094]
Aiming block [3.7094]
Bloque guía [3.7094]

40.8224.000



4	66	3.7095.604
6	86	3.7095.606

	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	3.1306.xxx	✓		✓			✓		3.5
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5
	3.1221.050	✓		✓	✓				



Palette für 5,0ChLP Platten 3.7094/3.7095 4x4 1/2H
Tray for 5.0ChLP plates 3.7094/3.7095 4x4 1/2H
Paleta para placas 5,0ChLP 3.7094/3.7095 4x4 1/2H

14.0205.422



Zielblock [3.7095]
Aiming block [3.7095]
Bloque guía [3.7095]


40.8225.000



		Len	Ti
4	56	3.7238.604	
5	66	3.7238.605	



		TiA	Co			VA		
	3.5200.xxx	✓		✓	✓		✓	3.5
	3.5205.xxx	✓		✓			✓	3.5
	3.5231.xxx	✓		✓	✓			3.9
	3.5225.xxx	✓		✓	✓		✓	2.4
	3.1306.xxx	✓		✓			✓	3.5
	4.5236.xxx		✓	✓	✓	✓	✓	3.5
	3.5207.006	✓		✓	✓			3.5
	3.1221.050	✓		✓	✓			





	Len	Ti
7	167	3.4033.507
9	197	3.4033.509

		TiA	Co			VA		
	3.5200.xxx	✓		✓	✓		✓	3.5
	3.5205.xxx	✓		✓			✓	3.5
	3.5231.xxx	✓		✓	✓		✓	3.9
	3.5225.xxx	✓		✓	✓		✓	2.4
	3.1306.xxx	✓		✓			✓	3.5
	4.5236.xxx		✓	✓	✓	✓	✓	3.5
	3.5207.006	✓		✓	✓			3.5
	3.1221.050	✓		✓	✓			





Palette für 5,0ChLP Platten 3.4033 4x4 1/2H
Tray for 5.0ChLP plates 3.4033 4x4 1/2H
Paleta para placas 5,0ChLP 3.4033 4x4 1/2H

14.0205.414



			Ti	
	Len	L		R
4	108	3.4039.604		3.4040.604
6	134	3.4039.606		3.4040.606
8	160	3.4039.608		3.4040.608
4-16	108 264			

		TiA	Co			VA		
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	3.5



Palette für 5,0ChLP Platten 3.4039/3.4040 4x4 1/2H
Tray for 5.0ChLP plates 3.4039/3.4040 4x4 1/2H
Paleta para placas 5,0ChLP 3.4039/3.4040 4x4 1/2H

14.0205.415

Zielblock [3.4039]
Aiming block [3.4039]
Bloque guía [3.4039]

Zielblock [3.4040]
Aiming block [3.4040]
Bloque guía [3.4040]

40.5726.100

Platte 3.4039.608 Probeimplantat
Plate 3.4039.608 trial
Plantilla para placa 3.4039.608

43.4039.608

40.5726.200


Platte 3.4040.608 Probeimplantat
Plate 3.4040.608 trial
Plantilla para placa 3.4040.608

43.4040.608



			Ti	
	Len	L		R
5	114	3.4051.605		3.4052.605
7	140	3.4051.607		3.4052.607
9	166	3.4051.609		3.4052.609
3-17	88 270			

		TiA	Co			VA		
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	3.5



Palette für 5,0ChLP Platten 3.4051/3.4052 4x4 1/2H
Tray for plates 5.0ChLP 3.4051/3.4052 4x4 1/2H
Paleta para placas 5,0ChLP 3.4051/3.4052 4x4 1/2H

14.0205.428



Zielblock L (3.4051.6xx)
Aiming block L [3.4051.6xx]
Bloque guía L (3.4051.6xx)

40.8235.000

Zielblock R (3.4052.6xx)
Aiming block R [3.4052.6xx]
Bloque guía R (3.4052.6xx)

40.8234.000



Platte 3.4051.607 Probeimplantat
Plate 3.4051.607 trial
Plantilla para placa 3.4051.607

43.4051.607



Platte 3.4052.607 Probeimplantat
Plate 3.4052.607 trial
Plantilla para placa 3.4052.607

43.4052.607



			Ti
	Len	L	R
4	85	3.7029.604	3.7030.604
5	95	3.7029.605	3.7030.605
6	105	3.7029.606	3.7030.606
8	125	3.7029.608	3.7030.608
10	145	3.7029.610	3.7030.610
4-12	85-165		

		TiA	Co			VA		
3.5200.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5205.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5231.xxx	✓	✓	✓	✓	✓	✓	✓	3.9
3.5225.xxx	✓	✓	✓	✓	✓	✓	✓	2.4
3.1306.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
4.5236.xxx	✓	✓	✓	✓	✓	✓	✓	3.5
3.5207.006	✓	✓	✓	✓	✓	✓	✓	3.5
3.1221.050	✓	✓	✓	✓	✓	✓	✓	3.5



Palette für 5,0ChLP Platten 3.7029/3.7030 4x4 1/2H
Tray for 5.0ChLP plates 3.7029/3.7030 4x4 1/2H
Paleta para placas 5,0ChLP 3.7029/3.7030 4x4 1/2H

14.0205.417




Platte 3.7029.605 Probeimplantat
Plate 3.7029.605 trial
Plantilla para placa 3.7029.605
Platte 3.7030.605 Probeimplantat
Plate 3.7030.605 trial
Plantilla para placa 3.7030.605




	Len	L	R
Klein Small Pequeña	59	3.7201.601	3.7200.601
Mittel Medium Media	63	3.7201.602	3.7200.602
Gross Large Grande	66	3.7201.603	3.7200.603

	TiA	Co			VA			
3.5200.xxx	✓		✓	✓		✓		3.5
3.5205.xxx	✓		✓			✓		3.5
3.5231.xxx	✓		✓	✓			✓	3.9
3.5225.xxx	✓		✓	✓		✓		2.4
4.5236.xxx		✓	✓	✓	✓	✓		3.5
3.5207.006	✓		✓	✓				3.5



Palette für 5,0ChLP Platten 3.7200/3.7201 4x4 1/2H
Tray for plates 5.0ChLP 3.7200/3.7201 4x4 1/2H
Paleta para placas 5,0ChLP 3.7200/3.7201 4x4 1/2H


14.0205.423



Ti	
L	R
3.4057.000	3.4058.000

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5







	Ti	
L		R
3.4059.000		3.4060.000

		TiA	Co			VA			
	3.5200.xxx	✓		✓	✓		✓		3.5
	3.5205.xxx	✓		✓			✓		3.5
	3.5231.xxx	✓		✓	✓			✓	3.9
	3.5225.xxx	✓		✓	✓		✓		2.4
	4.5236.xxx		✓	✓	✓	✓	✓		3.5
	3.5207.006	✓		✓	✓				3.5



VORSCHLAG ZUR ZUSAMMENSETZUNG DES 5,0 ChLP MODULAREN INSTRUMENTARIUMS
5.0ChLP MODULAR SET COMPLETION PROPOSAL
EJEMPLO DE COMPOSICIÓN EL SET MODULAR 5,0ChLP



VARIANTE I VARIANT I VARIANTE I

Stk. Pcs Uds.	Name Name Nombre	Kat.-Nr Catalogue no. Nº de catálogo
1	5,0ChLP Behälter 9x4H 5.0ChLP container 9x4 H Contenedor 5,0ChLP 9x4H	14.0205.101
2	5,0ChLP Behälterdeckel 9x4H 5.0ChLP container lid 9x4 H Tapa de contenedor 5,0ChLP 9x4H	14.0205.102
3	Instrumentarium für 5,0ChLP 4x1 1/2H Instrument set for 5.0ChLP 4x4 1/2H Set de instrumental para 5,0ChLP 4x1 1/2H	15.0205.206
4	Instrumentarium für 5,0ChLP 4x4 1/2H Instrument set for 5.0ChLP 4x4 1/2H Set de instrumental para 5,0ChLP 4x4 1/2H	15.0205.202
5	Behälter für Schrauben 5,0ChLP 4x2H Stand for 5.0ChLP screws 4x2 H Contenedor para tornillos 5,0ChLP 4x2H	14.0205.501
6	Behälter für Schrauben 5,0ChLP 4x2H Stand for 5.0ChLP screws 4x2 H Contenedor para tornillos 5,0ChLP 4x2H	14.0205.502
7	Palette 4x1 1/2H - x2szt Tray 4x1 1/2H - x2pcs Paleta 4x1 1/2H - x2Uds.	14.0000.201

VARIANTE II - ERWEITERT (KANÜLIERTE SCHRAUBEN) VARIANT II - EXTENDED (CANNULATED SCREWS) VARIANTE II - EXTENDIDO (TORNILLOS CANULADOS)

Stk. Pcs Uds.	Name Name Nombre	Kat.-Nr Catalogue no. Nº de catálogo
1	5,0ChLP Behälter 9x4H 5.0ChLP container 9x4 H Contenedor 5,0ChLP 9x4H	14.0205.101
2	5,0ChLP Behälterdeckel 9x4H 5.0ChLP container lid 9x4 H Tapa de contenedor 5,0ChLP 9x4H	14.0205.102
3	Instrumentarium für 5,0ChLP 4x1 1/2H Instrument set for 5.0ChLP 4x4 1/2H Set de instrumental para 5,0ChLP 4x1 1/2H	15.0205.206
4	Instrumentarium für 5,0ChLP 4x4 1/2H Instrument set for 5.0ChLP 4x4 1/2H Set de instrumental para 5,0ChLP 4x4 1/2H	15.0205.202
5	Instrumentarium für 5,0ChLP 4x1 1/2H Instrument set for 5.0ChLP 4x1 1/2H Set de instrumental para 5,0ChLP 4x1 1/2H	15.0205.203
6	Behälter für Schrauben 5,0ChLP 4x2H Stand for 5.0ChLP screws 4x2 H Contenedor para tornillos 5,0ChLP 4x2H	14.0205.503
7	Behälter für Schrauben 5,0ChLP 4x2H Stand for 5.0ChLP screws 4x2 H Contenedor para tornillos 5,0ChLP 4x2H	14.0205.504
8	Palette 4x1 1/2H Tray 4x1 1/2H Paleta 4x1 1/2H	14.0000.201



Stk. Pcs Uds.	Name Name Nombre	Kat.-Nr Catalogue no. Nº de catálogo	Stk. Pcs Uds.
1	5,0ChLP Behälterdeckel 9x4H 5.0ChLP container lid 9x4H Tapa de contenedor 5,0ChLP 9x4H	14.0205.102	1
2	5,0ChLP Behälter 9x4H 5.0ChLP container 9x4H Contenedor 5,0ChLP 9x4H	14.0205.101	2

INSTRUMENTARIUM FÜR 5,0ChLP 4X1 1/2H
INSTRUMENT SET FOR 5.0ChLP 4X4 1/2H
SET DE INSTRUMENTAL PARA 5,0ChLP 4X1 1/2H


15.0205.206	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 5,0ChLP 4x1 1/2H Tray for 5.0ChLP instrument set 4x4 1/2H Paleta para instrumental 5,0ChLP 4x1 1/2H	14.0205.206	1
	Kirschnerdraht 1,5/210 Kirschner wire 1.5/210 Aguja KIRSCHNER 1,5/210	40.4592.210	4
	Bohrer 1,8/210 Drill 1.8/210 Broca 1,8/210	40.2063.212	2
	Skalierter Bohrer 2,5/210 Drill with scale 2.5/210 Broca graduada 2,5/210	40.5912.212	2
	Skalierter Bohrer 2,8/210 Drill with scale 2.8/210 Broca graduada 2,8/210	40.5653.212	2
	Schraubendrehereinsatz T15 Screwdriver tip T15 Punta de destornillador T15	40.5677.000	1
	Drehmomentgriff mit Kupplung 2Nm Torque limiting ratchet handle 2Nm Mango dinamométrico con trinquete 2Nm	40.6652.000	1
	Knebelgriff mit Kupplung Handle ratchet device Mango con trinquete	40.6654.000	1
	Schutzhülse 7/5 Protective guide 7/5 Guía de protección 7/5	40.5672.000	2
	Führung VA 2,8 Guide VA 2.8 Guía VA 2,8	40.8206.028	1
	Kompressionsführung 2,5 Compression guide 2.5 Guía de compresión 2,5	40.4804.725	1
	Führungshülse 5,0/1,8 Guide sleeve 5.0/1.8 Guía de la broca 5,0/1,8	40.5673.718	2
	Führungshülse 5,0/2,8 Guide sleeve 5.0/2.8 Guía de la broca 5,0/2,8	40.5673.728	4
	Tiefenmesslehre Depth measure Medidor de profundidad	40.4639.550	1

INSTRUMENTARIUM 5,0ChLP 4x4 1/2H
INSTRUMENT SET FOR 5.0ChLP 4x4 1/2H
SET DE INSTRUMENTAL PARA 5.0ChLP 4x4 1/2H



15.0205.202	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 5,0ChLP 4x4 1/2H Tray for 5.0ChLP instrument set 4x4 1/2H Paleta para instrumental 5,0ChLP 4x4 1/2H	14.0205.202	1
	Plattenrepositionsschraube 2,8/180 Setting-compressing screw 2.8/180 Tornillo de ajuste-compresión 2,8/180	40.5674.728	1
	Messlehre für Schraubenlänge Screw length measure Medidor de tornillos	40.5675.500	1
	Biegegerät für Platten 5,0 Plates bender 5.0 Grifa para doblar placas 5,0	40.4643.500	2
	Dreiarmer Schraubendrehereinsatz 5,0ChLP Tripod screwdriver tip 5.0ChLP Punta de destornillador de tres brazos 5,0ChLP	40.6271.500	1
	Einsatz T15 mit Greifer T15 screwdriver tip with holder Punta de destornillador T15 con sujeción	40.6254.000	1
	HA Gewindebohrer 3,5 mit Handgriff Cortical tap HA 3.5 with handle Macho de roscar cortical HA 3,5 con mango	40.2548.200	1
	Gewindebohrer 5,0ChLP-3,5 Tap 5.0ChLP-3.5 Macho de roscar 5,0ChLP - 3,5	40.5661.000	1
Optionales Instrument Optional instrument Instrumento opcional			
	Dynamometrischer Verbinder 2Nm Torque connector 2Nm Conector dinamométrico 2Nm	40.5927.020	1

INSTRUMENTARIUM FÜR 5,0ChLP 4x4 1/2H
INSTRUMENT SET FOR 5.0ChLP 4x4 1/2H
SET DE INSTRUMENTAL PARA PLACAS 5,0ChLP 4x4 1/2H

15.0205.208	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 5,0ChLP 4x4 1/2H Tray for 5.0ChLP instrument set 4x4 1/2H Paleta para instrumental 5,0ChLP 4x4 1/2H	14.0205.208	1
	Kirschnerdraht 1,5/210 Kirschner wire 1.5/210 Aguja de Kirschner 1,5/210	40.4592.210	4
	Bohrer 1,8/210 Drill 1.8/210 Broca 1,8/210	40.2063.212	2
	Skalierter Bohrer 2,8/210 Drill with scale 2.8/210 Broca graduada 2,8/210	40.5653.212	2
	Plattenrepositionsschraube 2,8/180 Setting-compressing screw 2.8/180 Tornillo de ajuste-compresión 2,8/180	40.5674.728	1
	Drehmomentgriff mit Kupplung 2Nm Torque limiting ratchet handle 2Nm Mango dinamométrico de carraca 2Nm	40.6652.000	1
	Führung VA 2,8 Guide VA 2.8 Guía VA 2,8	40.8206.028	1
	Führungshülse 5,0/2,8 Guide sleeve 5.0/2.8 Guía de la broca 5,0/2,8	40.5673.728	3
	Führungshülse 5,0/1,8 Guide sleeve 5.0/1.8 Guía de la broca 5,0/1,8	40.5673.718	2
	Schutzhülse 7/5 Protective guide 7/5 Guía de protección 7/5	40.5672.000	2
	Schraubendrehereinsatz T15 Screwdriver tip T15 Punta de destornillador T15	40.5677.000	1
	Medidor de profundidad Depth measure Tiefenmesslehre	40.4639.550	1
Optionales Instrument Optional instrument Instrumento opcional			
	Dynamometrischer Verbinder 2Nm Torque connector 2Nm Conector dinamométrico 2Nm	40.5927.020	0

INSTRUMENTARIUM FÜR 5,0ChLP 4X1 1/2H
INSTRUMENT SET FOR 5.0ChLP 4X1 1/2H
SET DE INSTRUMENTAL PARA 5,0ChLP 4X1 1/2H

15.0205.203


Erweiterung auf 5,0ChLP kanülierte Schrauben Extension for cannulated screws 5.0ChLP Extensión para tornillos canulados de 5.0ChLP	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
--	------------------------	---------------------	---





Palette für Instrumentarium 5,0ChLP 4x1 1/2H
Tray for 5.0ChLP instrument set 4x1 1/2H
Paleta para instrumental 5,0ChLP 4x1 1/2H

1 14.0205.203

	Skalierter, Kanülierter Bohrer 3,2/1,5/210 Cannulated drill with scale 3.2/1.5/210 Broca canulada graduada 3,2/1,5/210	1	40.5990.212
	Kanülierter Schraubendreher T15 Cannulated screwdriver tip T15 Punta de destornillador canulada T15	1	40.5992.000
	Führungshülse 3,2/1,5 Guide sleeve 3.2/1.5 Guía de la broca 3,2/1,5	1	40.5994.715
	Führungshülse 5,0/3,2 Guide sleeve 5.0/3.2 Guía de la broca 5,0/3,2	1	40.5994.732




BEHÄLTER FÜR SCHRAUBEN 5,0ChLP 4X2H STAND FOR 5.0ChLP SCREWS 4X2 H BANDEJA PARA TORNILLOS 5,0ChLP 4X2H		14.0205.501																	
	5,0ChLP Selbstschneidende Schraube 3,5 5.0ChLP self-tapping screw 3.5 5,0ChLP tornillo auto-roscante 3,5																		
																		L	
	L[mm]	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40			
	Stk. Pcs Uds.	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4		
	L[mm]	42	44	46	48	50	52	54	56	58	60	65	70	75	80	85			
Stk. Pcs Uds.	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4			
5,0ChLP Schraube für Cerclage 5.0ChLP cerclage screw 5,0ChLP tortillo para cerclaje																			
Stk. Pcs Uds.	3																		
* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes																			

* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

BEHÄLTER FÜR SCHRAUBEN 5,0CHLP 4X2H STAND FOR 5.0CHLP SCREWS 4X2 H BANDEJA PARA TORNILLOS 5,0CHLP 4X2H	14.0205.506
	5,0ChLP Selbstschneidende Schraube 3,5 5.0ChLP self-tapping screw 3.5 5,0ChLP tornillo auto-roscante 3,5
	
	L[mm]121416182022242628303234363840
	Stk. Pcs Uds.4444444444444444
	L[mm]424446485052545658606570758085
Stk. Pcs Uds.4444444444444444	
5,0ChLP VA Schraube 3,5 5.0ChLP screw VA 3.5 5,0ChLP tornillo VA 3,5	
	
L[mm]121416182022242628303234363840	
Stk. Pcs Uds.3333333333333333	
L[mm]424446485052545658606570758085	
Stk. Pcs Uds.3333333333333333	

* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

BEHÄLTER FÜR SCHRAUBEN 5,0CHLP 4X2H STAND FOR 5.0CHLP SCREWS 4X2 H BANDEJA PARA TORNILLOS 5,0CHLP 4X2H	14.0205.502																
	Selbstschneidende Kortikalisschraube 3,5 Cortical self-tapping screw 3.5 Tornillo cortical auto-roscante 3,5																
	L[mm]	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40	
	Stk. Pcs Uds.	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	
	L[mm]	45	50	55	60	65	70	75	80	85							
	Stk. Pcs Uds.	2	2	2	2	2	2	2	2	2							
	4,5ChLP Schraube 2,4 4.5ChLP screw 2.4 4,5ChLP tornillo 2,4																
	L[mm]	10	12	14	16	18	20	22	24	26	28	30	32	34	36	38	
	Stk. Pcs Uds.	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3	
	L[mm]	40															
	Stk. Pcs Uds.	3															
5,0ChLP Spongiaschraube 3,9 5.0ChLP cancellous screw 3.9 5,0ChLP tornillo esponjosa 3,9																	
L[mm]	35	40	45	50	55	60	65	70	75	80							
Stk. Pcs Uds.	2	2	2	2	2	2	2	2	2	2							
* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes																	


* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

BEHÄLTER FÜR SCHRAUBEN 5,0ChLP 4X2H
STAND FOR 5.0ChLP SCREWS 4X2 H
BANDEJA PARA TORNILLOS 5,0ChLP 4X2H

14.0205.503




5,0ChLP Selbstschneidende Schraube 3,5
5.0ChLP self-tapping screw 3.5
5,0ChLP tornillo auto-roscante 3,5




L[mm]	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40
Stk. Pcs Uds.	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4
L[mm]	42	44	46	48											
Stk. Pcs Uds.	4	4	4	4											

Selbstschneidende Kortikalisschraube 3,5
Cortical self-tapping screw 3.5
Tornillo cortical auto-roscante 3,5




L[mm]	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40
Stk. Pcs Uds.	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
L[mm]	45														
Stk. Pcs Uds.	2														

4,5ChLP Schraube 2,4
4.5ChLP screw 2.4
4,5ChLP tornillo 2,4




L[mm]	12	14	16	18	20	22	24	26	28	30	32	34	36	38	40
Stk. Pcs Uds.	3	3	3	3	3	3	3	3	3	3	3	3	3	3	3

5,0ChLP Spongiaschraube 3,9
5.0ChLP cancellous screw 3.9
5,0ChLP tornillo esponjosa 3,9




L[mm]	35	40	45												
Stk. Pcs Uds.	3	3	3												

5,0ChLP Kanülierte Schraube 3,9
5.0ChLP cannulated screw 3.9
5,0ChLP tornillo canulado 3,9




L[mm]	35	40	45												
Stk. Pcs Uds.	3	3	3												

5,0ChLP Kanülierte Kortikalisschraube 3,9
5.0ChLP cortical cannulated screw 3.9x16
5,0ChLP tornillo cortical canulado 3,9x16



L[mm]	35	40	45												
Stk. Pcs Uds.	3	3	3												

5,0ChLP Konische kanülierte Schraube 3,9x16
5.0ChLP conical cannulated screw 3.9x16
5,0ChLP tornillo cónico canulado 3,9x16



L[mm]	35	40	45												
Stk. Pcs Uds.	3	3	3												

5,0ChLP Schraube für Cerclage
5.0ChLP cerclage screw
5,0ChLP tornillo para cerclaje



Stk. Pcs Uds.	3														
---------------	---	--	--	--	--	--	--	--	--	--	--	--	--	--	--

5,0ChLP Distanzschraube
5.0ChLP distance screw
5,0ChLP tornillo de distancia



Stk. Pcs Uds.	6														
---------------	---	--	--	--	--	--	--	--	--	--	--	--	--	--	--


* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

BEHÄLTER FÜR SCHRAUBEN 5,0ChLP 4X2H
STAND FOR 5.0ChLP SCREWS 4X2 H
BANDEJA PARA TORNILLOS 5,0ChLP 4X2H

14.0205.504




5,0ChLP Selbstschneidende Schraube 3,5
5.0ChLP self-tapping screw 3.5
5,0ChLP tornillo auto-roscante 3,5




L[mm]	50	52	54	56	58	60	65	70	75	80	85	
Stk. Pcs Uds.	4	4	4	4	4	4	4	4	4	4	4	

Selbstschneidende Kortikalisschraube 3,5
Cortical self-tapping screw 3.5
Tornillo cortical auto-roscante 3,5




L[mm]	50	55	60	65	70	75	80	85	
Stk. Pcs Uds.	2	2	2	2	2	2	2	2	

5,0ChLP Spongiaschraube 3,9
5.0ChLP cancellous screw 3.9
5,0ChLP tornillo esponjosa 3,9




L[mm]	50	55	60	65	70	75	80	
Stk. Pcs Uds.	4	4	4	4	4	4	4	

5,0ChLP Kanülierte Schraube 3,9
5.0ChLP cannulated screw 3.9
5,0ChLP tornillo canulado 3,9




L[mm]	50	55	60	65	70	75	80	
Stk. Pcs Uds.	4	4	4	4	4	4	4	

5,0ChLP Kanülierte Kortikalisschraube 3,9
5.0ChLP cortical cannulated screw 3.9x16
5,0ChLP tornillo cortical canulado 3,9x16



L[mm]	50	55	60	65	70	75	80	
Stk. Pcs Uds.	4	4	4	4	4	4	4	

5,0ChLP Konische kanülierte Schraube 3,9
5.0ChLP cortical cannulated screw 3.9x16
5,0ChLP tornillo cortical canulado 3,9x16



L[mm]	50	55	60	65	70	75	80	
Stk. Pcs Uds.	4	4	4	4	4	4	4	

* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

	3.5200.xxx	12÷95	✓		✓		✓		✓		3.5
	4.5236.xxx	12÷95		✓	✓			✓	✓		3.5
	3.5225.xxx	16÷40	✓		✓		✓		✓		2.4
	3.1306.xxx	10÷110	✓		✓				✓		3.5
	3.1220.xxx	6÷65	✓		✓				✓		2.7
	3.5205.xxx	18÷90	✓		✓				✓		3.5
	3.5231.xxx	30÷90	✓		✓		✓			✓	3.9
	3.5207.006	6	✓		✓		✓		✓		3.5
	3.1221.050		✓				✓				



5,0ChLP SELBSTSCHNEIDENDE SCHRAUBE 3,5
5,0ChLP SELF-TAPPING SCREW 3.5
5,0ChLP TORNILLO AUTO-ROSCANTE 3,5



12	3.5200.012
14	3.5200.014
16	3.5200.016
18	3.5200.018
20	3.5200.020
22	3.5200.022
24	3.5200.024
26	3.5200.026
28	3.5200.028
30	3.5200.030
32	3.5200.032
34	3.5200.034
36	3.5200.036
38	3.5200.038
40	3.5200.040
42	3.5200.042
44	3.5200.044
46	3.5200.046
48	3.5200.048
50	3.5200.050
52	3.5200.052
54	3.5200.054
56	3.5200.056
58	3.5200.058
60	3.5200.060
65	3.5200.065
70	3.5200.070
75	3.5200.075
80	3.5200.080
85	3.5200.085

Kerndurchmesser core diameter diámetro del núcleo	2.8	
Skalierter Bohrer drill with scale broca graduada	2.8	40.5653.212
Führungshülse guide sleeve guía de la broca	5.0/2.8	40.5673.728
Schraubendrehereinsatz screwdriver tip punta de destornillador	T15	40.5677.000

4,5ChLP SCHRAUBE 2,4
4.5ChLP SCREW 2.4
4,5ChLP TORNILLO 2,4



10	3.5225.010
12	3.5225.012
14	3.5225.014
16	3.5225.016
18	3.5225.018
20	3.5225.020
22	3.5225.022
24	3.5225.024
26	3.5225.026
28	3.5225.028
30	3.5225.030
32	3.5225.032
34	3.5225.034
36	3.5225.036
38	3.5225.038
40	3.5225.040

Kerndurchmesser core diameter diámetro del núcleo	1.8	
Bohrer drill broca	1.8	40.2063.222
Führungshülse guide sleeve guía de la broca	5.0/1.8	40.5673.718
Schraubendrehereinsatz screwdriver tip punta de destornillador	T15	40.5677.000




SELBSTSCHNEIDENDE KORTIKALISSCHRAUBE 3,5
CORTICAL SELF-TAPPING SCREW 3.5
TORNILLO CORTICAL AUTO-ROSCANTE 3,5



12	3.1306.012
14	3.1306.014
16	3.1306.016
18	3.1306.018
20	3.1306.020
22	3.1306.022
24	3.1306.024
26	3.1306.026
28	3.1306.028
30	3.1306.030
32	3.1306.032
34	3.1306.034
36	3.1306.036
38	3.1306.038
40	3.1306.040
45	3.1306.045
50	3.1306.050
55	3.1306.055
60	3.1306.060
65	3.1306.065
70	3.1306.070
75	3.1306.075
80	3.1306.080
85	3.1306.085

Kerndurchmesser core diameter diámetro del núcleo	2.4	
Skalierter Bohrer drill with scale broca graduada	2.5	40.5912.212
Kompressionsführung compression guide guía de compresión	2.5/2.8	40.4804.700
Schraubendrehereinsatz screwdriver tip punta de destornillador	T15	40.5677.000

SELBSTSCHNEIDENDE KORTIKALISSCHRAUBE 2,7
CORTICAL SELF-TAPPING SCREW 2.7
TORNILLO CORTICAL AUTO-ROSCANTE 2,7



10	3.1220.010
12	3.1220.012
14	3.1220.014
16	3.1220.016
18	3.1220.018
20	3.1220.020
22	3.1220.022
24	3.1220.024
26	3.1220.026
28	3.1220.028
30	3.1220.030
32	3.1220.032
34	3.1220.034
36	3.1220.036
38	3.1220.038
40	3.1220.040

Kerndurchmesser core diameter diámetro del núcleo	1.8	
Bohrer drill broca	1.8	40.2063.181
Kompressionsführung compression guide guía de compresión	1.8	40.4897.018
Sechskant-Schraubendreher star screwdriver destornillador de estrella	T8	40.0669.000



5,0ChLP KONISCHE SELBSTSCHNEIDENDE SCHRAUBE 3,5
5.0ChLP CONICAL SELF-TAPPING SCREW 3.5
5,0ChLP TORNILLO CONICO AUTO-ROSCANTE 3,5



18	3.5205.018
20	3.5205.020
22	3.5205.022
24	3.5205.024
26	3.5205.026
28	3.5205.028
30	3.5205.030
35	3.5205.035
40	3.5205.040
45	3.5205.045
50	3.5205.050
55	3.5205.055
60	3.5205.060
65	3.5205.065
70	3.5205.070
75	3.5205.075
80	3.5205.080
85	3.5205.085
90	3.5205.090

5,0ChLP SPONGIOSASCHRAUBE 3,9
5.0ChLP CANCELLOUS SCREW 3.9
5,0ChLP TORNILLO ESPONJOSA 3,9



35	3.5231.035
40	3.5231.040
45	3.5231.045
50	3.5231.050
55	3.5231.055
60	3.5231.060
65	3.5231.065
70	3.5231.070
75	3.5231.075
80	3.5231.080

Kerndurchmesser core diameter diámetro del núcleo	2.8	
Skalierter Bohrer drill with scale broca graduada	2.8	40.5653.212
Führungshülse guide sleeve guía de la broca	5.0/2.8	40.5673.728
Schraubendrehereinsatz screwdriver tip punta de destornillador	T15	40.5677.000

Kerndurchmesser core diameter diámetro del núcleo	2.5	
Skalierter Bohrer drill with scale broca graduada	2.5	40.5912.212
Führungshülse guide sleeve guía de la broca	2.5	40.5673.725
Schraubendrehereinsatz screwdriver tip punta de destornillador	T15	40.5677.000



5,0ChLP VA SCHRAUBE 3,5
5.0ChLP SCREW VA 3.5
5,0ChLP TORNILLO VA 3,5



12	4.5236.012
14	4.5236.014
16	4.5236.016
18	4.5236.018
20	4.5236.020
22	4.5236.022
24	4.5236.024
26	4.5236.026
28	4.5236.028
30	4.5236.030
32	4.5236.032
34	4.5236.034
36	4.5236.036
38	4.5236.038
40	4.5236.040
42	4.5236.042
44	4.5236.044
46	4.5236.046
48	4.5236.048
50	4.5236.050
52	4.5236.052
54	4.5236.054
56	4.5236.056
58	4.5236.058
60	4.5236.060
65	4.5236.065
70	4.5236.070
75	4.5236.075
80	4.5236.080
85	4.5236.085
90	4.5236.090
95	4.5236.095

Kerndurchmesser
core diameter
diámetro del núcleo

2.8

Skalierter Bohrer
drill with scale
Broca graduada

2.8

40.5653.212

Führung VA
guide VA
Guía VA

2.8

40.8206.028

Schraubendrehereinsatz
screwdriver tip
punta de destornillador

T15

40.5677.000

Gewindebohrer
tap
Macho de roscar

3.5

40.5661.000

5,0ChLP WIDERSTANDSSCHRAUBE
5.0ChLP RESISTANCE SCREW
5,0ChLP TORNILLO DE DISTANCIA



6 3.5207.006

5,0ChLP SCHRAUBE FÜR CERCLAGE
5.0ChLP CERCLAGE SCREW
5,0ChLP TORNILLO PARA CERCLAJE




3.1221.050

*Maximaler Drahtdurchmesser 1,5mm
Wire diameter max. 1.5 mm
Diámetro máximo de alambre 1,5 mm

Palette für 5,0ChLP Platten 3.3157-3.3159 4x4 1/2H
Tray for plates 5.0ChLP 3.3157-3.3159 4x4 1/2H
Paleta para placas 5,0ChLP 3.3157-3.3159 4x4 1/2H



14.0205.404

	5,0 ChLP Schmale Kompressionsplatte Low-Contact 5.0ChLP Narrow compression plate with limited contact 5,0ChLP placa estrecha de compresión de contacto limitado	
		3.3157.5xx
	Stk. Pcs Uds.	7
	oder/ or/ o	
	5,0ChLP Schmale Kompressionsplatte 5.0ChLP Narrow compression plate 5,0ChLP Placa estrecha de compresión	
		3.3158.5xx
	Stk. Pcs Uds.	7
	oder/ or/ o	
	5,0ChLP Schmale Platte 5.0ChLP Narrow plate 5,0ChLP Placa estrecha	
		3.3159.0xx
	Stk. Pcs Uds.	7

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4022/3.4080 4x4 1/2H
Tray for plates 5.0ChLP 3.4022/3.4080 4x4 1/2H
Paleta para placas 5,0ChLP 3.4022/3.4080 4x4 1/2H


14.0205.402

	4,5ChLP Gerade Rekonstruktionsplatte 4.5ChLP Reconstruction straight plate 4,5ChLP Placa de reconstrucción recta	
		3.4022.8xx
	Stk. Pcs Uds.	7
	oder/ or/ o	
	4,5ChLP Schmale gerade Rekonstruktionsplatte 4.5ChLP Reconstruction straight plate narrow 4,5ChLP Placa de reconstrucción recta estrecha	
		3.4080.5xx
	Stk. Pcs Uds.	7

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4025-3.4030 4x4 1/2H
Tray for plates 5.0ChLP 3.4025-3.4030 4x4 1/2H
Paleta para placas 5,0ChLP 3.4025-3.4030 4x4 1/2H

14.0205.403



4,5ChLP Klavikula-Hakenplatte
4.5ChLP clavicular hook plate
4,5ChLP Placa clavicular con gancho



Stk. Pcs Uds.

3.4025.5xx ÷ 3.4029.5xx

7

4,5ChLP Klavikula-Hakenplatte
4.5ChLP clavicular hook plate
4,5ChLP Placa clavicular con gancho



Stk. Pcs Uds.


3.4026.5xx ÷ 3.4030.5xx

7

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7014/3.7015 4x4 1/2H
Tray for plates 5.0ChLP 3.7014/3.7015 4x4 1/2H
Paleta para placas 5,0ChLP 3.7014/3.7015 4x4 1/2H

14.0205.405



5,0ChLP Klavikula S-Platte
5.0ChLP Clavicle S plate
5,0ChLP Placa clavicular S



Stk. Pcs Uds.

3.7014.5xx

6

5,0ChLP Klavikula S-Platte
5.0ChLP Clavicle S plate
5,0ChLP Placa clavicular S



Stk. Pcs Uds.


3.7015.5xx

6

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7048/3.7049 4x4 1/2H
Tray for plates 5.0ChLP 3.7048/3.7049 4x4 1/2H
Paleta para placas 5,0ChLP 3.7048/3.7049 4x4 1/2H

14.0205.406



5,0ChLP S-platte für Klavikulaschaft
5,0ChLP Clavicle shaft s plate
5,0ChLP Placa en S para el eje de la clavícula



3.7048.5xx

Stk. Pcs Uds.

4

5,0ChLP S-platte für Klavikulaschaft
5,0ChLP Clavicle shaft s plate
5,0ChLP Placa en S para el eje de la clavícula



3.7049.5xx


Stk. Pcs Uds.

4

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7207/3.7212 4x4 1/2H
Tray for plates 5.0ChLP 3.7207-3.7212 4x4 1/2H
Paleta para placas 5,0ChLP 3.7207/3.7212 4x4 1/2H

14.0205.426



5,0ChLP Mediale Schulterblattplatte
5,0ChLP Scapular medial plate
5,0ChLP placa medial de escápula



3.7207.6xx

Stk. Pcs Uds.

2

5,0ChLP Mediale Schulterblattplatte
5,0ChLP Scapular medial plate
5,0ChLP placa medial de escápula



3.7208.6xx

Stk. Pcs Uds.

2

5,0ChLP Laterale Schulterblattplatte
5,0ChLP Scapular lateral plate
5,0ChLP placa lateral de escápula



3.7209.6xx

Stk. Pcs Uds.

2

5,0ChLP Laterale Schulterblattplatte
5,0ChLP Scapular lateral plate
5,0ChLP placa lateral de escápula



3.7210.6xx

Stk. Pcs Uds.

2

5,0ChLP Pfannen-Schulterblattplatte
5,0ChLP Scapular glenoid plate
5,0ChLP placa de glenoide



3.7211.600

Stk. Pcs Uds.

1

5,0ChLP Pfannen-Schulterblattplatte
5,0ChLP Scapular glenoid plate
5,0ChLP placa de glenoide



3.7212.600

Stk. Pcs Uds.


1

* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4034 4x4 1/2H
Tray for plates 5.0ChLP 3.4034 4x4 1/2H
Paleta para placas 5,0ChLP 3.4034 4x4 1/2H

14.0205.401

15.0205.401



5,0ChLP Humerusplatte
5.0ChLP Humerus plate
5,0ChLP Placa humeral



3.4034.6xx

Stk. Pcs Uds.

6

Führungsblock
Aiming block
Bloque guía

40.5671.000




Palette für 5,0ChLP Platten 3.4034 4x4 1/2H
Tray for plates 5.0ChLP 3.4034 4x4 1/2H
Paleta para placas 5,0ChLP 3.4034 4x4 1/2H


* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7243/3.7244 4x4 1/2H
Tray for plates 5.0ChLP 3.7243/3.7244 4x4 1/2H
Paleta para placas 5,0ChLP 3.7243/3.7244 4x4 1/2H

14.0205.429



5,0ChLP Distale Humerusplatte Y
5.0ChLP Distal humerus y plate
5,0ChLP Placa dist.del húmer. en Y




3.7243.6xx

Szt. Pcs Utr.

4

5,0ChLP Distale Humerusplatte Y
5.0ChLP Distal humerus y plate
5,0ChLP Placa dist.del húmer. en Y



3.7244.605

Szt. Pcs Utr.


4

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 4x4 1/2H
Tray for plates 5.0ChLP 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.419

15.0205.419



5,0ChLP Distale laterale Humerusplatte 5,0ChLP Distal lateral humerus plate 5,0ChLP Placa distal lateral del húmero	
	3.7037.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale posteriorlaterale Humerusplatte 5,0ChLP Distal posterior lateral humerus plate 5,0ChLP Placa distal posterolateral del húmero	
	3.7039.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale posteriormediale Humerusplatte 5,0ChLP Distal posterior medial humerus plate 5,0ChLP Placa distal posteromedial del húmero	
	3.7041.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale mediale Humerusplatte 5,0ChLP Distal medial humerus plate 5,0ChLP Placa distal medial del húmero	
	3.7071.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale dorsolaterale Humerusplatte 5,0ChLP Distal dorsolateral humerus plate 5,0ChLP Placa distal dorsolateral del húmero	
	3.7073.6xx
Stk. Pcs Uds.	3

Führungsblock [3.7037] Aiming block [3.7037] Bloque guía [3.7037]		40.8211.000	1
Führungsblock [3.7039] Aiming block [3.7039] Bloque guía [3.7039]		40.8213.000	1
Führungsblock [3.7041] Aiming block [3.7041] Bloque guía [3.7041]		40.8215.000	1
Führungsblock [3.7071] Aiming block [3.7071] Bloque guía [3.7071]		40.8217.000	1
Führungsblock [3.7073] Aiming block [3.7073] Bloque guía [3.7073]		40.8219.000	1


Palette für 5,0ChLP Platten 4x4 1/2H
Tray for plates 5.0ChLP 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 4x4 1/2H
Tray for plates 5.0ChLP 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

14.0205.420

15.0205.420



5,0ChLP Distale laterale Humerusplatte 5,0ChLP Distal lateral humerus plate 5,0ChLP Placa distal lateral del húmero	
	3.7036.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale posteriorlaterale Humerusplatte 5,0ChLP Distal posterior lateral humerus plate 5,0ChLP Placa distal posterolateral del húmero	
	3.7038.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale posteriormediale Humerusplatte 5,0ChLP Distal posterior medial humerus plate 5,0ChLP Placa distal posteromedial del húmero	
	3.7040.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale mediale Humerusplatte 5,0ChLP Distal medial humerus plate 5,0ChLP Placa distal medial del húmero	
	3.7072.6xx
Stk. Pcs Uds.	3
5,0ChLP Distale dorsolaterale Humerusplatte 5,0ChLP Distal dorsolateral humerus plate 5,0ChLP Placa distal dorsolateral del húmero	
	3.7074.6xx
Stk. Pcs Uds.	3

Führungsblock [3.7036] Aiming block [3.7036] Bloque guía [3.7036]		40.8210.000	1
Führungsblock [3.7038] Aiming block [3.7038] Bloque guía [3.7038]		40.8212.000	1
Führungsblock [3.7040] Aiming block [3.7040] Bloque guía [3.7040]		40.8214.000	1
Führungsblock [3.7072] Aiming block [3.7072] Bloque guía [3.7072]		40.8216.000	1
Führungsblock [3.7074] Aiming block [3.7074] Bloque guía [3.7074]		40.8218.000	1


Palette für 5,0ChLP Platten 4x4 1/2H
Tray for plates 5.0ChLP 4x4 1/2H
Paleta para placas 5,0ChLP 4x4 1/2H

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4097/3.4098 4x4 1/2H
Tray for plates 5.0ChLP 3.4097/3.4098 4x4 1/2H
Paleta para placas 5,0ChLP 3.4097/3.4098 4x4 1/2H

14.0205.409

15.0205.409



5,0ChLP Olekranonplatte
5.0ChLP Olecranon Plate
5,0ChLP Placa de olécranon



3.4097.5xx

Stk. Pcs Uds.

5

5,0ChLP Olekranonplatte
5.0ChLP Olecranon Plate
5,0ChLP Placa de olécranon



3.4098.5xx

Stk. Pcs Uds.

5

Führungsblock [3.4097]
Aiming block [3.4097]
Bloque guía [3.4097]



40.5731.100

1

Führungsblock [3.4098]
Aiming block [3.4098]
Bloque guía [3.4098]



40.5731.200


1

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4097/3.4098 4x4 1/2H
Tray for plates 5.0ChLP 3.4097/3.4098 4x4 1/2H
Paleta para placas 5,0ChLP 3.4097/3.4098 4x4 1/2H

Palette für 5,0ChLP Platten 3.7033-3.7035 4x4 1/2H
Tray for 5.0ChLP plates 3.7033-3.7035 4x4 1/2H
Paleta para placas 5,0ChLP 3.7033-3.7035 4x4 1/2H

14.0205.410



5,0ChLP Handgelenkplatte für Arthrodesis, gerade
5.0ChLP Wrist arthrodesis plate - straight
5,0ChLP Placa para artrodesis de muñeca - recta



3.7033.004

Stk. Pcs Uds.

2

5,0ChLP Handgelenkplatte f.Arthrodesis kurz. gebog.
5.0ChLP Wrist arthrodesis plate - short bend
5,0ChLP Placa para artrod.de muñe.- doblado corto



3.7034.004

Stk. Pcs Uds.

2

5,0ChLP Handgelenkplatte f.Arthrodesis stan. gebog.
5.0ChLP Wrist arthrodesis plate - standard bend
5,0ChLP Placa para artrod.de muñe.- doblado están.



3.7035.004

Stk. Pcs Uds.


2

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7067/3.7068 4x4 1/2H
Tray for plates 5.0ChLP 3.7067/3.7068 4x4 1/2H
Paleta para placas 5,0ChLP 3.7067/3.7068 4x4 1/2H

14.0205.411

15.0205.411



5,0ChLP Volare Radiusplatte
5.0ChLP Radius volar plate
5,0ChLP Placa radial de mano



3.7067.6xx

Stk. Pcs Uds.

5

5,0ChLP Volare Radiusplatte
5.0ChLP Radius volar plate
5,0ChLP Placa radial de mano



3.7068.6xx

Stk. Pcs Uds.

5

Führungsblock [3.7067]
Aiming block [3.7067]
Bloque guía [3.7067]



40.5742.100

1

Führungsblock [3.7068]
Aiming block [3.7068]
Bloque guía [3.7068]



40.5742.200

1


* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7067/3.7068 4x4 1/2H
Tray for plates 5.0ChLP 3.7067/3.7068 4x4 1/2H
Paleta para placas 5,0ChLP 3.7067/3.7068 4x4 1/2H

Palette für 5,0ChLP Platten 3.7091/3.7092 4x4 1/2H
Tray for plates 5.0ChLP 3.7091/3.7092 4x4 1/2H
Paleta para placas 5,0ChLP 3.7091/3.7092 4x4 1/2H

14.0205.421

15.0205.421



5,0ChLP proximale Tibiplatte T
5.0ChLP proximal tibial T plate
5,0ChLP placa tibial proximal T



3.7091.6xx

Stk. Pcs Uds.

5

5,0ChLP proximale Tibiplatte T
5.0ChLP proximal tibial T plate
5,0ChLP placa tibial proximal T



3.7092.6xx

Stk. Pcs Uds.

5

Zielblock L [3.7091]
Aiming block [3.7091]
Bloque guía L [3.7091]



40.8223.000

1

Zielblock R [3.7092]
Aiming block [3.7092]
Bloque guía R [3.7092]



40.8222.000

1


* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7091/3.7092 4x4 1/2H
Tray for plates 5.0ChLP 3.7091/3.7092 4x4 1/2H
Paleta para placas 5,0ChLP 3.7091/3.7092 4x4 1/2H

Palette für 5,0ChLP Platten 3.7202-3.7203 4x4 1/2H
Tray for plates 5.0ChLP 3.7202/3.7203 4x4 1/2H
Paleta para placas 5,0ChLP 3.7202/3.7203 4x4 1/2H

14.0205.424

15.0205.424



5,0ChLP Proximale laterale Tibiaplatte
5.0ChLP proximal lateral tibia plate
5,0ChLP placa de la tibia proximal lateral



3.7203.6xx

Stk. Pcs Uds.

3

5,0ChLP Proximale laterale Tibiaplatte
5.0ChLP proximal lateral tibia plate
5,0ChLP placa de la tibia proximal lateral



3.7202.6xx

Stk. Pcs Uds.

3

Zielblock L [3.7203]
Aiming block L [3.7203]
Bloque guía L [3.7203]



40.8227.000

1

Zielblock R [3.7202]
Aiming block R [3.7202]
Bloque guía R [3.7202]



40.8226.000

1


* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7202-3.7203 4x4 1/2H
Tray for 5,0ChLP plates 3.7202/3.7203 4x4 1/2H
Paleta para placas 5,0ChLP 3.7202/3.7203 4x4 1/2H

Palette für 5,0ChLP Plat.3.4011/4012.6xx 4x4 1/2H
Tray for plates 5.0ChLP 3.4011/4012.6xx 4x4 1/2H
Paleta para placas 5,0ChLP 3.4011/4012.6xx 4x4 1/2H

14.0205.427

15.0205.427



5,0ChLP Schmale Tibiaplatte
5.0ChLP proximal lateral tibia plate
5,0ChLP Placa proximal lateral de tibia



3.4011.6xx

Stk. Pcs Uds.

4

5,0ChLP Schmale Tibiaplatte
5.0ChLP proximal lateral tibia plate
5,0ChLP Placa proximal lateral de tibia



3.4012.6xx

Stk. Pcs Uds.

4

Zielblock L [3.4011.6xx]
Aiming block L [3.4011.6xx]
Bloque guía L (3.4011.6xx)



40.8233.000

1

Zielblock R [3.4012.6xx]
Aiming block R [3.4012.6xx]
Bloque guía R (3.4012.6xx)



40.8232.000

1


* Behälter ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Plat.3.4011/4012.6xx 4x4 1/2H
Tray for 5,0ChLP plates 3.4011/4012.6xx 4x4 1/2H
Paleta para placas 5,0ChLP 3.4011/4012.6xx 4x4 1/2H

Palette für 5,0ChLP Platten 3.7094/3.7095 4x4 1/2H
Tray for plates 5.0ChLP 3.7094/3.7095 4x4 1/2H
Paleta para placas 5,0ChLP 3.7094/3.7095 4x4 1/2H

14.0205.422

15.0205.422



5,0 ChLP schmale posteriore Tibiaplatt
5.0ChLP posterior narrow tibia plate
5,0ChLP placa tibial posterior estrecha



3.7094.6xx

Stk. Pcs Uds.

4

5,0 ChLP breite posteriore Tibiaplatt
5.0ChLP posterior wide tibia plate
5,0ChLP placa tibial posterior ancha



3.7095.6xx

Stk. Pcs Uds.

4

Führungsblock [3.7094]
Aiming block [3.7094]
Bloque guía [3.7094]



40.8224.000

1

Führungsblock [3.7095]
Aiming block [3.7095]
Bloque guía [3.7095]



40.8225.000


1

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes


Palette für 5,0ChLP Platten 3.7094/3.7095 4x4 1/2H
Tray for plates 5.0ChLP 3.7094/3.7095 4x4 1/2H
Paleta para placas 5,0ChLP 3.7094/3.7095 4x4 1/2H

Palette für 5,0ChLP Platten 3.4033 4x4 1/2H
Tray for plates 5.0ChLP 3.4033 4x4 1/2H
Paleta para placas 5,0ChLP 3.4033 4x4 1/2H

14.0205.414



5,0ChLP Distale mediale Tibiaplatt
5.0ChLP Distal medial tibial plate
5,0ChLP Placa tibial distal medial



3.4033.5xx

Stk. Pcs Uds.


2

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4039/3.4040 4x4 1/2H
Tray for plates 5.0ChLP 3.4039/3.4040 4x4 1/2H
Paleta para placas 5,0ChLP 3.4039/3.4040 4x4 1/2H

14.0205.415

15.0205.415



5,0ChLP Distale mediale Tibiaplatte
5.0ChLP Distal medial tibial plate
5,0ChLP Placa tibial d/m izquierda



3.4039.5xx

Stk. Pcs Uds.

4

5,0ChLP Distale mediale Tibiaplatte
5.0ChLP Distal medial tibial plate
5,0ChLP Placa tibial d/m izquierda



3.4040.5xx

Stk. Pcs Uds.

4

Führungsblock [3.4039]
Aiming block [3.4039]
Bloque guía [3.4039]



40.5726.100

1

Führungsblock [3.4040]
Aiming block [3.4040]
Bloque guía [3.4040]



40.5726.200

1


* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4039/3.4040 4x4 1/2H
Tray for plates 5.0ChLP 3.4039/3.4040 4x4 1/2H
Paleta para placas 5,0ChLP 3.4039/3.4040 4x4 1/2H

Palette für 5,0ChLP Platten 3.4051/3.4052 4x4 1/2H
Tray for plates 5.0ChLP 3.4051/3.4052 4x4 1/2H
Paleta para placas 5,0ChLP 3.4051/3.4052 4x4 1/2H

14.0205.428

15.0205.428



5,0ChLP Distale L-Tibiaplatte
5.0ChLP Distal tibial L plate
5,0ChLP Placa tibial distal L



3.4051.6xx

Stk. Pcs Uds.

4

5,0ChLP Distale L-Tibiaplatte
5.0ChLP Distal tibial L plate
5,0ChLP Placa tibial distal L



3.4052.6xx

Stk. Pcs Uds.

4

Führungsblock [3.4051.6xx]
Aiming block [3.4051.6xx]
Bloque guía [3.4051.6xx]



40.8235.000

1

Führungsblock [3.4052.6xx]
Aiming block [3.4052.6xx]
Bloque guía [3.4052.6xx]



40.8234.000


1

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.4051/3.4052 4x4 1/2H
Tray for plates 5.0ChLP 3.4051/3.4052 4x4 1/2H
Paleta para placas 5,0ChLP 3.4051/3.4052 4x4 1/2H

Palette für 5,0ChLP Platten 3.7029/3.7030 4x4 1/2H
Tray for plates 5.0ChLP 3.7029/3.7030 4x4 1/2H
Paleta para placas 5,0ChLP 3.7029/3.7030 4x4 1/2H

14.0205.417



5,0 ChLP Distale laterale Fibulaplatte
5.0ChLP Distal lateral fibula plate
5,0ChLP Placa distal lateral de peroné



3.7029.5xx /6xx

Stk. Pcs Uds.

3

5,0 ChLP Distale laterale Fibulaplatte
5.0ChLP Distal lateral fibula plate
5,0ChLP Placa distal lateral de peroné



3.7030.5xx /6xx


Stk. Pcs Uds.

3

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes

Palette für 5,0ChLP Platten 3.7200/3.7201 4x4 1/2H
Tray for plates 5.0ChLP 3.7200/3.7201 4x4 1/2H
Paleta para placas 5,0ChLP 3.7200/3.7201 4x4 1/2H

14.0205.423



5,0ChLP Kalkaneusplatte
5.0ChLP Calcaneal plate
5.0ChLP Placa para calcáneo



3.7200.6xx

Stk. Pcs Uds.

6

5,0ChLP Kalkaneusplatte
5.0ChLP Calcaneal plate
5.0ChLP Placa para calcáneo



3.7201.6xx

Stk. Pcs Uds.

6

* Palette ohne Implantaten/ Tray does not include implants/ Bandeja no incluye implantes



3.3155.606



	Len	Ti
6	131	3.3155.606
8	173	3.3155.608
10	215	3.3155.610
12	257	3.3155.612
14	299	3.3155.614



		TiA	Co			VA			
	3.5210.xxx	✓		✓	✓		✓		5.0
	3.5232.xxx	✓		✓	✓			✓	5.4
	3.5216.xxx	✓		✓			✓		5.0
	3.1471.xxx	✓		✓			✓		4.5
	4.5244.xxx		✓	✓	✓	✓	✓		5.0
	3.5219.008	✓		✓	✓				5.0
	3.1221.070	✓			✓				





Palette für 7,0ChLP gerade Platten
Palette for 7.0ChLP plates - straight
Paleta para placas 7,0ChLP rectas

40.5704.550



			Ti
	Len	L	R
4	138	3.4023.604	3.4024.604
6	180	3.4023.606	3.4024.606
8	221	3.4023.608	3.4024.608
10	263	3.4023.610	3.4024.610
4-16	138 387		

		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.5221.xxx	✓		✓	✓			✓	✓	7.3
	3.5224.xxx	✓		✓				✓	✓	7.3
	3.5228.xxx	✓		✓	✓			✓		6.5
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					



Palette für 7,0ChLP Platten 3.4023/4024
Palette for 7.0ChLP plates - 3.4023/3.4024
Paleta para placas 7,0ChLP 3.4023/4024

40.5704.510



Zielblock L [3.4023]
Aiming block L [3.4023]
Bloque guía L (3.4023)

40.5725.100



Zielblock R [3.4024]
Aiming block R [3.4024]
Bloque guía R (3.4024)

40.5725.200



Platte 3.4023.606 Probeimplantat
Plate 3.4023.606 trial
Plantilla para placa 3.4023.606

43.4023.606

Platte 3.4024.606 Probeimplantat
Plate 3.4024.606 trial
Plantilla para placa 3.4024.606

43.4024.606



			Ti
	Len	L	R
4	129	3.4013.604	3.4014.604
5	150	3.4013.605	3.4014.605
6	171	3.4013.606	3.4014.606
8	213	3.4013.608	3.4014.608
4-12	129-297		



	3.5210.xxx	✓	✓	✓	✓		5.0
	3.5232.xxx	✓	✓	✓		✓	5.4
	3.5216.xxx	✓	✓		✓		5.0
	3.1471.xxx	✓	✓		✓		4.5
	4.5244.xxx		✓	✓	✓	✓	5.0
	3.5219.008	✓	✓	✓			5.0
	3.1221.070	✓		✓			



Palette für 7,0ChLP Platten 3.4013/3.4014 mit Instrum.
Palette for 7.0ChLP plates - 3.4013/3.4014;
Paleta para placas 7,0ChLP 3.4013/3.4014 (con instrumentos)

40.5704.520

Zielblock L [3.4013]
Aiming block L [3.4013]
Bloque guía L (3.4013)

40.5724.100

Zielblock R [3.4014]
Aiming block R [3.4014]
Bloque guía R (3.4014)


40.5724.200

Platte 3.4013.605 Probeimplantat
Plate 3.4013.605 trial
Plantilla para placa 3.4013.605

43.4013.605


Platte 3.4014.605 Probeimplantat
Plate 3.4014.605 trial
Plantilla para placa 3.4014.605

43.4014.605



3.4009.605

			Ti
	Len	L	R
4	116	3.4009.604	3.4010.604
6	158	3.4009.606	3.4010.606
8	200	3.4009.608	3.4010.608
10	242	3.4009.610	3.4010.610
4-14	116 - 326		



		TiA	Co			VA			
	3.5210.xxx	✓		✓	✓		✓		5.0
	3.5232.xxx	✓		✓	✓			✓	5.4
	3.5216.xxx	✓		✓			✓		5.0
	3.1471.xxx	✓		✓			✓		4.5
	4.5244.xxx		✓	✓	✓	✓	✓		5.0
	3.5219.008	✓		✓	✓				5.0
	3.1221.070	✓			✓				



Palette für 7,0ChLP Platten 3.4009/4010
Palette for 7.0ChLP plates - 3.4009/3.4010
Paleta para placas 7,0ChLP 3.4009/4010

40.5704.540



Platte 3.4009.606 Probeimplantat
Plate 3.4009.606 trial
Plantilla para placa 3.4009.606

43.4009.606

Platte 3.4010.606 Probeimplantat
Plate 3.4010.606 trial
Plantilla para placa 3.4010.606

43.4010.606



		Ti	
	Len	L	R
3	131	3.4089.603	3.4090.603
4	152	3.4089.604	3.4090.604
6	194	3.4089.606	3.4090.606
8	236	3.4089.608	3.4090.608
3-10	131 - 278		



		TiA	Co			VA			
	3.5210.xxx	✓		✓	✓		✓		5.0
	3.5232.xxx	✓		✓	✓			✓	5.4
	3.5216.xxx	✓		✓			✓		5.0
	3.1471.xxx	✓		✓			✓		4.5
	4.5244.xxx		✓	✓	✓	✓	✓		5.0
	3.5219.008	✓		✓	✓				5.0
	3.1221.070	✓			✓				



Palette für 7,0ChLP Platten 3.4089/3.4090 mit Instr.
Palette for 7.0ChLP plates - 3.4089/3.4090
Paleta para placas 7,0ChLP 3.4089/3.4090 (con instrumentos)

40.5704.560



Zielblock L [3.4089]
Aiming block L [3.4089]
Bloque guía L [3.4089]

40.5709.100

Zielblock R [3.4090]
Aiming block R [3.4090]
Bloque guía R [3.4090]

40.5709.200




Platte 3.4089.604 Probeimplantat
Plate 3.4089.604 trial
Plantilla para placa 3.4089.604


43.4089.604

Platte 3.4090.604 Probeimplantat
Plate 3.4089.604 trial
Plantilla para placa 3.4090.604

43.4090.604



			Ti
	Len	L	R
4	134	3.7055.604	3.7054.604
6	176	3.7055.606	3.7054.606
8	218	3.7055.608	3.7054.608
10	260	3.7055.610	3.7054.610
3-13	113 323		



		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					



Palette für 7,0ChLP Platten 3.7054/3.7055 mit Instrum.
Palette for 7.0ChLP plates - 3.7054/3.7055
Paleta para placas 7,0ChLP 3.7054/3.7055 (con
instrumentos)

40.6292.000

Zielblock L [3.7055]
Aiming block L [3.7055]
Bloque guía L [3.7055]

40.8221.000


Zielblock R [3.7054]
Aiming block R [3.7054]
Bloque guía R [3.7054]

40.8220.000




			
1	84	3.7253.601	3.7254.601
4	138	3.7253.604	3.7254.604
7	192	3.7253.607	3.7254.607

		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.1471.xxx	✓		✓			✓			4.5
	3.5242.xxx	✓		✓				✓		5.4
	4.5244.xxx		✓	✓	✓	✓	✓			5.0



Palette für Platte 7,0ChLP 5x4 1/2H
 Tray for plates 7.0ChLP 5x4 1/2H
 Paleta para placas 7,0ChLP 5x4 1/2H

14.0207.412




Palette für Platte 7,0ChLP 5x4 1/2H
 Tray for plates 7.0ChLP 5x4 1/2H
 Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413



			Ti
	Len	L	R
1	126	3.7255.601	3.7256.601
4	180	3.7255.604	3.7256.604




		TiA	Co			VA	C		
	3.5210.xxx	✓		✓	✓		✓		5.0
	3.1471.xxx	✓		✓			✓		4.5
	3.5242.xxx	✓		✓				✓	5.4
	4.5244.xxx		✓	✓	✓	✓	✓		5.0



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H

14.0207.412



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413



		Ti	
	Len	L	R
1	92	3.7257.601	3.7258.601
4	146	3.7257.604	3.7258.604




		TiA	Co		VA			
	3.5210.xxx	✓	✓	✓	✓			5.0
	3.1471.xxx	✓	✓		✓			4.5
	3.5242.xxx	✓	✓			✓		5.4
	4.5244.xxx		✓	✓	✓	✓		5.0



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H

14.0207.412



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413



			Ti
	Len	L	R
1	101	3.7259.601	3.7260.601
4	155	3.7259.604	3.7260.604




		TiA	Co			VA	C		
	3.5210.xxx	✓		✓	✓		✓		5.0
	3.1471.xxx	✓		✓			✓		4.5
	3.5242.xxx	✓		✓				✓	5.4
	4.5244.xxx		✓	✓	✓	✓	✓		5.0



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7,0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H

14.0207.412



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7,0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413



		Ti		
	Len	L	R	
1	93	3.7261.601	3.7262.601	
4	147	3.7261.604	3.7262.604	




		TiA	Co		VA				
	3.5210.xxx	✓	✓	✓	✓				5.0
	3.1471.xxx	✓	✓		✓				4.5
	3.5242.xxx	✓	✓			✓			5.4
	4.5244.xxx		✓	✓	✓	✓			5.0



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H

14.0207.412



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413



		Ti
	Len	L R
1	49	3.7263.601
4	103	3.7263.604




		TiA	Co			VA	C			
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.1471.xxx	✓		✓			✓			4.5
	3.5242.xxx	✓		✓				✓		5.4
	4.5244.xxx		✓	✓	✓	✓	✓			5.0





Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.412



Palette für Platte 7,0ChLP 5x4 1/2H
Tray for plates 7.0ChLP 5x4 1/2H
Paleta para placas 7,0ChLP 5x4 1/2H


14.0207.413

15.0207.207	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 7,0ChLP 4x4 1/2H Tray for 7.0ChLP instrument set 4x4 1/2H Paleta para instrumental 7,0ChLP 4x4 1/2H	14.0207.207	1
	Kirschnerdraht 2,0/210 Kirschner wire 2.0/210 Aguja de Kirschner 2,0/210	40.4815.210	4
	Skalierter Bohrer 3,2/210 Drill with scale 3.2/210 Broca graduada 3,2/210	40.5650.212	2
	Skalierter Bohrer 4,0/210 Drill with scale 4.0/210 Broca graduada 4,0/210	40.5651.212	2
	Plattenrepositionsschraube 4,0/180 Setting-compressing screw 4.0/180 Tornillo de ajuste-compresión 4,0/180	40.5706.740	1
	Führung VA 4,0 Guide VA 4.0 Guía VA 4,0	40.8207.040	1
	Führungshülse 7,0/3,2 Guide sleeve 7.0/3.2 Guía de la broca 7,0/3,2	40.5705.732	1
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	40.5705.740	2
	Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Llave dinamométrica con trinquete 4Nm	40.6660.000	1
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta de destornillador T25-1/4	40.5684.200	1
	Tiefenmesslehre Depth measure Medidor de profundidad	40.4639.550	1
	Plattenrepositionsschraube 2,0/75 Setting-compressing screw 2.0/75 Tornillo de ajuste-compresión 2,0/075	40.8209.075	4
	Dynamometrischer Verbinder 4Nm Torque connector 4Nm Conector dinamométrico 4Nm	40.5927.040	1

15.0207.110	Name Name Nombre	Kat.-Nr., Catalogue no. Número decatálogo	Stk. Pcs Uds.
	7,0ChLP Behälter 9x4 1/2H 7.0ChLP container 9x4 1/2H Contenedor 7,0ChLP 9x4 1/2H	14.0207.110	1
	Zielbügel für Platte 3.7254 Plate 3.7254 targeter Guía de la placa 3.7254	40.8270.000	1
	Zielbügel für Platte 3.7253 Plate 3.7253 targeter Guía de la placa 3.7253	40.8271.000	1
	Schutzhülse 9,0/7,0 Protective guide 9.0/7.0 Guía de protección 9,0/7,0	40.5693.770	2
	Trokar 7,0 Trocar 7.0 Trocar 7,0	40.5695.770	1
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	40.5690.740	2
	Schutzhülse 10,0/8,0 Protective guide 10.0/8.0 Guía de protección 10,0/8,0	40.5694.780	1
	Trokar 8,0 Trocar 8.0 Trocar 8,0	40.5696.780	1
	Führungshülse 8,0/3,2 Guide sleeve 8.0/3.2 Guía de la broca 8,0/3,2	40.5691.732	1
	Skalierter Bohrer 3,2/300 Drill with scale 3.2/300 Broca graduada 3,2/300	40.5650.302	1
	Skalierter Bohrer 4,0/300 Drill with scale 4.0/300 Broca graduada 4,0/300	40.5651.302	1
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T 25-1/4 Punta de destornillador T25-1/4	40.5684.000	1
	Plattenrepositionsschraube 4,0 Setting-compressing screw 4.0 Tornillo de ajuste-compresión 4,0	40.5698.740	1

BEHÄLTER FÜR SCHRAUBEN 7,0ChLP 4X2H
STAND FOR 7.0ChLP SCREWS 4X2H
BANDEJA PARA TORNILLOS 7,0ChLP 4X2H

14.0207.506

	7,0ChLP Selbstschneidende Schraube 5,0 7,0ChLP self-tapping screw 5,0 7,0ChLP Tornillo Auto-Roscante 5,0 																							
	L [mm]	16	18	20	22	24	26	28	30	32	34	36	38	40	42	44	46	48	50	52	54	56		
	Stk. Pcs Uds.	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4		
	L [mm]	58	60	65	70																			
	Stk. Pcs Uds.	4	4	4	4																			

* Behälter ohne Implantaten/ Stand does not include implants/ Bandeja no incluye implantes

BEHÄLTER FÜR SCHRAUBEN 7,0ChLP 4X2H
STAND FOR 7.0ChLP SCREWS 4X2H
BANDEJA PARA TORNILLOS 7,0ChLP 4X2H

14.0207.507

Kortikalisschraube Ø4,5
Cortical screw Ø4.5
Tornillo cortical Ø4,5


L [mm]	16	18	20	22	24	26	28	30	32	34	36	38	40	42	44	46	48	50	52	54	56
Stk. Pcs Uds.	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2
L [mm]	58	60	65	70																	
Stk. Pcs Uds.	2	2	2	2																	

Konische Spongiaschraube ø 5,4
Conical cancellous screw Ø5.4
Tornillo cónico esponjosa ø5,4


L [mm]	30	35	40	45	50	55	60	65	70	75	80	85	90	95							
Szt.Pcs Utr.	2	2	2	2	2	2	2	2	2	2	2	2	2	2							

* Behälter ohne Implantaten/ Stand does not include implants/ Bandeja no incluye implantes

VORSCHLAG ZUR ZUSAMMENSETZUNG DES 7,0 ChLP FÜR ARTHRODESE DES SPRUNGGELENKS
7.0 ChLP SET COMPLETATION PROPOSAL FOR TALUS ARTHRODESIS
EJEMPLO DE COMPOSICIÓN PARA EL SET DE 7,0ChLP PARA ARTRODESIS DE TALUS



Nº	Name Name Nombre	Catalogue no. Kat.-Nr. Nº de catálogo	Pcs Stk. Pcs.
1	7,0ChLP Behälterdeckel 9x4H 7.0ChLP container lid 9x4H Tapa de contenedor 7,0ChLP 9x4H	14.0207.102	1
2	7,0ChLP Behälter 9x4H 7.0ChLP container 9x4H Contenedor 7,0ChLP 9x4H	14.0207.101	1
3	Palette für Platte 7,0ChLP 5x4 1/2H Tray for 7.0ChLP plates 5x4 1/2H Paleta para placas 7,0ChLP 5x4 1/2H	14.0207.413	1
4	Behälter für Schrauben 7,0ChLP 4x2H Stand for 7.0ChLP screws 4x2H Contenedor para tornillos 7,0ChLP 4x2H	14.0207.412	1
5	Behälter für Schrauben 7,0ChLP 4x2H Stand for 7.0ChLP screws 4x2H Contenedor para tornillos 7,0ChLP 4x2H	14.0207.506	1
6	Behälter für Schrauben 7,0ChLP 4x2H Stand for 7.0ChLP screws 4x2H Contenedor para tornillos 7,0ChLP 4x2H	14.0207.507	1
7	Instrumentarium für 7,0ChLP 9x4 1/2H Instrument set for 7.0ChLP 9x4 1/2H Set de instrumental para 7,0ChLP 9x4 1/2H	15.0207.110	1
8	Instrumentarium für 7,0ChLP 4x4 1/2H Instrument set for 7.0ChLP 4x4 1/2H Set de instrumental para 7,0ChLP 4x4 1/2H	15.0207.207	1



1	Len	L	R	Ti	KURZ SHORT CORTO
6	222	3.7221.606	3.7220.606		
8	274	3.7221.608	3.7220.608		
10	326	3.7221.610	3.7220.610		
12	378	3.7221.612	3.7220.612		

2	Len	L	R	Ti	LANG LONG LARGO
6	229	3.7223.606	3.7222.606		
8	281	3.7223.608	3.7222.608		
10	333	3.7223.610	3.7222.610		
12	384	3.7223.612	3.7222.612		

* Anzahl der Löcher in der Mitte des Einsatzschafes

* Number of holes in the middle shaft part of the plate

* Número de agujeros en el centro del eje de inserción



		TiA	Co				VA			
	3.5210.xxx	✓		✓	✓	✓				5.0
	3.1471.xxx	✓		✓		✓				4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0




Palette für 7,0ChLP Platten - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225
 Palette for 7.0ChLP plates - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225
 Paleta para placas 7,0ChLP - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225

40.6299.100



		Ti	KURZ SHORT CORTO
①	Len		
2	42	3.7224.600	
②			LANG LONG LARGO
2	50	3.7225.600	



		TiA	Co				VA			
	3.5210.xxx	✓		✓	✓	✓				5.0
	3.1471.xxx	✓		✓		✓				4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0



Palette für 7,0ChLP Platten - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225
Palette for 7.0ChLP plates - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225
Paleta para placas 7,0ChLP - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225

40.6299.100



CERKLAGEDRAHT MIT KLEMME - GEFLEHT - 2,0X600
WIRE FOR CERCLAGE WITH CLAMP-WEAVE-2.0X600
ALAMBRE PARA CERCLAJE CON GRAPA - TEJIDO - 2.0X600






Len

600

Co

4.5250.600S

15.0506.101		Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
		Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta de destornillador T25-1/4	40.5684.220	1
		Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Mango dinamométrico T con trinquete 4Nm	40.6660.000	1
		Führung 60 Guide 60 Guía 60	40.8265.060	1
		Führung 76 Guide 76 Guía 76	40.8265.076	1
		Griff Handle Mango	40.8266.000	1
		Schneidezange Cutting pliers Alicates para cortar	40.8267.000	1
		Instrument für Draht Wire instrument Instrumento para alambre	40.8268.000	1


15.0506.101	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Zange für Plattenschneiden Plate cutting pliers Alicates para cortar placas	40.4650.000	1
	Führung Guide Guía	40.8272.000	1
	Behälter - Cerklage 9x4 H Container - cerclage 9x4H Contenedor - cerclaje 9x4H	14.0506.101	1
	Behälterdeckel - Cerklage 9x4H Container lid - cerclage 9x4H Tapa de contenedor - cerclaje 9x4H	14.0506.102	1



3.7023.605

			Ti
	Len	L	R
2	132	3.7023.602	3.7022.602
4	174	3.7023.604	3.7022.604
6	216	3.7023.606	3.7022.606
8	258	3.7023.608	3.7022.608
2-16	132 - 426		

		TiA	Co			VA	C			
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.5221.xxx	✓		✓	✓			✓	✓	7.3
	3.5224.xxx	✓		✓				✓	✓	7.3
	3.5228.xxx	✓		✓	✓			✓		6.5
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					



Palette für 7,0ChLP Platten 3.7022/7023 mit Instrum.
Palette for 7.0ChLP plates - 3.7022/3.7023;
Paleta para placas 7,0ChLP 3.7022/7023 (con
instrumentos)

40.5704.590



Zielblock L [3.7023]
Aiming block L [3.7023]
Bloque guía L (3.7023)

40.5732.100


Zielblock R [3.7022]
Aiming block R [3.7022]
Bloque guía R (3.7022)

40.5732.200



Platte 3.7022.604/3.7023.604
Probeimplantat
Plate 3.7022.604/ 3.7023.604 trial
Plantilla para placa 3.7022.604/3.7023.604

43.7022.604



3.7062.612



	Len	Ti
10	209	3.7062.610
12	251	3.7062.612
14	292	3.7062.614
16	333	3.7062.616




		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					



Palette für 7,0ChLP Platten 3.7062
Palette for 7.0ChLP plates - 3.7062;
Paleta para placas 7,0ChLP 3.7062

40.5704.900

7,0 ChLP SCHMALE KOMPRESSIONSPLATTE LOW CONTACT 7.0ChLP NARROW COMPRESSION PLATE WITH LIMITED CONTACT 7,0ChLP PLACA ESTRECHA DE COMPRESIÓN DE CONTACTO LIMITADO




		
5	90	3.7093.605
6	108	3.7093.606
8	144	3.7093.608
10	180	3.7093.610
12	216	3.7093.612
14	252	3.7093.614

										
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					



Palette für 7,0ChLP Platten 3.7093
Palette for 7.0ChLP plates - 3.7093
Paleta para placas 7,0ChLP 3.7093


40.6296.000




Nº	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
①	7,0ChLP Behälterdeckel 9x4H/ 7.0ChLP container lid 9x4H/ Tapa de contenedor 7,0ChLP 9x4H	14.0207.102	1
②	7,0ChLP Behälter 9x4 H/ 7.0ChLP container 9x4 H/ Contenedor 7,0ChLP 9x4 H	14.0207.101	1
③	Instrumentarium für 7,0ChLP 4x4 1/2H/ Instrument set for 7.0ChLP 4x4 1/2H/ Set de instrumental para 7,0ChLP 4x4 1/2H	15.0207.201	1
④	Instrumentarium für 7,0ChLP 4x4 1/2H/ Instrument set for 7.0ChLP 4x4 1/2H/ Set de instrumental para 7,0ChLP 4x4 1/2H	15.0207.202	1
⑤	Behälter für Schrauben 7,0ChLP 4x2H/ Stand for 7.0ChLP screws 4x2H/ Contenedor para tornillos 7,0ChLP 4x2H	14.0207.501	1
⑥	Behälter für Schrauben 7,0ChLP 4x2H / Stand for 7.0ChLP screws 4x2H/ Contenedor para tornillos 7,0ChLP 4x2H	14.0207.502	1
⑦	Palette 4x1 1/2H/ Tray 4x1 1/2H/ Paleta 4x1 1/2H	14.0000.201	2

Nº	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
①	7,0ChLP Behälterdeckel 9x4H 7.0ChLP container lid 9x4H Tapa de contenedor 7,0ChLP 9x4H	14.0207.102	1
②	7,0ChLP Behälter 9x4 H 7.0ChLP container 9x4 H Contenedor 7,0ChLP 9x4 H	14.0207.101	2

INSTRUMENTARIUM FÜR 7,0ChLP 4x4 1/2H
INSTRUMENT SET FOR 7.0ChLP 4x4 1/2H
SET DE INSTRUMENTAL PARA 7,0ChLP 4x4 1/2H

15.0207.201	Name Name Nombre	Kat.-Nr Catalogue no. Nº de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 7,0ChLP 4x4 1/2H Tray for 7.0ChLP instrument set 4x4 1/2H Paleta para instrumental 7,0ChLP 4x4 1/2H	14.0207.201	1
	Kirschnerdraht 2,0/210 Kirschner wire 2.0/210 Aguja KIRSCHNER 2,0/210	40.4815.210	4
	Skalierter Bohrer 3,2/210 Drill with scale 3.2/210 Broca graduada 3,2/210	40.5650.212	2
	Skalierter Bohrer 4,0/210 Drill with scale 4.0/210 Broca graduada 4,0/210	40.5651.212	2
	Kanülierter skalierter Bohrer 5,0/2.2/210 Cannulated drill with scale 5.0/2.2/210 Broca canulada graduada 5,0/2.2/210	40.5652.212	1
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta de destornillador T25-1/4	40.5684.200	1
	Kanül. Schraubendrehereinsatz T30-1/4 Cannulated screwdriver tip T30-1/4 Punta de destornillador canulada T30-1/4	40.5685.200	1
	Kompressionsführung 3,2 Compression guide 3.2 Guía de compresión 3,2	40.4802.732	1
	Führungshülse 9,0/5,0 Guide sleeve 9.0/5.0 Guía de la broca 9,0/5,0	40.5654.750	1
	Führungshülse 9,0/3,2 Guide sleeve 9.0/3.2 Guía de la broca 9,0/3,2	40.5707.732	1
	Führungshülse 5,0/2,0 Guide sleeve 5.0/2.0 Guía de la broca 5,0/2,0	40.5654.120	1
	Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Llave dinamométrica con trinquete en T 4Nm	40.6660.000	1
	Schutzhülse 9/7 Protective guide 9/7 Guía de protección 9/7	40.5708.000	2
	Führungshülse 7,0/3,2 Guide sleeve 7.0/3.2 Guía de la broca 7,0/3,2	40.5705.732	2
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	40.5705.740	4
	Tiefenmesslehre Depth measure Medidor de profundidad	40.4639.550	1


INSTRUMENTARIUM FÜR 7,0ChLP 4x4 1/2H
INSTRUMENT SET FOR 7.0ChLP 4x4 1/2H
SET DE INSTRUMENTAL PARA 7,0ChLP 4x4 1/2H

15.0207.202	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 7,0ChLP 4x4 1/2H Tray for 7.0ChLP instrument set 4x4 1/2H Paleta para instrumental 7,0ChLP 4x4 1/2H	14.0207.202	1
	Biegegerät für Platten 7,0 Plates bender 7.0 Grifa para doblar placas 7,0	40.4643.700	2
	Plattenrepositionsschraube 4,0/180 Setting-compressing screw 4.0/180 Tornillo de ajuste-compresión 4,0/180	40.5706.740	2
	Dreiarmer Schraubendrehereinsatz 7,0ChLP Tripod screwdriver tip 7.0ChLP Punta de destornillador de tres brazos 7,0ChLP	40.6271.700	1
	Einsatz T25 mit Greifer Screwdriver tip T25 with holder Punta de destornillador T25 con sujeción	40.6259.200	1
	Gewindebohrer 7,0ChLP - 5,0 Tap 7.0ChLP - 5.0 Macho de roscar 7,0ChLP - 5,0	40.5646.200	1
	Gewindebohrer HA 4,5 Cortical tap HA 4.5 Macho de roscar HA4,5	40.5647.200	1
	Messlehre für Schraubenlänge Screw length measure Medidor de tornillos	40.5675.500	1
OPTIONALES INSTRUMENT / OPTIONAL INSTRUMENT / INSTRUMENTO OPCIONAL			
	Dynamometrischer Verbinder 4Nm Torque connector 4Nm Conector dinamométrico 4Nm	40.5927.040	1

INSTRUMENTARIUM FÜR 7,0ChLP 4x4 1/2H
INSTRUMENT SET FOR 7.0ChLP 4x4 1/2H
SET DE INSTRUMENTAL PARA PLACAS 7,0ChLP 4x4 1/2H





15.0207.208	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Palette für Instrumentarium 7,0ChLP 4x4 1/2H Tray for 7.0ChLP instrument set 4x4 1/2H Paleta para instrumental 7,0ChLP 4x4 1/2H	14.0207.208	1
	Kirschnerdraht 2,0/210 Kirschner wire 2.0/210 Aguja de Kirschner 2,0/210	40.4815.210	4
	Skalierter Bohrer 3,2/210 Drill with scale 3.2/210 Broca graduada 3,2/210	40.5650.212	2
	Skalierter Bohrer 4,0/210 Drill with scale 4.0/210 Broca graduada 4,0/210	40.5651.212	2
	Kanülierter skalierter Bohrer 5,0/2,2/210 Cannulated drill with scale 5.0/2.2/210 Broca canulada graduada 5,0/2,2/210	40.5652.212	1
	Plattenrepositionsschraube 4,0/180 Setting-compressing screw 4.0/180 Tornillo de ajuste-compresión 4,0/180	40.5706.740	1
	Führung VA 4,0 Guide VA 4.0 Guía VA 4,0	40.8207.040	1
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	40.5705.740	3
	Führungshülse 7,0/3,2 Guide sleeve 7.0/3.2 Guía de la broca 7,0/3,2	40.5705.732	2
	Führungshülse 9,0/5,0 Guide sleeve 9.0/5.0 Guía de la broca 9,0/5,0	40.5654.750	1
	Führungshülse 5,0/2,0 Guide sleeve 5.0/2.0 Guía de la broca 5,0/2,0	40.5654.120	1
	Schutzhülse 9/7 Protective guide 9/7 Guía de protección 9/7	40.5708.000	2

15.0207.208	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo	Stk. Pcs Uds.
	Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Mango dinamométrico T con embrague 4Nm	40.6660.000	1
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta de destornillador T25-1/4	40.5684.200	1
	Kanül. Schraubendrehereinsatz T30-1/4 Cannulated screwdriver tip T30-1/4 Punta de destornillador canulada T30-1/4	40.5685.200	1
	Tiefenmesslehre Depth measure Medidor de profundidad	40.4639.550	1
Optionales Instrument Optional instrument Instrumento opcional			
	Dynamometrischer Verbinder 4Nm Torque connector 4Nm Conector dinamométrico 4Nm	40.5927.040	0





BEHÄLTER FÜR SCHRAUBEN 7,0ChLP 4X2H
STAND FOR 7.0ChLP SCREWS 4X2 H
BANDEJA PARA TORNILLOS 7,0ChLP 4X2H

14.0207.501

	7,0ChLP Schraube 5,0 7.0ChLP Screw 5.0 7,0ChLP tornillo 5,0			
	L [mm]	16 18 20 22 24 26 28 30 32 34 36 38 40 42 44 46 48 50		
	Stk. Pcs Uds.	4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4 4		
	L [mm]	52 54 56 58 60 65 70 75 80 85 90 95		
	7,0ChLP Spongiaschraube 5,4 7.0ChLP Cancellous screw 5,4 7,0ChLP tornillo esponjosa 5,4			
	L [mm]	35 40 45 50 55 60 65 70 75 80		
	Stk. Pcs Uds.	2 2 2 2 2 2 2 2 2 2		
	7,0ChLP Schraube für Cerclage 7.0ChLP cerclage screw 7,0ChLP tornillo para cerclaje			
		Stk. Pcs Uds.	3	

BEHÄLTER FÜR SCHRAUBEN 7,0ChLP 4X2H
STAND FOR 7.0ChLP SCREWS 4X2 H
BANDEJA PARA TORNILLOS 7,0ChLP 4X2H

14.0207.502

	Kortikalisschraube 4,5 Cortical screw 4.5 Tornillo cortical 4,5			
	L [mm]	20 22 24 26 28 30 32 34 36 38 40 42 44 46 48 50 52 54		
	Stk. Pcs Uds.	2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2		
	L [mm]	56 58 60 65 70 75 80 85 90		
	7,0ChLP Spongiaschraube 6,5 7.0ChLP Cancellous screw 6.5 7,0ChLP tornillo esponjosa 6,5			
	L [mm]	35 40 45 50 55 60 65 70 75 80		
	Stk. Pcs Uds.	2 2 2 2 2 2 2 2 2 2		
	7,0ChLP Kanülierte Schraube 7,3 7.0ChLP Cannulated screw 7.3 7,0ChLP tornillo canulado 7,3			
	7,0ChLP Kanül. konische Schraube 7,3 7.0ChLP Conical cannulated screw 7.3 7,0ChLP tornillo cónico canulado 7,3			
	L [mm]	50 55 60 65 70 75 80 85 90 95		
	Stk. Pcs Uds.	1 1 1 1 1 1 1 1 1 1		
	L [mm]	50 55 60 65 70 75 80 85 90 95		
		Stk. Pcs Uds.	1 1 1 1 1 1 1 1 1 1	






	Len	TiA
2	59	3.5314.002
3	80	3.5314.003



	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5316.xxx	✓								7.3
	4.5244.xxx		✓	✓	✓	✓	✓			5.0
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					


40.6600.000	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Zielbügel Targeter arm Brazo guía	1	40.6601.000
	Distaler Zielbügelauflauf Distal targeter Guía distal	1	40.6602.000
	Mehrloch-Führung Polyhole guide Guía multiagujero	1	40.6612.000
	Schutzführung 9/2 Protective guide 9/2 Guía de protección 9/2	3	40.6603.000
	Führungsdraht 2,8/350 Guide rod 2.8/350 Varilla guía 2,8/350	2	40.6610.000
	Kirschnerdraht 2,0/385 Kirschner wire 2.0/385 Aguja de Kirschner 2,0/385	3	40.6604.000
	Messlehre für Schrauben Screw length measure Medidor de tornillos	1	40.6605.000
	Bohrer 7,3/5 Drill 7.3/5 Broca 7,3/5	1	40.6606.000
	Schlüssel für Schrauben 7,3 Screw wrench 7.3 Llave para tornillo 7,3	1	40.6607.000
	Verbinder Connector Conector	1	40.6608.000
	Schlüssel für Schraubenentfernung 7,3 Wrench for 7.3 screw removal Llave para destornillar tornillos 7,3	1	40.6609.000
	Trokar 7,0 Trocar 7.0 Trocar 7,0	1	40.5695.570

40.6600.000	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue no. Número de catálogo
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	2	40.5690.540
	Schutzhülse 9,0/7,0 Protective guide 9.0/7.0 Guía de protección 9,0/7,0	2	40.5693.570
	Skalierter Bohrer 3,2/300 Drill with scale 3.2/300 Broca graduada 3,2/300	1	40.5650.301
	Skalierter Bohrer 4,0/300 Drill with scale 4.0/300 Broca graduada 4,0/300	2	40.5651.301
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta del destornillador T25-1/4	1	40.5684.000
	Schraubendrehereinsatz T30-1/4 Screwdriver tip T30-1/4 Punta del destornillador T30-1/4	1	40.6611.000
	Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Mango dinamométrico T con embrague 4Nm	1	40.6660.000
	Führungshülse 8,0/3,2 Guide sleeve 8.0/3.2 Guía de la broca 8,0/3,2	1	40.5691.532
	Schutzhülse 10,0/8,0 Protective guide 10.0/8.0 Guía de protección 10,0/8,0	1	40.5694.580
	Trokar 8,0 Trocar 8,0 Trocar 8,0	1	40.5696.580
	Verschlusschraube für Zielbügelaufsatz Targeter end cap Tapón guía	15	40.5612.000
	Messlehre für Schraubenlänge Screw length measure Medidor de tornillos	1	40.5700.000
	AO-Verbindungsstück - 7,0ChLP Connector AO - 7.0ChLP Conector AO - 7,0ChLP	1	40.4898.070
	Gewindebohrer 7,0ChLP - 5,0 Tap 7.0ChLP - 5.0 Macho de roscar 7,0ChLP - 5,0	1	40.5646.000
	Gewindebohrer HA4,5 Cortical tap HA 4.5 Macho de roscar HA4,5	1	40.5647.000
	Behälter Stand Contenedor	1	40.6619.000



	TiA		
	H	L	R
3		3.7065.030	3.7066.030
5		3.7065.050	3.7066.050
7.5		3.7065.075	3.7066.075
9		3.7065.090	3.7066.090
10		3.7065.100	3.7066.100
11		3.7065.110	3.7066.110
12.5		3.7065.125	3.7066.125
15		3.7065.150	3.7066.150
17.5		3.7065.175	3.7066.175

		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.5221.xxx	✓		✓	✓			✓	✓	7.3
	3.5224.xxx	✓		✓				✓	✓	7.3
	3.5228.xxx	✓		✓	✓			✓		6.5
	4.5244.xxx		✓	✓	✓	✓	✓			5.0



Anwenden mit zusätzlichen Instrumenten


Use with additional instruments:

Utilizar con instrumentos adicionales:

40.5350.600 – Instrumentarium für Osteotomie


Instrument set for osteotomy

Set de instrumental para osteotomía



H	Len	
5		3.3538.505
7.5		3.3538.507
9		3.3538.509
10		3.3538.510
11	125	3.3538.511
12.5		3.3538.512
15		3.3538.515
17.5		3.3538.517

		TiA	Co			VA				
	3.5210.xxx	✓		✓	✓		✓			5.0
	3.5232.xxx	✓		✓	✓			✓		5.4
	3.5216.xxx	✓		✓			✓			5.0
	3.5221.xxx	✓		✓	✓			✓	✓	7.3
	3.5224.xxx	✓		✓				✓	✓	7.3
	3.5228.xxx	✓		✓	✓			✓		6.5
	3.1471.xxx	✓		✓			✓			4.5
	4.5244.xxx		✓	✓	✓	✓	✓			
	3.5219.008	✓		✓	✓					5.0
	3.1221.070	✓			✓					








Anwenden mit zusätzlichen Instrumenten

Use with additional instruments:

Utilizar con instrumentos adicionales:

40.5350.600 – Instrumentarium für Osteotomie
Instrument set for osteotomy
Set de instrumental para osteotomía

40.5350.600	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue No. Número de catálogo
	Führungsnadel mit Loch Guide pin with eyelet Perno guía con ojo	1	40.3943.100
	Zielgerät für Osteotomie Targeter for osteotomy Guía para osteotomía	1	40.3974.200
	Instrumentarium für Osteotomie Instrument set for osteotomy Set para osteotomía	1	40.5360.000
	Räumer 20° Reamer 20° Escariador 20°	2	40.5357.020
	Räumer 16° Reamer 16° Escariador 16°	2	40.5357.016
	Tiefenmesslehre Depth measure Medidor de profundidad	1	40.4639.700
	Kompressionsführung 3,2 Compression guide 3.2 Guía de compresión 3,2	1	40.4802.732
	Kirschnerdraht 2,0/220 Kirschner wire 2.0/220 Aguja de Kirschner 2,0/220	2	40.4815.220
	Drehmomentgriff T mit Kupplung 4Nm Torque limiting ratchet handle T 4Nm Llave dinamométrica con trinquete en T 4Nm	1	40.6660.000
	Führungsdraht Guide rod Varilla guía	2	40.5353.000
	Osteotom 25 Osteotome 25 Osteótomo 25	1	40.5361.000
	Osteotom 35 Osteotome 35 Osteótomo 35	1	40.5362.000
	Applikator Applicator Aplicador	1	40.5369.000
	Skalierter Bohrer 3,2/220 Drill with scale 3.2/220 Broca graduada 3,2/220	1	40.5650.222
	Skalierter Bohrer 4,0/220 Drill with scale 4.0/220 Broca graduada 4,0/220	1	40.5651.222

40.5350.600	Name Name Nombre	Stk. Pcs Uds.	Kat.-Nr. Catalogue No. Número de catálogo
	Messlehre für Schraubenlänge Screw length measure Medidor de tornillos	1	40.5675.100
	Schraubendrehereinsatz T25-1/4 Screwdriver tip T25-1/4 Punta del destornillador T25-1/4	1	40.5684.000
	Kanül. Schraubendrehereinsatz T30-1/4 Cannulated screwdriver tip T30-1/4 Punta del destornillador canulada T30-1/4	1	40.5685.000
	Schraubendrehereinsatz S3,5-1/4 Screwdriver tip S3.5-1/4 Punta del destornillador S3,5-1/4	1	40.5686.000
	Kanülierter Schraubendrehereinsatz S5-1/4 Cannulated screwdriver tip S5-1/4 Punta del destornillador canulada S5-1/4	1	40.5687.000
	Führungshülse 7,0/4,0 Guide sleeve 7.0/4.0 Guía de la broca 7,0/4,0	1	40.5705.740
	Führungshülse 9,0/3,2 Guide sleeve 9.0/3.2 Guía de la broca 9,0/3,2	1	40.5707.732
	Behälter für Implantate für Osteotomie Stand for implants for osteotomy Contenedor para implantes de osteotomía	1	40.5358.000
	Instrumentariumbehälter für Osteotomie Stand for instrument set of osteotomy Contenedor para el set de instrumental de osteotomía	1	40.5359.600
	Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm Perforated aluminum lid 1/1 595x275x15mm Gray Tapa de aluminio perforada 1/1 595x275x15mm Gris	1	12.0750.200
	Behälter mit solidem Boden 1/1 595x275x86mm Container with solid bottom 1/1 595x275x86mm Contenedor de fondo sólido 1/1 595x275x86mm	1	12.0750.100

		Len	TiA	Co			VA				
	3.5210.xxx	16÷110	✓		✓	✓		✓			5.0
	3.5216.xxx	30÷90	✓		✓			✓			5.0
	3.5232.xxx	30÷90	✓		✓	✓			✓		5.4
	3.5221.xxx	30÷100	✓		✓	✓			✓	✓	7.3
	3.5222.xxx	30÷100	✓		✓	✓			✓	✓	7.3
	3.5223.xxx	30÷95	✓		✓				✓	✓	7.3
	3.5224.xxx	30÷95	✓		✓				✓	✓	7.3
	3.5242.xxx		✓		✓				✓		5.4
	3.5228.xxx	30÷95	✓		✓	✓			✓		6.5
	3.1471.xxx	16÷110	✓		✓			✓			4.5
	4.5244.xxx	16÷110		✓	✓	✓	✓	✓			5.0
	3.5219.008	8	✓		✓	✓					5.0
	3.1221.070		✓			✓					



7,0ChLP SELBSTSCHNEIDENDE SCHRAUBE 5,0
7,0ChLP SELF-TAPPING SCREW 5.0
7,0ChLP TORNILLO AUTO-ROSCANTE 5,0



16	3.5210.016
18	3.5210.018
20	3.5210.020
22	3.5210.022
24	3.5210.024
26	3.5210.026
28	3.5210.028
30	3.5210.030
32	3.5210.032
34	3.5210.034
36	3.5210.036
38	3.5210.038
40	3.5210.040
42	3.5210.042
44	3.5210.044
46	3.5210.046
48	3.5210.048
50	3.5210.050
52	3.5210.052
54	3.5210.054
56	3.5210.056
58	3.5210.058
60	3.5210.060
65	3.5210.065
70	3.5210.070
75	3.5210.075
80	3.5210.080
85	3.5210.085
90	3.5210.090
95	3.5210.095

Kerndurchmesser core diameter diámetro del núcleo	4.0	
Skalierter Bohrer drill with scale broca graduada	4.0	40.5651.212
Führungshülse guide sleeve guía de la broca	7.0/4.0	40.5705.740
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.200

SELBSTSCHNEIDENDE KORTIKALISSCHRAUBE 4,5
CORTICAL SELF-TAPPING SCREW 4.5
TORNILLO CORTICAL AUTO-ROSCANTE 4,5



20	3.1471.020
22	3.1471.022
24	3.1471.024
26	3.1471.026
28	3.1471.028
30	3.1471.030
32	3.1471.032
34	3.1471.034
36	3.1471.036
38	3.1471.038
40	3.1471.040
42	3.1471.042
44	3.1471.044
46	3.1471.046
48	3.1471.048
50	3.1471.050
52	3.1471.052
54	3.1471.054
56	3.1471.056
58	3.1471.058
60	3.1471.060
65	3.1471.065
70	3.1471.070
75	3.1471.075
80	3.1471.080
85	3.1471.085
90	3.1471.090
95	3.1471.095

Kerndurchmesser core diameter diámetro del núcleo	3.0	
Skalierter Bohrer drill with scale broca graduada	3.2	40.5650.212
Kompressionsführung compression guide guía de compresión	3.2	40.4802.732
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.200



7,0ChLP SPONGIOSASCHRAUBE 5,4
7,0ChLP CANCELLOUS SCREW 5.4
7,0ChLP TORNILLO ESPONJOSA 5,4



35	3.5232.035
40	3.5232.040
45	3.5232.045
50	3.5232.050
55	3.5232.055
60	3.5232.060
65	3.5232.065
70	3.5232.070
75	3.5232.075
80	3.5232.080

7,0ChLP KANÜLIERTE KONISCHE SCHRAUBE 7,3X16
7,0ChLP CONICAL CANNULATED SCREW 7.3X16
7,0ChLP TORNILLO CÓNICO CANULADO 7,3X16



50	3.5224.050
55	3.5224.055
60	3.5224.060
65	3.5224.065
70	3.5224.070
75	3.5224.075
80	3.5224.080
85	3.5224.085
90	3.5224.090
95	3.5224.095

Kerndurchmesser core diameter diámetro del núcleo	3.2	
Skalierter Bohrer drill with scale broca graduada	3.2	40.5650.212
Führungshülse guide sleeve guía de la broca	3.2	40.5705.732
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.200

Kerndurchmesser core diameter diámetro del núcleo	5.2	
Kirschnerdraht Kirschner wire aguja de Kirschner	2.0	40.4815.210
Führungshülse guide sleeve guía de la broca	5.0/2.0 9.0/5.0	40.5654.120 40.5654.750
kanülierter skalierter Bohrer cannulated drill with scale broca canulada graduada	5.0/2.2	40.5652.212
kanülierter Schraubendrehereinsatz cannulated screwdriver tip punta de destornillador canulada	T30	40.5685.200



7,0ChLP SELBSTSCHN. KONISCHE SCHRAUBE 5,0
7,0ChLP CONICAL SELF - TAPPING SCREW 5.0
7,0ChLP TORNILLO CÓNICO AUTO-ROSCANTE 5,0



30	3.5216.030
35	3.5216.035
40	3.5216.040
45	3.5216.045
50	3.5216.050
55	3.5216.055
60	3.5216.060
65	3.5216.065
70	3.5216.070
75	3.5216.075
80	3.5216.080
85	3.5216.085
90	3.5216.090

7,0ChLP SPONGIOSASCHRAUBE 6,5
7,0ChLP CANCELLOUS SCREW 6.5
7,0ChLP TORNILLO ESPONJOSA 6,5



30	3.5228.030
35	3.5228.035
40	3.5228.040
45	3.5228.045
50	3.5228.050
55	3.5228.055
60	3.5228.060
65	3.5228.065
70	3.5228.070
75	3.5228.075
80	3.5228.080
85	3.5228.085
90	3.5228.090
95	3.5228.095

Kerndurchmesser core diameter diámetro del núcleo	4.0	
Skalierter Bohrer drill with scale broca graduada	4.0	40.5651.212
Führungshülse guide sleeve guía de la broca	7.0/4.0	40.5705.740
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.200

Kerndurchmesser core diameter diámetro del núcleo	3.2	
Skalierter Bohrer drill with scale broca graduada	3.2	40.5650.212
Führungshülse guide sleeve guía de la broca	9.0/3.2	40.5707.732
kanülierter Schraubendrehereinsatz cannulated screwdriver tip punta de destornillador canulada	T30	40.5685.200




7,0ChLP KANÜLIERTE SCHRAUBE 7,3
7,0ChLP CANNULATED SCREW 7.3
7,0ChLP TORNILLO CANULADO 7,3



L	A
30	3.5221.030
35	3.5221.035
40	3.5221.040
45	3.5221.045
50	3.5221.050
55	3.5221.055
60	3.5221.060
65	3.5221.065
70	3.5221.070
75	3.5221.075
80	3.5221.080
85	3.5221.085
90	3.5221.090
95	3.5221.095
100	3.5221.100

7,0ChLP KONISCHE SPONGIOSASCHRAUBE Ø5,4
7,0 ChLP CONICAL CANCELLOUS SCREW Ø5.4
7,0 ChLP TORNILLO CÓNICO ESPONJOSA Ø5,4



L	A	
30	16	3.5242.030
35		3.5242.035
40		3.5242.040
45		3.5242.045
50		3.5242.050
55	32	3.5242.055
60		3.5242.060
65		3.5242.065
70		3.5242.070
75		3.5242.075
80		3.5242.080
85		3.5242.085
90		3.5242.090
95		3.5242.095

Kerndurchmesser core diameter diámetro del núcleo	5.2	
Kirschnerdraht kirschner wire aguja de Kirschner	2.0	40.4815.210
Führungshülse guide sleeve guía de la broca	5.0/2.0 9.0/5.0	40.5654.120 40.5654.750
kanülierter skaliert Bohrer cannulated drill with scale broca canulada graduada	5.0/2.2	40.5652.210
kanülierter Schraubendrehereinsatz cannulated screwdriver tip punta de destornillador canulada	T30	40.5685.200

Kerndurchmesser core diameter diámetro del núcleo	3.2	
Skaliert Bohrer drill with scale broca graduada	3.2	40.5650.222
Führungshülse guide sleeve guía de la broca	3.2	40.5705.732
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.000



7,0ChLP VA SCHRAUBE 5,0
7,0ChLP SCREW VA 5.0
7,0ChLP TORNILLO VA 5,0



16	4.5244.016
18	4.5244.018
20	4.5244.020
22	4.5244.022
24	4.5244.024
26	4.5244.026
28	4.5244.028
30	4.5244.030
32	4.5244.032
34	4.5244.034
36	4.5244.036
38	4.5244.038
40	4.5244.040
42	4.5244.042
44	4.5244.044
46	4.5244.046
48	4.5244.048
50	4.5244.050
52	4.5244.052
54	4.5244.054
56	4.5244.056
58	4.5244.058
60	4.5244.060
65	4.5244.065
70	4.5244.070
75	4.5244.075
80	4.5244.080
85	4.5244.085
90	4.5244.090
95	4.5244.095
100	4.5244.100
105	4.5244.105
110	4.5244.110

Kerndurchmesser core diameter diámetro del núcleo	4.0	
Skaliert Bohrer drill with scale broca graduada	4.0	40.5651.212
Führung VA guide VA Guía VA	4.0	40.8207.040
Schraubendrehereinsatz screwdriver tip punta de destornillador	T25	40.5684.200

7,0ChLP WIDERSTANDSSCHRAUBE
7,0ChLP DISTANCE SCREW
7,0ChLP TORNILLO DE DISTANCIA



8	3.5219.008
---	------------

7,0ChLP SCHRAUBE FÜR CERCLAGE
7,0ChLP CERCLAGE SCREW
7,0ChLP TORTILLO PARA CERCLAJE



	3.1221.070
--	------------


*Maximaler Durchmesser 1,5mm
Wire diameter max. 1.5 mm
Diámetro máximo de alambre 1,5 mm

Dreiarmliger Schraubendreher screwdriver Destornillador de tres brazos	40.4746.000
--	-------------

TELESKOPSCHRAUBE 7,3
TELESCOPIC SCREW 7.3
TORNILLO TELESCÓPICO 7,3



70	3.5316.070
75	3.5316.075
80	3.5316.080
85	3.5316.085
90	3.5316.090
95	3.5316.095
100	3.5316.100
105	3.5316.105
110	3.5316.110
115	3.5316.115
120	3.5316.120

Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo			Stk. Pcs, Uds.
Palette für 7,0ChLP gerade Platten Palette for 7.0ChLP straight plates Paleta para placas 7,0ChLP rectas				
	Palette	Palette	Paleta	40.5704.350 1
	Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100 1
	Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200 1
				40.5704.550

ohne Implantaten; implants not included; implantes no incluidos

Palette für 7,0ChLP Platten 3.4023/3.4024 (mit Instrum.)
Palette for 7.0ChLP plates - 3.4023/3.4024 (w/instr.)
Paleta para placas 7,0ChLP 3.4023/3.4024 (con instrumentos)

	Zielblock L [3.4023]	Aiming block [3.4023]	Bloque guía L (3.4023)	40.5725.100 1	40.5704.310
	Zielblock R [3.4024]	Aiming block R [3.4024]	Bloque guía R (3.4024)	40.5725.200 1	
	Schutzhülse 9/7	Protective guide 9/7	Guía de protección 9/7	40.5708.000 2	
	Palette	Palette	Paleta	40.5704.410 1	40.5704.510
	Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100 1	
	Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200 1	

ohne Implantaten, mit zusätzlichen Instrumenten
implants not included; with additional instruments
implantes no incluidos, con instrumentos adicionales

Palette für 7,0ChLP Platten 3.4013/3.4014 (mit Instrum.)
Palette for 7.0ChLP plates - 3.4013/3.4014 (w/instr.)
Paleta para placas 7,0ChLP 3.4013/3.4014 (con instrumentos)

	Zielblock L [3.4013]	Aiming block L [3.4013]	Bloque guía L (3.4013)	40.5724.100 1	40.5704.320
	Zielblock R [3.4014]	Aiming block R [3.4014]	Bloque guía R (3.4014)	40.5724.200 1	
	Schutzhülse 9/7	Protective guide 9/7	Guía de protección 9/7	40.5708.000 2	
	Palette	Palette	Paleta	40.5704.420 1	40.5704.520
	Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100 1	
	Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200 1	

ohne Implantaten, mit zusätzlichen Instrumenten
implants not included; with additional instruments
implantes no incluidos, con instrumentos adicionales

Palette für 7,0ChLP Platten 3.4009/4010
Palette for 7.0ChLP plates - 3.4009/3.4010
Paleta para placas 7,0ChLP 3.4009/4010



Palette	Palette	Paleta	40.5704.340	1
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1

40.5704.540

ohne Implantaten; implants not included; implantes no incluidos

Palette für 7,0ChLP Platten 3.4089/3.4090 (mit Instrum.)
Palette for 7.0ChLP plates - 3.4089/3.4090 (w/instr.)
Paleta para placas 7,0ChLP 3.4089/3.4090 (con instrumentos)



Zielblock L [3.4089]	Aiming block L [3.4089]	Bloque guía L (3.4089)	40.5709.100	1
Zielblock R [3.4090]	Aiming block R [3.4090]	Bloque guía R (3.4090)	40.5709.200	1
Schutzhülse 9/7	Protective guide 9/7	Guía de protección 9/7	40.5708.000	2
Palette	Palette	Paleta	40.5704.460	1
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1

40.5704.360

40.5704.560

ohne Implantaten, mit zusätzlichen Instrumenten
implants not included; with additional instruments
implantes no incluidos, con instrumentos adicionales

Palette für 7,0ChLP Platten 3.7054/3.7055 (mit Instrum.)
Palette for 7.0ChLP plates - 3.7054/3.7055 (w/instr.)
Paleta para placas 7,0ChLP 3.7054/3.7055 (con instrumentos)



Zielblock L [3.7055]	Aiming block L [3.7055]	Bloque guía L (3.7055)	40.8221.000	1
Zielblock R [3.7054]	Aiming block R [3.7054]	Bloque guía R (3.7054)	40.8220.000	1
Schutzhülse 9/7	Protective guide 9/7	Guía de protección 9/7	40.5708.000	2
Palette	Palette	Paleta	40.6293.000	1
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1

40.6292.000

ohne Implantaten, mit zusätzlichen Instrumenten
implants not included; with additional instruments
implantes no incluidos, con instrumentos adicionales

Paleta na płytki 7,0ChLP - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225;
Palette for 7.0ChLP plates - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225;
Поддон для пластин 7,0ChLP - 3.7220/3.7221, 3.7222/3.7223, 3.7224, 3.7225



Palette	Palette	Paleta	40.6299.000	1
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1

40.6299.100

ohne Implantaten; implants not included; implantes no incluidos

Palette für 7,0ChLP Platten 3.7022/7023 (mit Instrum.)
Palette for 7.0ChLP plates - 3.7022/3.7023 (w/instr.)
Paleta para placas 7,0ChLP 3.7022/7023 (con instrumentos)



Zielblock L [3.7023]	Aiming block L [3.7023]	Bloque guía L (3.7023)	40.5732.100	1	40.5704.390
Zielblock R [3.7022]	Aiming block R [3.7022]	Bloque guía R (3.7022)	40.5732.200	1	
Schutzhülse 9/7	Protective guide 9/7	Guía de protección 9/7	40.5708.000	2	
Palette	Palette	Paleta	40.5704.490	1	
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1	40.5704.590
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1	

ohne Implantaten, mit zusätzlichen Instrumenten
implants not included; with additional instruments
implantes no incluidos, con instrumentos adicionales

Palette für 7,0ChLP Platten 3.7062
Palette for 7.0ChLP plates - 3.7062
Paleta para placas 7,0ChLP 3.7062



Palette	Palette	Paleta	40.5704.800	1	40.5704.900
Behälter mit solidem Boden 1/1 595x275x86mm	Container with solid bottom 1/1 595x275x86mm	Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100	1	
Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm	Perforated aluminum lid 1/1 595x275x15mm Gray	Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200	1	

ohne Implantaten; implants not included; implantes no incluidos

Palette für 7,0ChLP Platten 3.7093
Palette for 7.0ChLP plates - 3.7093
Paleta para placas 7,0ChLP 3.7093



Palette	Palette	Paleta	40.6296.100	1	40.6296.000
Behälter mit solidem Boden 1/2 306x272x85mm	Container with solid bottom 1/2 306x272x85mm	Contenedor de fondo sólido 1/2 306x272x85mm	12.0751.100	1	
Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm	Perforated aluminum lid 1/2 306x272x15mm Gray	Tapa de aluminio perforada 1/2 306x272x15mm Gris	12.0751.200	1	


ohne Implantaten; implants not included; implantes no incluidos

Behälter für Schrauben 7,0ChLP
Stand for screws 7.0ChLP
Contenedor para tornillos 7,0ChLP



Behälter für Schrauben 7,0ChLP	Stand for screws 7.0ChLP	Contenedor para tornillos 7,0ChLP	40.5749.600	1	40.5749.700
Behälter mit solidem Boden 1/2 306x272x184mm	Container with solid bottom 1/2 306x272x135mm	Contenedor de fondo sólido 1/2 306x272x184mm	12.0751.103	1	
Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm	Perforated aluminum cover 1/2 306x272x15mm Gray	Tapa de aluminio perforada 1/2 306x272x15mm Gris	12.0751.200	1	

ohne Implantaten; implants not included; implantes no incluidos


	Name Name Nombre	Links/Rechts Left/Right Izq./Der.	Catalogue No. Kat.-Nr. Número de catálogo
	Platte 3.4031.030 Probeimplantat Plate 3.4031.030 trial Plantilla para placa 3.4031.030	L	43.4031.030
	Platte 3.4032.030 Probeimplantat Plate 3.4032.030 trial Plantilla para placa 3.4032.030	R	43.4032.030
	Platte 3.4037.603 Probeimplantat Plate 3.4037.603 trial Plantilla para placa 3.4037.603	L	43.4037.603
	Platte 3.4038.603 Probeimplantat Plate 3.4038.603 trial Plantilla para placa 3.4038.603	R	43.4038.603
	Platte 3.4041.603 Probeimplantat Plate 3.4041.603 trial Plantilla para placa 3.4041.603	L	43.4041.603
	Platte 3.4042.603 Probeimplantat Plate 3.4042.603 trial Plantilla para placa 3.4042.603	R	43.4042.603
	Platte 3.4043.603 Probeimplantat Plate 3.4043.603 trial Plantilla para placa 3.4043.603	L	43.4043.603
	Platte 3.4044.603 Probeimplantat Plate 3.4044.603 trial Plantilla para placa 3.4044.603	R	43.4044.603
	Platte 3.4045.603 Probeimplantat Plate 3.4045.603 trial Plantilla para placa 3.4045.603	L	43.4045.603
	Platte 3.4046.603 Probeimplantat Plate 3.4046.603 trial Plantilla para placa 3.4046.603	R	43.4046.603
	Platte 3.4047.603 Probeimplantat Plate 3.4047.603 trial Plantilla para placa 3.4047.603	L/R	43.4047.603
	Platte 3.4048.603 Probeimplantat Plate 3.4048.603 trial Plantilla para placa 3.4048.603	L/R	43.4048.603

	Name Name Nombre	Links/Rechts Left/Right Izq./Der.	Catalogue No. Kat.-Nr. Número de catálogo
	Platte 3.7024.502 Probeimplantat Plate 3.7024.502 trial Plantilla para placa 3.7024.502	R	43.7024.502
	Platte 3.7025.502 Probeimplantat Plate 3.7025.502 trial Plantilla para placa 3.7025.502	L	43.7025.502
	Platte 3.7026.502 Probeimplantat Plate 3.7026.502 trial Plantilla para placa 3.7026.502	R	43.7026.502
	Platte 3.7027.502 Probeimplantat Plate 3.7027.502 trial Plantilla para placa 3.7027.502	L	43.7027.502
	Platte 3.7028.502 Probeimplantat Plate 3.7028.502 trial Plantilla para placa 3.7028.502	L/R	43.7028.502
	Platte 3.7229.606 Probeimplantat Plate 3.7229.606 trial Plantilla para placa 3.7229.606	L	43.7229.606
	Platte 3.7228.606 Probeimplantat Plate 3.7228.606 trial Plantilla para placa 3.7228.606	R	43.7228.606
	Platte 3.7235.606 Probeimplantat Plate 3.7235.606 trial Plantilla para placa 3.7235.606	L	43.7235.606
	Platte 3.7234.606 Probeimplantat Plate 3.7234.606 trial Plantilla para placa 3.7234.606	R	43.7234.606

	Name Name Nombre	Links/Rechts Left/Right Izq./Der.	Catalogue No. Kat.-Nr. Número de catálogo
	Platte 3.4011.605 Probeimplantat Plate 3.4011.605 trial Plantilla para placa 3.4011.605	L	43.4011.605
	Przymiar płytki 3.4012.605 Plate 3.4012.605 trial Plantilla para placa 3.4012.605	R	43.4012.605
	Platte 3.4039.608 Probeimplantat Plate 3.4039.608 trial Plantilla para placa 3.4039.608	L	43.4039.608
	Platte 3.4040.608 Probeimplantat Plate 3.4040.608 trial Plantilla para placa 3.4040.608	R	43.4040.608
	Platte 3.4051.607 Probeimplantat Plate 3.4051.607 trial Plantilla para placa 3.4051.607	L	43.4051.607
	Platte 3.4052.607 Probeimplantat Plate 3.4052.607 trial Plantilla para placa 3.4052.607	R	43.4052.607
	Platte 3.4097.504 Probeimplantat Plate 3.4097.504 trial Plantilla para placa 3.4097.504	L	43.4097.504
	Platte 3.4098.504 Probeimplantat Plate 3.4098.504 trial Plantilla para placa 3.4098.504	R	43.4098.504
	Platte 3.7014.505 Probeimplantat Plate 3.7014.505 trial Plantilla para placa 3.7014.505	R	43.7014.505
	Platte 3.7015.505 Probeimplantat Plate 3.7015.505 trial Plantilla para placa 3.7015.505	L	43.7015.505
	Platte 3.7029.605 Probeimplantat Plate 3.7029.605 trial Plantilla para placa 3.7029.605	L	43.7029.605
	Platte 3.7030.605 Probeimplantat Plate 3.7030.605 trial Plantilla para placa 3.7030.605	R	43.7030.605
	Platte 3.7033.004 Probeimplantat Plate 3.7033.004 trial Plantilla para placa 3.7033.004	L/R	43.7033.004
	Platte 3.7203.606 Probeimplantat Plate 3.7203.606 trial Plantilla para placa 3.7203.606	L	43.7203.606
	Platte 3.7202.606 Probeimplantat Plate 3.7202.606 trial Plantilla para placa 3.7202.606	R	43.7202.606
	Platte 3.7243.608 Probeimplantat Plate 3.7243.608 trial Plantilla para placa 3.7243.608	L	43.7243.608
	Platte 3.7244.608 Probeimplantat Plate 3.7244.608 trial Plantilla para placa 3.7244.608	R	43.7244.608


	Name Name Nombre	Links/Rechts Left/Right Izq./Der.	Catalogue No. Kat.-Nr. Número de catálogo
	Platte 3.7034.004 Probeimplantat Plate 3.7034.004 trial Plantilla para placa 3.7034.004	L/R	43.7034.004
	Platte 3.7035.004 Probeimplantat Plate 3.7035.004 trial Plantilla para placa 3.7035.004	L/R	43.7035.004
	Platte 3.7036.606 Probeimplantat Plate 3.7036.606 trial Plantilla para placa 3.7036.606	R	43.7036.606
	Platte 3.7037.606 Probeimplantat Plate 3.7037.606 trial Plantilla para placa 3.7037.606	L	43.7037.606
	Platte 3.7038.606 Probeimplantat Plate 3.7038.606 trial Plantilla para placa 3.7038.606	R	43.7038.606
	Platte 3.7039.606 Probeimplantat Plate 3.7039.606 trial Plantilla para placa 3.7039.606	L	43.7039.606
	Platte 3.7040.606 Probeimplantat Plate 3.7040.606 trial Plantilla para placa 3.7040.606	R	43.7040.606
	Platte 3.7041.606 Probeimplantat Plate 3.7041.606 Plantilla para placa 3.7041.606	L	43.7041.606
	Platte 3.7048.508 Probeimplantat Plate 3.7048.508 trial Plantilla para placa 3.7048.508	R	43.7048.508
	Platte 3.7049.508 Probeimplantat Plate 3.7049.508 trial Plantilla para placa 3.7049.508	L	43.7049.508
	Platte 3.7071.606 Probeimplantat Plate 3.7071.606 trial Plantilla para placa 3.7071.606	L	43.7071.606
	Platte 3.7072.606 Probeimplantat Plate 3.7072.606 trial Plantilla para placa 3.7072.606	R	43.7072.606
	Platte 3.7073.606 Probeimplantat Plate 3.7073.606 trial Plantilla para placa 3.7073.606	L	43.7073.606
	Platte 3.7074.606 Probeimplantat Plate 3.7074.606 trial Plantilla para placa 3.7074.606	R	43.7074.606



PLATTE 3.4025 ÷ 3.4028 PROBEIMPLANTATE
PLATE 3.4025 ÷ 3.4028 TRIAL
PLANTILLA PARA PLACA 3.4025 ÷ 3.4028



			
12	43.4025.505	43.4026.505	
15	43.4027.505	43.4028.505	

	Name Name Nombre	Links/Rechts Left/Right Izq./Der.	Catalogue No. Kat.-Nr. Número de catálogo
	Platte 3.4009.606 Probeimplantat Plate 3.4009.606 trial Plantilla para placa 3.4009.606	L	43.4009.606
	Platte 3.4010.606 Probeimplantat Plate 3.4010.606 trial Plantilla para placa 3.4010.606	R	43.4010.606
	Platte 3.4013.605 Probeimplantat Plate 3.4013.605 trial Plantilla para placa 3.4013.605	L	43.4013.605
	Platte 3.4014.605 Probeimplantat Plate 3.4014.605 trial Plantilla para placa 3.4014.605	R	43.4014.605
	Platte 3.4023.606 Probeimplantat Plate 3.4023.606 trial Plantilla para placa 3.4023.606	L	43.4023.606
	Platte 3.4024.606 Probeimplantat Plate 3.4024.606 trial Plantilla para placa 3.4024.606	R	43.4024.606
	Platte 3.4089.604 Probeimplantat Plate 3.4089.604 trial Plantilla para placa 3.4089.604	L	43.4089.604
	Platte 3.4090.604 Probeimplantat Plate 3.4090.604 trial Plantilla para placa 3.4090.604	R	43.4090.604
	Platte 3.7022.604/3.7023.604 Probeimplantat Plate 3.7022.604/3.7023.604 trial Plantilla para placa 3.7022.604/3.7023.604	L/R	43.7022.604



	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Raspatorium 6mm, 180mm Raspator 6 mm, 180 mm Periostótomo 6mm, 180mm	30.4060
	Raspatorium 6mm, 190mm Raspator 6 mm, 190 mm Periostótomo 6mm, 190mm	30.4063
	Raspatorium 13mm, 200mm Raspator 13 mm, 200 mm Periostótomo 13mm, 200mm	30.4064
	Raspatorium 13mm, 200mm Raspator 13 mm, 200 mm Periostótomo 13mm, 200mm	30.4065
	Buck-Gramcko Hebel 7,5 BUCK-GRAMCKO elevator 7.5 mm Elevador BUCK-GRAMCKO 7,5	40.2185
	Mini-Hohmann Hebel 6 MINI-HOHMANN elevator 6 Elevador MINI-HOHMANN 6	40.2186
	Mini-Hohmann Hebel 8 MINI-HOHMANN elevator 8 Elevador MINI-HOHMANN 8	40.2187
	Hebel 15 Elevator 15 Elevador 15	40.2188
	Hohmann-Hebel 17 HOHMANN elevator 17 Elevador HOHMANN 17	40.2195
	Repositionszange 120mm Bone holding forceps 120 mm Pinzas para huesos 120mm	40.4149
	Repositionszange 135mm Bone holding forceps 135 mm Pinzas para huesos 135mm	40.4150

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Repositionszange 140mm Bone holding forceps 140 mm Pinzas para huesos 140mm	40.4153
	Knochenrepositionszange 180mm Bone holding forceps 180 mm Pinzas para huesos 180mm	40.7262
	Knochenrepositionszange 170mm Bone holding forceps 170 mm Pinzas para huesos 170mm	40.7261
	Knochenrepositionszange 205mm Bone holding forceps 205 mm Pinzas para huesos 205mm	40.7263
	Repositionszange 210mm Bone holding forceps 210 mm Pinzas para huesos 210mm	40.4152

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Verbrugge-Baby Repositionsange 175mm VERBRUGGE-BABY bone holding forceps 175 mm Pinzas para huesos VERBRUGGE-BABY 175mm	40.4228.100
	Verbrugge-Baby Repositionsange 175mm VERBRUGGE-BABY bone holding forceps 175 mm Pinzas para huesos VERBRUGGE-BABY 175mm	40.4228
	Verbrugge Repositionsange 250mm VERBRUGGE bone holding forceps 250 mm Pinzas para huesos VERBRUGGE 250mm	40.4225
	Verbrugge Repositionsange 260mm VERBRUGGE bone holding forceps 260 mm Pinzas para huesos VERBRUGGE 260mm	40.4226
	Verbrugge Repositionsange 2,5x150mm VERBRUGGE bone holding forceps 2,5x150mm Pinzas para huesos VERBRUGGE 2,5x150mm	40.4229
	Verbrugge Repositionsange 6,5x190mm VERBRUGGE bone holding forceps 6,5x190 mm Pinzas para huesos VERBRUGGE 6,5x190mm	40.4230
	Verbrugge Repositionsange 9,5x240mm VERBRUGGE bone holding forceps 9,5x240 mm Pinzas para huesos VERBRUGGE 9,5x240mm	40.4231
	Repositionszange 160mm Bone holding forceps 160 mm Pinzas para huesos 160mm	40.4157
	Repositionszange 240mm Bone holding forceps 240 mm Pinzas para huesos 240mm	40.4158
	Repositionszange Nr. 0 Bone holding forceps no. 0 Pinzas para huesos n.º 0	40.0100
	Repositionszange Nr. 1 Bone holding forceps no. 1 Pinzas para huesos n.º 1	40.0101
	Repositionszange Nr. 2 Bone holding forceps no. 2 Pinzas para huesos n.º 2	40.0102
	Repositionszange Nr. 3 Bone holding forceps no. 3 Pinzas para huesos n.º 3	40.0103

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Luer Hohlmeisselzange 150mm LUER bone rongeur 150 mm Pinza gubia LUER 150mm	40.3820
	Ruskin Hohlmeisselzange 180/4mm RUSKIN bone rongeur 180/4 mm Pinza gubia RUSKIN 180/4mm	40.3318
	Liston Mini Knochenzange 150mm Mini LISTON bone cutting forceps 150 mm Pinzas para cortar hueso LISTON Mini 150mm	40.3900
	Drahtbiegezange, scharfe Spitze 14cm Wire bending forceps with sharp ends 14 cm Pinzas para doblar alambre de punta aguda 14cm	40.4444
	Drahtbiegezange 17cm Wire bending pliers 17 cm Alicates para doblar alambre 17cm	40.3185

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Inge Sperrer 6 mm, 17cm INGE reamer 6mm, 17 cm Separador INGE 6mm, 17cm	40.1392
	Inge Sperrer 6mm, 17cm INGE reamer 6mm, 17 cm Escariador INGE 6mm, 17cm	40.1394
	Lambotte Mini-Osteotom 6x125mm LAMBOTTE osteotome mini 6x125 mm Osteótomo LAMBOTTE mini 6x125mm	40.0223
	Lambotte Osteotom 13x250mm LAMBOTTE osteotome 13x250 mm Osteótomo LAMBOTTE 13x250mm	40.0213
	 Volkman Knochenkürette 17cm, 00 VOLKMANN bone curette 17 cm, 00 Cucharilla para hueso VOLKMANN 17cm, 00	40.3985
	 Volkman Knochenkürette 17cm, 0 VOLKMANN bone curette 17 cm, 0 Cucharilla para hueso VOLKMANN 17cm, 0	40.3986
	Haken Retractor Gancho	40.3599
	Schraubengreifer 3,5mm - 6,5mm Holder for 3.5 mm - 6.5 mm screws Sujetador para tornillos 3,5mm - 6,5mm	40.7271
	Schraubengreifer 1,5mm - 2,7mm Holder for 1.5 mm - 2.7 mm screws Sujetador para tornillos 1,5mm - 2,7mm	40.7272

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Tablett aus Stahldraht 255x250x50mm Steel wire tray 255x250x50mm Bandeja de alambre de acero 255x250x50mm	12.0321.050
	Tablett aus Stahldraht 255x250x70mm Steel wire tray 255x250x70 mm Bandeja de alambre de acero 255x250x70mm	12.0321.070
	Tablett aus Stahldraht 540x250x50mm Steel wire tray 540x250x50 mm Bandeja de alambre de acero 540x250x50mm	12.0318.050
	Tablett aus Stahldraht 540x250x70mm Steel wire tray 540x250x70mm Bandeja de alambre de acero 540x250x70mm	12.0318.070
	Tablett aus Stahldraht 480x250x50mm Steel wire tray 480x250x50mm Bandeja de alambre de acero 480x250x50mm	12.0320.050
	Tablett aus Stahldraht 480x250x70mm Steel wire tray 480x250x70mm Bandeja de alambre de acero 480x250x70mm	12.0320.070
	Silikonmatte 250x235mm Silicone mat 250x235mm Colchoneta de silicona 250x235mm	12.0719
	Silikonmatte 470x230mm Silicone mat 470x230 mm Colchoneta de silicona 470x230mm	12.0725
	Silikonmatte 520x230mm Silicone mat 520x230mm Colchoneta de silicona 520x230mm	12.0723
	Perforierter Aluminiumdeckel, grau 1/2 306x272x15mm Perforated aluminum cover 1/2 306x272x15mm Gray Tapa de aluminio perforada 1/2 306x272x15mm Gris	12.0751.200
	Behälter mit solidem Boden 1/2 306x272x85mm Container with solid bottom 1/2 306x272x85mm Contenedor de fondo sólido 1/2 306x272x85mm	12.0751.100

	Name Name Nombre	Kat.-Nr., Catalogue no. Número de catálogo
	Perforierter Aluminiumdeckel, grau, 1/1 595x275x15mm Perforated aluminum lid 1/1 595x275x15mm Gray Tapa de aluminio perforada 1/1 595x275x15mm Gris	12.0750.200
	Behälter mit solidem Boden 1/1 595x275x86mm Container with solid bottom 1/1 595x275x86mm Contenedor de fondo sólido 1/1 595x275x86mm	12.0750.100

Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página
3.1220 ... 9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 28, 31, 38, 39, 40, 41, 56, 58, 81, 101, 103, 144		3.4041.60315, 155		36, 37, 38, 39, 40, 41, 44, 45, 46, 47, 49, 50, 56, 57		3.7027.50227, 156	
3.1221.050 ..59, 60, 61, 62, 63, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 101, 105		3.404213, 15		3.5200 ..59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 101, 102, 104		3.702828	
3.1221.070 ..118, 119, 120, 121, 122, 123, 126, 127, 128, 132, 133, 136, 137, 138, 144, 147, 151		3.4042.60315, 155		3.5205 ..59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 101		3.7028.50228, 156	
3.1306 ..56, 59, 60, 61, 63, 64, 66, 67, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 101, 103		3.404313, 19		3.5207.006 ..59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 101, 105		3.702992, 117	
3.1379144		3.4043.60319, 155		3.5210 .118, 119, 120, 121, 122, 123, 124, 125, 126, 127, 128, 132, 133, 136, 137, 138, 142, 143, 144, 147, 148		3.7029.60592, 157	
3.1443124, 125		3.404413, 19		3.5216 .118, 119, 120, 121, 122, 123, 126, 127, 128, 142, 143, 144, 147, 150		3.703092	
3.1471 .118, 119, 120, 121, 122, 123, 126, 127, 128, 133, 136, 137, 138, 144, 147, 148		3.4044.60319, 155		3.5219151		3.7030.60592, 157	
3.3155118		3.404513, 17		3.5219.008 ..118, 119, 120, 121, 122, 123, 126, 127, 128, 132, 133, 136, 137, 138, 144, 147		3.703381	
3.315761, 106		3.4045.60317, 155		3.5221 .119, 126, 142, 143, 144, 147, 151		3.7033.00481, 112, 157	
3.3158106		3.404613, 17		3.5222 ..119, 126, 142, 143, 147, 149		3.703481	
3.3159106		3.4046.60317, 155		3.5225 ..59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 101, 102		3.7034.00481, 112, 158	
3.3538144		3.404713, 14		3.5228 .119, 126, 142, 143, 144, 147, 150		3.703581	
3.4008117		3.4047.60314, 155		3.5314132, 133		3.7035.00481, 112, 158	
3.4009121		3.404813, 18		3.5316132, 133, 151		3.703674, 111	
3.4010121		3.4048.60318, 155		3.701467, 107		3.7036.606158	
3.401185, 114		3.405191, 116, 157		3.7014.50567, 157		3.703774, 110	
3.4011.605157		3.405291, 116, 157		3.701567, 107		3.7037.60674, 158	
3.401285, 114		3.405794		3.7015.50567, 157		3.703875, 111	
3.4012.605157		3.405894		3.702044		3.7038.606158	
3.4013120		3.405995		3.7021.04130		3.703975, 110	
3.4014120		3.406095		3.7022126		3.7039.606158	
3.402260, 106		3.406131		3.7023126		3.704076, 111	
3.4023119		3.406247		3.702427		3.7040.606158	
3.4024119		3.406332		3.7024.50227, 156		3.704176, 110	
3.402566, 107, 158		3.406432		3.702527		3.7041.606158	
3.402666, 107		3.406546		3.702627		3.704868, 108	
3.402766		3.406633		3.7026.50227, 156		3.7048.50868, 158	
3.402866		3.406734		3.702727		3.7049108	
3.4029107		3.406834		3.7027.50227, 156		3.7049.50868, 158	
3.4030107		3.406935		3.702827		3.7054123	
3.403145		3.407036		3.7028.50227, 156		3.7055123	
3.4031.030155		3.407136		3.702927		3.705649	
3.403245		3.407241		3.703027		3.705750	
3.4032.030155		3.407333		3.7030.60592, 157		3.706052	
3.403389, 115		3.407440		3.703381		3.706152	
3.403472, 109		3.407540		3.703481		3.7062127	
3.403713, 16		3.408059, 106		3.7034.00481, 112, 158		3.7065142, 143	
3.4037.60316, 155		3.408162, 63		3.703581		3.7066142, 143	
3.403813, 16		3.408439		3.703674, 111		3.706782, 113	
3.4038.60316, 155		3.408537		3.7036.606158		3.706882, 113	
3.403990, 116		3.408637		3.703774, 110		3.707177, 110	
3.4039.506157		3.408738		3.7037.60674, 158		3.7071.606158	
3.404090, 116		3.408838		3.703875, 111		3.707277, 111	
3.4040.506157		3.4089122		3.7038.606158		3.7072.606158	
3.404113, 15		3.4090122		3.703975, 110		3.707378, 110	
		3.409780, 112		3.7039.606158		3.707478, 111	
		3.409880, 112		3.704076, 111		3.7074.606158	
		3.4098.50480, 157		3.7040.606158		3.707512	
		3.409920		3.704176, 110		3.707612	
		3.4159138		3.7041.606158		3.707712	
		3.4160138		3.704868, 108		3.707812	
		3.51649, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 24, 28, 30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 44, 45, 46, 47, 49, 50, 52, 56, 57		3.7048.50868, 158		3.707912	
		3.51659, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 28, 30, 31, 32, 33, 34, 35,		3.7049108		3.708012	
				3.7049.50868, 158		3.708112	
				3.7054123		3.708212	
				3.7055123		3.708311	
				3.705649			
				3.705750			
				3.706052			
				3.706152			
				3.7062127			
				3.7065142, 143			
				3.7066142, 143			
				3.706782, 113			
				3.706882, 113			
				3.707177, 110			
				3.7071.606158			
				3.707277, 111			
				3.7072.606158			
				3.707378, 110			
				3.7073.606158			
				3.707478, 111			
				3.7074.606158			
				3.707512			
				3.707612			
				3.707712			
				3.707812			
				3.707912			
				3.708012			
				3.708112			
				3.708212			
				3.708311			

Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página
3.7084	11	12.0321.050	165	14.0207.601	139	40.4146	48
3.7085	11	12.0321.070	165	15.0204.001	53	40.4146.000	51
3.7086	11	12.0719	165	15.0204.201	53, 54	40.4149	160
3.7091	83, 113	12.0723	165	15.0204.202	53, 55	40.4150	160
3.7092	83, 113	12.0725	165	15.0204.602	25	40.4152	161
3.7093	128	12.0750.100	146, 152, 153, 154, 166	15.0205.106	141	40.4153	161
3.7094	86, 115	12.0750.200	146, 152, 153, 166	15.0205.203	97	40.4157	162
3.7095	87, 115	12.0751.100	21, 22, 26, 29, 42, 48, 51, 165	15.0205.208	96	40.4158	162
3.7098	136	12.0751.103	154	15.0205.401	109	40.4225	162
3.7099	137	12.0751.200	21, 22, 26, 29, 42, 48, 51, 152, 153, 154, 165	15.0205.409	112	40.4226	162
3.7200	93, 117	14.0204.103	53	15.0205.411	113	40.4228	162
3.7201	93, 117	14.0204.104	53	15.0205.415	116	40.4228.100	162
3.7202	84, 114	14.0204.201	54	15.0205.419	110	40.4229	162
3.7203	84, 114	14.0204.601	53	15.0205.420	111	40.4230	162
3.7206.014	24	14.0204.602	25	15.0205.421	113	40.4231	162
3.7206.021	24	14.0205.106	141	15.0205.422	115	40.4444	163
3.7207	69	14.0205.107	141	15.0205.424	84, 114	40.4592.130	141
3.7208	69	14.0205.203	97	15.0205.427	114	40.4592.180	25
3.7209	70	14.0205.208	96	15.0205.430	65	40.4592.210	96
3.7210	70	14.0205.401	72, 109	15.0205.601	65	40.4639.550	96, 130, 140, 141
3.7211	71	14.0205.402	60, 106	15.0207.105	140	40.4639.700	145
3.7212.600	71	14.0205.403	66, 107	15.0207.208	129	40.4640.000	23, 26, 29, 43, 48, 51, 54
3.7220	124	14.0205.404	61, 106	30.3303.000	23, 26, 29, 43, 48, 51	40.4643.000	23, 26, 29, 43
3.7221	124	14.0205.405	67, 107	30.4060	160	40.4746.000	151
3.7222	124	14.0205.406	68, 108	30.4063	160	40.4802.732	140, 145, 148
3.7223	124	14.0205.409	80, 112	30.4064	160	40.4804.700	103
3.7224	125	14.0205.410	81, 112	30.4065	160	40.4804.725	141
3.7225	125	14.0205.411	82, 113	40.0100	162	40.4814.000	23, 26, 29, 43, 48, 51
3.7226	9	14.0205.414	89, 115	40.0101	162	40.4815.210	129, 140, 151
3.7227	9	14.0205.415	90, 116	40.0102	162	40.4815.220	145, 149
3.7228	9	14.0205.417	92, 117	40.0103	162	40.4818.000	43
3.7229	9	14.0205.419	74, 75, 76, 77, 78, 110	40.0213	164	40.4818.100	23, 26, 29, 54
3.7230	9	14.0205.420	74, 75, 76, 77, 78, 111	40.0223	164	40.4896.018	23, 26, 29, 43, 48, 51, 54, 57
3.7231	9	14.0205.421	83, 113	40.0669.000	23, 48, 58, 103	40.4897.018	23, 26, 29, 43, 58, 103
3.7232	10	14.0205.422	86, 87, 115	40.0669.100	26, 29, 43, 51	40.4898.070	135
3.7233	10	14.0205.423	93, 117	40.1392	164	40.5350	142, 144
3.7234	10	14.0205.424	114	40.1394	164	40.5350.600	142, 143, 144, 145
3.7235	10	14.0205.426	69, 70, 71, 108	40.2063.181	23, 26, 29, 43, 48, 51, 54, 57, 58, 103	40.5353	145
3.7236	10	14.0205.427	85, 114	40.2063.212	96	40.5357.016	145
3.7237	10	14.0205.428	91, 116	40.2063.222	102	40.5357.020	145
3.7238	88	14.0205.429	73, 109	40.2185	160	40.5358	146
3.7240	79, 80	14.0205.430	65	40.2185.000	54	40.5359.600	146
3.7242	64	14.0205.501	98	40.2186	160	40.5360	145
3.7243	73, 109	14.0205.502	98	40.2187	160	40.5361	145
3.7244	73, 109	14.0205.503	99	40.2188	160	40.5362	145
4.5235	9, 10, 11, 12, 14, 15, 16, 17, 18, 19, 20, 24, 28, 44, 45, 46, 47, 49, 50, 52, 58	14.0205.504	100	40.2195	160	40.5369	145
4.5236	59, 60, 61, 62, 63, 64, 66, 67, 68, 69, 70, 71, 72, 73, 74, 75, 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93, 94, 95, 105	14.0205.506	98	40.3185	163	40.5612.000	135
12.0318.050	165	14.0205.601	65	40.3318	163	40.5646.000	135
12.0318.070	165	14.0205.602	139	40.3599	164	40.5647.000	135
12.0320.050	165	14.0207.105	140	40.3820	163	40.5650.212	129, 140, 150
12.0320.070	165	14.0207.106	140	40.3900	163	40.5650.222	145, 148, 149
		14.0207.208	129	40.3943.100	145	40.5650.301	135
		14.0207.501	131	40.3974.200	145	40.5651.212	129, 140, 150
		14.0207.502	131	40.3985	164	40.5651.222	145, 148
				40.3986	164		

Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página
40.5651.301	135	40.5707.732	146, 150	40.6297.000	51	40.8224.000	86
40.5652.210	151	40.5708.000	129, 152, 153, 154	40.6298.000	49, 50, 51	40.8225	115
40.5652.212	129	40.5709.100	153	40.6299.000	153	40.8225.000	87
40.5652.222	149	40.5709.200	153	40.6299.100	124, 125, 153	40.8226	114
40.5653.212	65, 96, 104, 105	40.5710.020	42	40.6600.000	134	40.8226.000	84
40.5653.222	102	40.5710.141	21	40.6601.000	134	40.8227.000	84, 114
40.5654.120	129, 149, 151	40.5711.000	42, 43	40.6602.000	134	40.8232.000	85, 114
40.5654.750	129, 149, 151	40.5711.200	21, 23, 26	40.6603.000	134	40.8233.000	85, 114
40.5661.000	105	40.5711.300	26	40.6604.000	134	40.8234.000	91, 116
40.5671	109	40.5712.000	43	40.6605.000	134	40.8235.000	91, 116
40.5671.000	72	40.5712.100	23, 26	40.6606.000	134	40.8236	9, 10
40.5672.000	96	40.5713.040	14, 15, 16, 17, 18, 19	40.6607.000	134	40.8237	9, 10
40.5673.718	96, 102	40.5713.050	22	40.6608.000	134	40.8238	9, 10
40.5673.725	104	40.5713.150	11, 12, 22	40.6609.000	134	40.8239	9, 10
40.5673.728	65, 96, 102, 104	40.5713.170	9, 10, 21	40.6610.000	134	40.8240	10
40.5675.100	146	40.5713.230	29	40.6611.000	135	40.8240.	9
40.5677.000	65, 96, 102, 103, 104, 105, 141	40.5713.330	27, 28, 29	40.6612.000	134	40.8241	9, 10
40.5678	23	40.5714.020	30, 31, 32, 33, 34, 35, 36, 37, 38, 39, 40, 41, 42	40.6619.000	135	40.8262.015	25
40.5678.000	26	40.5724.100	120, 152	40.6650.000	23, 26, 29, 43, 48, 51, 54	40.8262.022	25
40.5682.000	23, 26, 29, 43, 48, 51, 54, 57, 58	40.5724.200	120, 152	40.6652.000	96	40.8263.000	65
40.5684.000	135, 146, 148, 149	40.5725.100	119, 152	40.6654.000	65, 141	43.4009.606	121, 159
40.5684.200	130, 140, 150	40.5725.200	119, 152	40.6660.000	130, 140	43.4010.606	121, 159
40.5685.000	146, 149	40.5726.100	90, 116	40.7261	161	43.4011.605	85, 157
40.5685.200	130, 150, 151	40.5726.200	90, 116	40.7262	161	43.4012.605	85, 157
40.5686.000	146	40.5731.100	80, 112	40.7263	161	43.4013.605	120, 159
40.5687.000	146	40.5731.200	80, 112	40.7271	164	43.4014.605	120, 159
40.5690.540	135	40.5732.100	126, 154	40.7272	164	43.4023.606	119, 159
40.5691.532	135	40.5732.200	126, 154	40.8206.028	96, 105	43.4024.606	119, 159
40.5693.570	135	40.5741.200	48	40.8207.040	129	43.4025.505	66, 158
40.5694.580	135	40.5741.300	44, 45, 46, 47, 48	40.8210	111	43.4026.505	66, 158
40.5695.570	134	40.5742.100	82, 113	40.8210.000	74	43.4027.505	66, 158
40.5696.580	135	40.5742.200	82, 113	40.8211	110	43.4028.505	66, 158
40.5700.000	135	40.5749.600	154	40.8211.000	74	43.4031.030	155
40.5704.310	152	40.5749.700	154	40.8212	111	43.4032.030	155
40.5704.320	152	40.5912.212	104, 141	40.8212.000	75	43.4037.603	16, 155
40.5704.340	153	40.5912.222	103	40.8213	110	43.4038.603	16, 155
40.5704.350	152	40.5927.020	96	40.8213.000	75	43.4039.506	157
40.5704.360	153	40.5927.040	130	40.8214	111	43.4039.608	90
40.5704.390	154	40.5928.018	23, 26, 29, 48, 51, 54, 58	40.8214.000	76	43.4040.506	157
40.5704.410	152	40.5987.024	23, 26	40.8215	110	43.4040.608	90
40.5704.420	152	40.5988.000	23, 26	40.8215.000	76	43.4041.603	15, 155
40.5704.460	153	40.5989.000	23, 26, 54	40.8216	111	43.4042.603	15, 155
40.5704.490	154	40.5990.212	97	40.8216.000	77	43.4043.603	19, 155
40.5704.510	119, 152	40.5992.000	97	40.8217	110	43.4044.603	19, 155
40.5704.520	120, 152	40.5994.715	97	40.8217.000	77	43.4045.603	17, 155
40.5704.540	121, 153	40.5994.732	97	40.8218	111	43.4046.603	17, 155
40.5704.550	118, 152	40.6281.000	55	40.8218.000	78	43.4047.603	14, 155
40.5704.560	122, 123, 153	40.6282.000	55	40.8219	110	43.4048.603	18, 155
40.5704.590	126, 153, 154	40.6283.000	54	40.8219.000	78	43.4051	157
40.5704.800	154	40.6284.000	54	40.8220.000	123, 153	43.4051.611	91
40.5704.900	127, 128, 154	40.6285.000	55	40.8221.000	123, 153	43.4052	157
40.5705.732	129, 149	40.6293.000	153	40.8222	113	43.4052.611	91
40.5705.740	129, 140, 146, 148, 150	40.6296.000	154	40.8223	113	43.4089.604	122, 159
40.5706.740	129	40.6296.100	154	40.8223.000	83	43.4090.604	122, 159
				40.8224	115	43.4097.504	80, 157
						43.4098.504	80, 157

Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página	Kat.-Nr., Catalogue no. Número de catálogo	Seite Page Página
43.7014.505	67, 157						
43.7015.505	67, 157						
43.7022.604	126, 159						
43.7024.502	27, 156						
43.7025.502	27, 156						
43.7026.502	27, 156						
43.7027.502	27, 156						
43.7028.502	28, 156						
43.7029.605	92, 157						
43.7030.605	92, 157						
43.7033.004	81, 157						
43.7034.004	81, 158						
43.7035.004	81, 158						
43.7036.606	74, 158						
43.7037.606	74, 158						
43.7038.606	75, 158						
43.7039.606	75						
43.7040.606	76, 158						
43.7041.606	76						
43.7048.508	68, 158						
43.7049.508	68, 158						
43.7071.606	77, 158						
43.7072.606	77, 158						
43.7073.606	78, 158						
43.7074.606	78						
43.7202.606	84, 157						
43.7203.606	84, 157						
43.7228.606	9, 156						
43.7229.606	9, 156						
43.7234.606	10, 156						
43.7235.606	156						
43.7243.608	73, 157						
43.7244.608	73, 157						

ChM Ltd.
Lewickie 3b
16-061 Juchnowiec Kościelny
POLAND
tel.: +48 85 86 86 100
fax: +48 85 86 86 101
www.chm.eu